



COSIL
INSTRUMENTS
COMPANY

MANUFACTURAR & EXPORTER OF ALL
KIND OF ELECTRO SURGICAL INSTRUMENTS



ELECTROSURGICAL INSULATED BIPOLAR FORCEPS



ELECTROSURGICAL INSULATED BIPOLAR FORCEPS

Our range of ElectroSurgery Forceps is designed to provide precision and control during electroSurgical procedures. Featuring highly polished tips for optimal heat and electrical conductivity, these forceps ensure rapid heat dispersion, allowing for effective cutting and coagulation. Available in both standard and advanced models, our forceps are built to minimize tissue adhesion and enhance the surgeon's tactile feedback.



The middle area of non-stick forceps often features a textured or grooved surface to provide a secure and controlled grip for the user. This section is usually reinforced for added stability and balance, allowing for precise manipulation even during delicate procedures.



The tips of non-stick forceps are finely crafted for precision and coated with a non-stick material, such as Teflon, to prevent adhesion and ensure smooth handling. These tips are typically narrow and tapered, allowing for delicate manipulation of tissues or materials without sticking or tearing.



The Nonstick Forceps Plugin is an essential accessory for medical and laboratory equipment, designed to

Whether for single-use or reusable options, these instruments offer reliability and efficiency, ensuring safe and effective surgical outcomes.



Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM101A	CEFM101B	CEFM101C	CEFM101D	CEFM101E	0.50mm
	CEFM102A	CEFM102B	CEFM102C	CEFM102D	CEFM102E	1.00mm
	CEFM103A	CEFM103B	CEFM103C	CEFM103D	CEFM103E	2.00mm
Single-Use	CEFS101A	CEFS101B	CEFS101C	CEFS101D	CEFS101E	0.50mm
	CEFS102A	CEFS102B	CEFS102C	CEFS102D	CEFS102E	1.00mm
	CEFS103A	CEFS103B	CEFS103C	CEFS103D	CEFS103E	2.00mm

Box of 20 Units



Insulated Straight Angle Tip Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM104A	CEFM104B	CEFM104C	CEFM104D	CEFM104E	0.50mm
	CEFM105A	CEFM105B	CEFM105C	CEFM105D	CEFM105E	1.00mm
	CEFM106A	CEFM106B	CEFM106C	CEFM106D	CEFM106E	2.00mm
Single-Use	CEFS104A	CEFS104B	CEFS104C	CEFS104D	CEFS104E	0.50mm
	CEFS105A	CEFS105B	CEFS105C	CEFS105D	CEFS105E	1.00mm
	CEFS106A	CEFS106B	CEFS106C	CEFS106D	CEFS106E	2.00mm

Box of 20 Units



Insulated Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM107A	CEFM107B	CEFM107C	CEFM107D	CEFM107E	0.50mm
	CEFM108A	CEFM108B	CEFM108C	CEFM108D	CEFM108E	1.00mm
	CEFM109A	CEFM109B	CEFM109C	CEFM109D	CEFM109E	2.00mm
Single-Use	CEFS107A	CEFS107B	CEFS107C	CEFS107D	CEFS107E	0.50mm
	CEFS108A	CEFS108B	CEFS108C	CEFS108D	CEFS108E	1.00mm
	CEFS109A	CEFS109B	CEFS109C	CEFS109D	CEFS109E	2.00mm

Box of 20 Units



Insulated Bayonet Tip Down Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM110A	CEFM110B	CEFM110C	CEFM110D	CEFM110E	0.50mm
	CEFM111A	CEFM111B	CEFM111C	CEFM111D	CEFM111E	1.00mm
	CEFM112A	CEFM112B	CEFM112C	CEFM112D	CEFM112E	2.00mm
Single-Use	CEFS110A	CEFS110B	CEFS110C	CEFS110D	CEFS110E	0.50mm
	CEFS111A	CEFS111B	CEFS111C	CEFS111D	CEFS111E	1.00mm
	CEFS112A	CEFS112B	CEFS112C	CEFS112D	CEFS112E	2.00mm

Box of 20 Units



Insulated Bayonet Tip Up Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM113A	CEFM113B	CEFM113C	CEFM113D	CEFM113E	0.50mm
	CEFM114A	CEFM114B	CEFM114C	CEFM114D	CEFM114E	1.00mm
	CEFM115A	CEFM115B	CEFM115C	CEFM115D	CEFM115E	2.00mm
Single-Use	CEFS113A	CEFS113B	CEFS113C	CEFS113D	CEFS113E	0.50mm
	CEFS114A	CEFS114B	CEFS114C	CEFS114D	CEFS114E	1.00mm
	CEFS115A	CEFS115B	CEFS115C	CEFS115D	CEFS115E	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Bipolar Forceps

Flat Plug Connection



Mchpherson Bipolar Forceps

102mm	
Multi-Use	CEFM116T 0.50mm
	CEFM117T 1.00mm
	CEFM118T 2.00mm
Single-Use	CEFS116T 0.50mm
	CEFS117T 1.00mm
	CEFS118T 2.00mm

Box of 20 Units



Jeweler#7 Bipolar Forceps

115mm	
Multi-Use	CEFM119U 0.50mm
	CEFM120U 1.00mm
	CEFM121U 2.00mm
Single-Use	CEFS119U 0.50mm
	CEFS120U 1.00mm
	CEFS121U 2.00mm

Box of 20 Units



Jeweler#4 Ang. Tip Bipolar Forceps

108mm	
Multi-Use	CEFM122U 0.50mm
	CEFM123U 1.00mm
	CEFM124U 2.00mm
Single-Use	CEFS122U 0.50mm
	CEFS123U 1.00mm
	CEFS124U 2.00mm

Box of 20 Units



Jeweler#4 Str. Tip Bipolar Forceps

115mm	
Multi-Use	CEFM125V 0.50mm
	CEFM126V 1.00mm
	CEFM127V 2.00mm
Single-Use	CEFS125V 0.50mm
	CEFS126V 1.00mm
	CEFS127V 2.00mm

Box of 20 Units



Adson Bipolar Forceps

121mm		155mm		180mm		200mm	
Multi-Use	CEFM128W	CEFM128A	CEFM128B	CEFM128C	CEFM129W	CEFM129A	CEFM129B
	CEFM129W	CEFM129A	CEFM129B	CEFM129C	CEFM130W	CEFM130A	CEFM130B
	CEFM130W	CEFM130A	CEFM130B	CEFM130C			
Single-Use	CEFS128W	CEFS128A	CEFS128B	CEFS128C	CEFS129W	CEFS129A	CEFS129B
	CEFS129W	CEFS129A	CEFS129B	CEFS129C	CEFS130W	CEFS130A	CEFS130B
	CEFS130W	CEFS130A	CEFS130B	CEFS130C			

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available

Tip Sizes	0.25mm
	0.75mm
	1.25mm
	2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
 Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Bipolar Forceps

Two Pin Connection



Iris Str. Bipolar Forceps

108mm		
Multi-Use	CEFM167U	0.50mm
	CEFM168U	1.00mm
	CEFM169U	2.00mm
Single-Use	CEFS167U	0.50mm
	CEFS168U	1.00mm
	CEFS169U	2.00mm

Box of 20 Units



Iris Ang. Bipolar Forceps

108mm		
Multi-Use	CEFM170U	0.50mm
	CEFM171U	1.00mm
	CEFM172U	2.00mm
Single-Use	CEFS170U	0.50mm
	CEFS171U	1.00mm
	CEFS172U	2.00mm

Box of 20 Units



Jeweler#5 Str. Bipolar Forceps

115mm		
Multi-Use	CEFM173V	0.50mm
	CEFM174V	1.00mm
	CEFM175V	2.00mm
Single-Use	CEFS173V	0.50mm
	CEFS174V	1.00mm
	CEFS175V	2.00mm

Box of 20 Units



Jeweler#4 Str. Bipolar Forceps

115mm		
Multi-Use	CEFM176V	0.50mm
	CEFM177V	1.00mm
	CEFM178V	2.00mm
Single-Use	CEFS176V	0.50mm
	CEFS177V	1.00mm
	CEFS178V	2.00mm

Box of 20 Units



Jeweler#4 Pointed Bipolar Forceps

115mm		
Multi-Use	CEFM179V	0.50mm
	CEFM180V	1.00mm
	CEFM181V	2.00mm
Single-Use	CEFS179V	0.50mm
	CEFS180V	1.00mm
	CEFS181V	2.00mm

Box of 20 Units



Jeweler#7 Bipolar Forcep

115mm		
Multi-Use	CEFM182V	0.50mm
	CEFM183V	1.00mm
	CEFM184V	2.00mm
Single-Use	CEFS182V	0.50mm
	CEFS183V	1.00mm
	CEFS184V	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available

Tip Sizes	0.25mm
	0.75mm
	1.25mm
	2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
 Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Bipolar Forceps

Two Pin Connection



Coaptation Special Bipolar Forceps

115mm	
Multi-Use	CEFM185V 0.50mm
	CEFM186V 1.00mm
	CEFM187V 2.00mm
Single-Use	CEFS185V 0.50mm
	CEFS186V 1.00mm
	CEFS187V 2.00mm

Box of 20 Units



Jeweler#9 Bipolar Forceps

108mm	
Multi-Use	CEFM188U 0.50mm
	CEFM189U 1.00mm
	CEFM190U 2.00mm
Single-Use	CEFS188U 0.50mm
	CEFS189U 1.00mm
	CEFS190U 2.00mm

Box of 20 Units



Jeweler#4 Monopolar Forceps

115mm	
Multi-Use	CEFM191V 0.50mm
	CEFM192V 1.00mm
	CEFM193V 2.00mm
Single-Use	CEFS191V 0.50mm
	CEFS192V 1.00mm
	CEFS193V 2.00mm

Box of 20 Units



Jeweler#7 Monopolar Forceps

108mm	
Multi-Use	CEFM194U 0.50mm
	CEFM195U 1.00mm
	CEFM196U 2.00mm
Single-Use	CEFS194U 0.50mm
	CEFS195U 1.00mm
	CEFS196U 2.00mm

Box of 20 Units



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM197W	CEFM197A	CEFM197B	CEFM197C	0.50mm
	CEFM198W	CEFM198A	CEFM198B	CEFM198C	1.00mm
	CEFM199W	CEFM199A	CEFM199B	CEFM199C	2.00mm
Single-Use	CEFS197W	CEFS197A	CEFS197B	CEFS197C	0.50mm
	CEFS198W	CEFS198A	CEFS198B	CEFS198C	1.00mm
	CEFS199W	CEFS199A	CEFS199B	CEFS199C	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Regular	
Tip Sizes	0.25mm
	0.75mm
	1.25mm
	2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
 Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Bipolar Forceps

Two Pin Connection



Insulated Malis Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM131A	CEFM131B	CEFM131C	CEFM131D	CEFM131E	0.50mm
	CEFM132A	CEFM132B	CEFM132C	CEFM132D	CEFM132E	1.00mm
	CEFM133A	CEFM133B	CEFM133C	CEFM133D	CEFM133E	2.00mm
Single-Use	CEFS131A	CEFS131B	CEFS131C	CEFS131D	CEFS131E	0.50mm
	CEFS132A	CEFS132B	CEFS132C	CEFS132D	CEFS132E	1.00mm
	CEFS133A	CEFS133B	CEFS133C	CEFS133D	CEFS133E	2.00mm

Box of 20 Units



Insulated Hardy Angled Shape Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM140A	CEFM140B	CEFM140C	CEFM140D	CEFM140E	0.50mm
	CEFM141A	CEFM141B	CEFM141C	CEFM141D	CEFM141E	1.00mm
	CEFM142A	CEFM142B	CEFM142C	CEFM142D	CEFM142E	2.00mm
Single-Use	CEFS140A	CEFS140B	CEFS140C	CEFS140D	CEFS140E	0.50mm
	CEFS141A	CEFS141B	CEFS141C	CEFS141D	CEFS141E	1.00mm
	CEFS142A	CEFS142B	CEFS142C	CEFS142D	CEFS142E	2.00mm

Box of 20 Units



Insulated Jensen Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM134A	CEFM134B	CEFM134C	CEFM134D	CEFM134E	0.50mm
	CEFM135A	CEFM135B	CEFM135C	CEFM135D	CEFM135E	1.00mm
	CEFM136A	CEFM136B	CEFM136C	CEFM136D	CEFM136E	2.00mm
Single-Use	CEFS134A	CEFS134B	CEFS134C	CEFS134D	CEFS134E	0.50mm
	CEFS135A	CEFS135B	CEFS135C	CEFS135D	CEFS135E	1.00mm
	CEFS136A	CEFS136B	CEFS136C	CEFS136D	CEFS136E	2.00mm

Box of 20 Units



Insulated Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM143A	CEFM143B	CEFM143C	CEFM143D	CEFM143E	0.50mm
	CEFM144A	CEFM144B	CEFM144C	CEFM144D	CEFM144E	1.00mm
	CEFM145A	CEFM145B	CEFM145C	CEFM145D	CEFM145E	2.00mm
Single-Use	CEFS143A	CEFS143B	CEFS143C	CEFS143D	CEFS143E	0.50mm
	CEFS144A	CEFS144B	CEFS144C	CEFS144D	CEFS144E	1.00mm
	CEFS145A	CEFS145B	CEFS145C	CEFS145D	CEFS145E	2.00mm

Box of 20 Units



Insulated Hardy Bayonet Straight Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM137A	CEFM137B	CEFM137C	CEFM137D	CEFM137E	0.50mm
	CEFM138A	CEFM138B	CEFM138C	CEFM138D	CEFM138E	1.00mm
	CEFM139A	CEFM139B	CEFM139C	CEFM139D	CEFM139E	2.00mm
Single-Use	CEFS137A	CEFS137B	CEFS137C	CEFS137D	CEFS137E	0.50mm
	CEFS138A	CEFS138B	CEFS138C	CEFS138D	CEFS138E	1.00mm
	CEFS139A	CEFS139B	CEFS139C	CEFS139D	CEFS139E	2.00mm

Box of 20 Units



Insulated Cushing Crush Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM146A	CEFM146B	CEFM146C	CEFM146D	CEFM146E	0.50mm
	CEFM147A	CEFM147B	CEFM147C	CEFM147D	CEFM147E	1.00mm
	CEFM148A	CEFM148B	CEFM148C	CEFM148D	CEFM148E	2.00mm
Single-Use	CEFS146A	CEFS146B	CEFS146C	CEFS146D	CEFS146E	0.50mm
	CEFS147A	CEFS147B	CEFS147C	CEFS147D	CEFS147E	1.00mm
	CEFS148A	CEFS148B	CEFS148C	CEFS148D	CEFS148E	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
 Cosil Multi-Use/Single-Use Length



Insulated Gerald Straight Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM149A CEFM150A CEFM151A	CEFM149B CEFM150B CEFM151B	CEFM149C CEFM150C CEFM151C	CEFM149D CEFM150D CEFM151D	CEFM149E CEFM150E CEFM151E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS149A CEFS150A CEFS151A	CEFS149B CEFS150B CEFS151B	CEFS149C CEFS150C CEFS151C	CEFS149D CEFS150D CEFS151D	CEFS149E CEFS150E CEFS151E	0.50mm 1.00mm 2.00mm

Box of 20 Units



Adson Bipolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM158W CEFM159W CEFM160W	CEFM158A CEFM159A CEFM160A	CEFM158B CEFM159B CEFM160B	CEFM158C CEFM159C CEFM160C	0.50mm 1.00mm 2.00mm
Single-Use	CEFS158W CEFS159W CEFS160W	CEFS158A CEFS159A CEFS160A	CEFS158B CEFS159B CEFS160B	CEFS158C CEFS159C CEFS160C	0.50mm 1.00mm 2.00mm

Box of 20 Units



Insulated Gerald Straight Angled Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM152A CEFM153A CEFM154A	CEFM152B CEFM153B CEFM154B	CEFM152C CEFM153C CEFM154C	CEFM152D CEFM153D CEFM154D	CEFM152E CEFM153E CEFM154E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS152A CEFS153A CEFS154A	CEFS152B CEFS153B CEFS154B	CEFS152C CEFS153C CEFS154C	CEFS152D CEFS153D CEFS154D	CEFS152E CEFS153E CEFS154E	0.50mm 1.00mm 2.00mm

Box of 20 Units



Smart Iris Str. Bipolar Forceps

	102mm	
Multi-Use	CEFM161T CEFM162T CEFM163T	0.50mm 1.00mm 2.00mm
Single-Use	CEFS161T CEFS162T CEFS163T	0.50mm 1.00mm 2.00mm

Box of 20 Units



Insulated Gerald Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM155A CEFM156A CEFM157A	CEFM155B CEFM156B CEFM157B	CEFM155C CEFM156C CEFM157C	CEFM155D CEFM156D CEFM157D	CEFM155E CEFM156E CEFM157E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS155A CEFS156A CEFS157A	CEFS155B CEFS156B CEFS157B	CEFS155C CEFS156C CEFS157C	CEFS155D CEFS156D CEFS157D	CEFS155E CEFS156E CEFS157E	0.50mm 1.00mm 2.00mm

Box of 20 Units



Smart Iris Ang. Bipolar Forceps

	102mm	
Multi-Use	CEFM164T CEFM165T CEFM166T	0.50mm 1.00mm 2.00mm
Single-Use	CEFS164T CEFS165T CEFS166T	0.50mm 1.00mm 2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Irrigation Forceps

Two Pin, Flat Plug Connection



Irrigation Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM200A	CEFM200B	CEFM200C	CEFM200D	CEFM200E	0.50mm
	CEFM201A	CEFM201B	CEFM201C	CEFM201D	CEFM201E	1.00mm
	CEFM202A	CEFM202B	CEFM202C	CEFM202D	CEFM202E	2.00mm
Single-Use	CEFS200A	CEFS200B	CEFS200C	CEFS200D	CEFS200E	0.50mm
	CEFS201A	CEFS201B	CEFS201C	CEFS201D	CEFS201E	1.00mm
	CEFS202A	CEFS202B	CEFS202C	CEFS202D	CEFS202E	2.00mm

Box of 20 Units



Irrigation Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM203A	CEFM203B	CEFM203C	CEFM203D	CEFM203E	0.50mm
	CEFM204A	CEFM204B	CEFM204C	CEFM204D	CEFM204E	1.00mm
	CEFM205A	CEFM205B	CEFM205C	CEFM205D	CEFM205E	2.00mm
Single-Use	CEFS203A	CEFS203B	CEFS203C	CEFS203D	CEFS203E	0.50mm
	CEFS204A	CEFS204B	CEFS204C	CEFS204D	CEFS204E	1.00mm
	CEFS205A	CEFS205B	CEFS205C	CEFS205D	CEFS205E	2.00mm

Box of 20 Units



Irrigation Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM206A	CEFM206B	CEFM206C	CEFM206D	CEFM206E	0.50mm
	CEFM207A	CEFM207B	CEFM207C	CEFM207D	CEFM207E	1.00mm
	CEFM208A	CEFM208B	CEFM208C	CEFM208D	CEFM208E	2.00mm
Single-Use	CEFS206A	CEFS206B	CEFS206C	CEFS206D	CEFS206E	0.50mm
	CEFS207A	CEFS207B	CEFS207C	CEFS207D	CEFS207E	1.00mm
	CEFS208A	CEFS208B	CEFS208C	CEFS208D	CEFS208E	2.00mm

Box of 20 Units



Irrigation Jensen Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM209A	CEFM209B	CEFM209C	CEFM209D	CEFM209E	0.50mm
	CEFM210A	CEFM210B	CEFM210C	CEFM210D	CEFM210E	1.00mm
	CEFM211A	CEFM211B	CEFM211C	CEFM211D	CEFM211E	2.00mm
Single-Use	CEFS209A	CEFS209B	CEFS209C	CEFS209D	CEFS209E	0.50mm
	CEFS210A	CEFS210B	CEFS210C	CEFS210D	CEFS210E	1.00mm
	CEFS211A	CEFS211B	CEFS211C	CEFS211D	CEFS211E	2.00mm

Box of 20 Units



Irrigation Mals Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM212A	CEFM212B	CEFM212C	CEFM212D	CEFM212E	0.50mm
	CEFM213A	CEFM213B	CEFM213C	CEFM213D	CEFM213E	1.00mm
	CEFM214A	CEFM214B	CEFM214C	CEFM214D	CEFM214E	2.00mm
Single-Use	CEFS212A	CEFS212B	CEFS212C	CEFS212D	CEFS212E	0.50mm
	CEFS213A	CEFS213B	CEFS213C	CEFS213D	CEFS213E	1.00mm
	CEFS214A	CEFS214B	CEFS214C	CEFS214D	CEFS214E	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length



Bayonet Premium Tips Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM215A	CEFM215B	CEFM215C	CEFM215D	CEFM215E	0.50mm
	CEFM217A	CEFM217B	CEFM217C	CEFM217D	CEFM217E	1.00mm
	CEFM218A	CEFM218B	CEFM218C	CEFM218D	CEFM218E	2.00mm
Single-Use	CEFS215A	CEFS215B	CEFS215C	CEFS215D	CEFS215E	0.50mm
	CEFS217A	CEFS217B	CEFS217C	CEFS217D	CEFS217E	1.00mm
	CEFS218A	CEFS218B	CEFS218C	CEFS218D	CEFS218E	2.00mm




Bayonet Monopolar Angle Up Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM219A	CEFM219B	CEFM219C	CEFM219D	CEFM219E	0.50mm
	CEFM220A	CEFM220B	CEFM220C	CEFM220D	CEFM220E	1.00mm
	CEFM221A	CEFM221B	CEFM221C	CEFM221D	CEFM221E	2.00mm
Single-Use	CEFS219A	CEFS219B	CEFS219C	CEFS219D	CEFS219E	0.50mm
	CEFS220A	CEFS220B	CEFS220C	CEFS220D	CEFS220E	1.00mm
	CEFS221A	CEFS221B	CEFS221C	CEFS221D	CEFS221E	2.00mm



Insulated Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM222A	CEFM222B	CEFM222C	CEFM222D	CEFM222E	0.50mm
	CEFM223A	CEFM223B	CEFM223C	CEFM223D	CEFM223E	1.00mm
	CEFM224A	CEFM224B	CEFM224C	CEFM224D	CEFM224E	2.00mm
Single-Use	CEFS222A	CEFS222B	CEFS222C	CEFS222D	CEFS222E	0.50mm
	CEFS223A	CEFS223B	CEFS223C	CEFS223D	CEFS223E	1.00mm
	CEFS224A	CEFS224B	CEFS224C	CEFS224D	CEFS224E	2.00mm



Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units
UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM225W	CEFM225A	CEFM225B	CEFM225C	0.50mm
	CEFM226W	CEFM226A	CEFM226B	CEFM226C	1.00mm
	CEFM227W	CEFM227A	CEFM227B	CEFM227C	2.00mm
Single-Use	CEFS225W	CEFS225A	CEFS225B	CEFS225C	0.50mm
	CEFS226W	CEFS226A	CEFS226B	CEFS226C	1.00mm
	CEFS227W	CEFS227A	CEFS227B	CEFS227C	2.00mm



Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM234A	CEFM234B	CEFM234C	CEFM234D	CEFM234E	0.50mm
	CEFM235A	CEFM235B	CEFM235C	CEFM235D	CEFM235E	1.00mm
	CEFM236A	CEFM236B	CEFM236C	CEFM236D	CEFM236E	2.00mm
Single-Use	CEFS234A	CEFS234B	CEFS234C	CEFS234D	CEFS234E	0.50mm
	CEFS235A	CEFS235B	CEFS235C	CEFS235D	CEFS235E	1.00mm
	CEFS236A	CEFS236B	CEFS236C	CEFS236D	CEFS236E	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM228A	CEFM228B	0.50mm
	CEFM229A	CEFM229B	1.00mm
	CEFM230A	CEFM230B	2.00mm
Single-Use	CEFS228A	CEFS228B	0.50mm
	CEFS229A	CEFS229B	1.00mm
	CEFS230A	CEFS230B	2.00mm




Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM237A	CEFM237B	CEFM237C	0.50mm
	CEFM238A	CEFM238B	CEFM238C	1.00mm
	CEFM239A	CEFM239B	CEFM239C	2.00mm
Single-Use	CEFS237A	CEFS237B	CEFS237C	0.50mm
	CEFS238A	CEFS238B	CEFS238C	1.00mm
	CEFS239A	CEFS239B	CEFS239C	2.00mm



Gilies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM231A	CEFM231B	CEFM231C	0.50mm
	CEFM232A	CEFM232B	CEFM232C	1.00mm
	CEFM233A	CEFM233B	CEFM233C	2.00mm
Single-Use	CEFS231A	CEFS231B	CEFS231C	0.50mm
	CEFS232A	CEFS232B	CEFS232C	1.00mm
	CEFS233A	CEFS233B	CEFS233C	2.00mm



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM240A	CEFM240B	CEFM240C	CEFM240D	CEFM240E	0.50mm
	CEFM241A	CEFM241B	CEFM241C	CEFM241D	CEFM241E	1.00mm
	CEFM242A	CEFM242B	CEFM242C	CEFM242D	CEFM242E	2.00mm
Single-Use	CEFS240A	CEFS240B	CEFS240C	CEFS240D	CEFS240E	0.50mm
	CEFS241A	CEFS241B	CEFS241C	CEFS241D	CEFS241E	1.00mm
	CEFS242A	CEFS242B	CEFS242C	CEFS242D	CEFS242E	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Monopolar Forceps

Banana Plug Connection



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM243A	CEFM243B	CEFM243C	CEFM243D	CEFM243E	0.50mm
	CEFM244A	CEFM244B	CEFM244C	CEFM244D	CEFM244E	1.00mm
	CEFM245A	CEFM245B	CEFM245C	CEFM245D	CEFM245E	2.00mm
Single-Use	CEFS243A	CEFS243B	CEFS243C	CEFS243D	CEFS243E	0.50mm
	CEFS244A	CEFS244B	CEFS244C	CEFS244D	CEFS244E	1.00mm
	CEFS245A	CEFS245B	CEFS245C	CEFS245D	CEFS245E	2.00mm

S

Serrated

ST

Serrated With 1/2 Toothed

SH

Serrated With 1/2 Toothed

PT

Plain With 1/2 Toothed

FOR ORDER
Please mention tip code along with product code



Regminton Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM246C	CEFM246D	CEFM246E	0.50mm
	CEFM247C	CEFM247D	CEFM247E	1.00mm
	CEFM248C	CEFM248D	CEFM248E	2.00mm
Single-Use	CEFS246C	CEFS246D	CEFS246E	0.50mm
	CEFS247C	CEFS247D	CEFS247E	1.00mm
	CEFS248C	CEFS248D	CEFS248E	2.00mm

S

Serrated

ST

Serrated With 1/2 Toothed

SH

Serrated With 1/2 Toothed

PT

Plain With 1/2 Toothed

FOR ORDER
Please mention tip code along with product code



Monopolar Grasping Holder 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM248A	CEFM248B	CEFM248C	CEFM248D	CEFM248E
	Angled CEFM249A	CEFM249B	CEFM249C	CEFM249D	CEFM249E
Single-Use	Straight CEFS248A	CEFS248B	CEFS248C	CEFS248D	CEFS248E
	Angled CEFS249A	CEFS249B	CEFS249C	CEFS249D	CEFS249E



Monopolar Dressing Scissors 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM250A	CEFM250B	CEFM250C	CEFM250D	CEFM250E
	Angled CEFM251A	CEFM251B	CEFM251C	CEFM251D	CEFM251E
Single-Use	Straight CEFS250A	CEFS250B	CEFS250C	CEFS250D	CEFS250E
	Angled CEFS251A	CEFS251B	CEFS251C	CEFS251D	CEFS251E



Monopolar Riche Forceps 48mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM252A	CEFM252B	CEFM252C	CEFM252D	CEFM252E
	Angled CEFM253A	CEFM253B	CEFM253C	CEFM253D	CEFM253E
Single-Use	Straight CEFS252A	CEFS252B	CEFS252C	CEFS252D	CEFS252E
	Angled CEFS253A	CEFS253B	CEFS253C	CEFS253D	CEFS253E

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units
UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Electrosurgical Forceps Serial


C E F M / S X X X X

Cosil Multi-Use/Single-Use Length



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM254W	CEFM254A	CEFM254B	CEFM254C	0.50mm
	CEFM255W	CEFM255A	CEFM255B	CEFM255C	1.00mm
	CEFM256W	CEFM256A	CEFM256B	CEFM256C	2.00mm
Single-Use	CEFS254W	CEFS254A	CEFS254B	CEFS254C	0.50mm
	CEFS255W	CEFS255A	CEFS255B	CEFS255C	1.00mm
	CEFS256W	CEFS256A	CEFS256B	CEFS256C	2.00mm



Dissecting Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM263A	CEFM263B	CEFM263C	CEFM263D	CEFM263E	0.50mm
	CEFM264A	CEFM264B	CEFM264C	CEFM264D	CEFM264E	1.00mm
	CEFM265A	CEFM265B	CEFM265C	CEFM265D	CEFM265E	2.00mm
Single-Use	CEFS263A	CEFS263B	CEFS263C	CEFS263D	CEFS263E	0.50mm
	CEFS264A	CEFS264B	CEFS264C	CEFS264D	CEFS264E	1.00mm
	CEFS265A	CEFS265B	CEFS265C	CEFS265D	CEFS265E	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM257A	CEFM257B	0.50mm
	CEFM258A	CEFM258B	1.00mm
	CEFM259A	CEFM259B	2.00mm
Single-Use	CEFS257A	CEFS257B	0.50mm
	CEFS258A	CEFS258B	1.00mm
	CEFS259A	CEFS259B	2.00mm




Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM266A	CEFM266B	CEFM266C	0.50mm
	CEFM267A	CEFM267B	CEFM267C	1.00mm
	CEFM268A	CEFM268B	CEFM268C	2.00mm
Single-Use	CEFS266A	CEFS266B	CEFS266C	0.50mm
	CEFS267A	CEFS267B	CEFS267C	1.00mm
	CEFS268A	CEFS268B	CEFS268C	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM260A	CEFM260B	CEFM260C	0.50mm
	CEFM261A	CEFM261B	CEFM261C	1.00mm
	CEFM262A	CEFM262B	CEFM262C	2.00mm
Single-Use	CEFS260A	CEFS260B	CEFS260C	0.50mm
	CEFS261A	CEFS261B	CEFS261C	1.00mm
	CEFS262A	CEFS262B	CEFS262C	2.00mm



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM269A	CEFM269B	CEFM269C	CEFM269D	CEFM269E	0.50mm
	CEFM270A	CEFM270B	CEFM270C	CEFM270D	CEFM270E	1.00mm
	CEFM271A	CEFM271B	CEFM271C	CEFM271D	CEFM271E	2.00mm
Single-Use	CEFS269A	CEFS269B	CEFS269C	CEFS269D	CEFS269E	0.50mm
	CEFS270A	CEFS270B	CEFS270C	CEFS270D	CEFS270E	1.00mm
	CEFS271A	CEFS271B	CEFS271C	CEFS271D	CEFS271E	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
CEFM/S XXXX
Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Monopolar Forceps

Without Connector, Female Plug Connection



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM272A	CEFM272B	CEFM272C	CEFM272D	CEFM272E	0.50mm
	CEFM273A	CEFM273B	CEFM273C	CEFM273D	CEFM273E	1.00mm
	CEFM274A	CEFM274B	CEFM274C	CEFM274D	CEFM274E	2.00mm
Single-Use	CEFS272A	CEFS272B	CEFS272C	CEFS272D	CEFS272E	0.50mm
	CEFS273A	CEFS273B	CEFS273C	CEFS273D	CEFS273E	1.00mm
	CEFS274A	CEFS274B	CEFS274C	CEFS274D	CEFS274E	2.00mm



Regminton Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM275C	CEFM275D	CEFM275E	0.50mm
	CEFM276C	CEFM276D	CEFM276E	1.00mm
	CEFM277C	CEFM277D	CEFM277E	2.00mm
Single-Use	CEFS275C	CEFS275D	CEFS275E	0.50mm
	CEFS276C	CEFS276D	CEFS276E	1.00mm
	CEFS277C	CEFS277D	CEFS277E	2.00mm



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM278W	CEFM278A	CEFM278B	CEFM278C	0.50mm
	CEFM279W	CEFM279A	CEFM279B	CEFM279C	1.00mm
	CEFM280W	CEFM280A	CEFM280B	CEFM280C	2.00mm
Single-Use	CEFS278W	CEFS278A	CEFS278B	CEFS278C	0.50mm
	CEFS279W	CEFS279A	CEFS279B	CEFS279C	1.00mm
	CEFS280W	CEFS280A	CEFS280B	CEFS280C	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM281A	CEFM281B	0.50mm
	CEFM282A	CEFM282B	1.00mm
	CEFM283A	CEFM283B	2.00mm
Single-Use	CEFS281A	CEFS281B	0.50mm
	CEFS282A	CEFS282B	1.00mm
	CEFS283A	CEFS283B	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM284A	CEFM284B	CEFM284C	0.50mm
	CEFM285A	CEFM285B	CEFM285C	1.00mm
	CEFM286A	CEFM286B	CEFM286C	2.00mm
Single-Use	CEFS284A	CEFS284B	CEFS284C	0.50mm
	CEFS285A	CEFS285B	CEFS285C	1.00mm
	CEFS286A	CEFS286B	CEFS286C	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

Box of 20 Units
UMAX: 0.5KVP

Also Available
 0.25mm
 0.75mm
 1.25mm
 2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
 Cosil Multi-Use/Single-Use Length



Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM287A	CEFM287B	CEFM287C	CEFM287D	CEFM287E	0.50mm
	CEFM288A	CEFM288B	CEFM288C	CEFM288D	CEFM288E	1.00mm
	CEFM289A	CEFM289B	CEFM289C	CEFM289D	CEFM289E	2.00mm
Single-Use	CEFS287A	CEFS287B	CEFS287C	CEFS287D	CEFS287E	0.50mm
	CEFS288A	CEFS288B	CEFS288C	CEFS288D	CEFS288E	1.00mm
	CEFS289A	CEFS289B	CEFS289C	CEFS289D	CEFS289E	2.00mm



FOR ORDER
Please mention tip code along with product code



Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM290A	CEFM290B	CEFM290C	0.50mm
	CEFM291A	CEFM291B	CEFM291C	1.00mm
	CEFM292A	CEFM292B	CEFM292C	2.00mm
Single-Use	CEFS290A	CEFS290B	CEFS290C	0.50mm
	CEFS291A	CEFS291B	CEFS291C	1.00mm
	CEFS292A	CEFS292B	CEFS292C	2.00mm




FOR ORDER
Please mention tip code along with product code



Dressing s Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM293A	CEFM293B	CEFM293C	CEFM293D	CEFM293E	0.50mm
	CEFM294A	CEFM294B	CEFM294C	CEFM294D	CEFM294E	1.00mm
	CEFM295A	CEFM295B	CEFM295C	CEFM295D	CEFM295E	2.00mm
Single-Use	CEFS293A	CEFS293B	CEFS293C	CEFS293D	CEFS293E	0.50mm
	CEFS294A	CEFS294B	CEFS294C	CEFS294D	CEFS294E	1.00mm
	CEFS295A	CEFS295B	CEFS295C	CEFS295D	CEFS295E	2.00mm



FOR ORDER
Please mention tip code along with product code



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM296A	CEFM296B	CEFM296C	CEFM296D	CEFM296E	0.50mm
	CEFM297A	CEFM297B	CEFM297C	CEFM297D	CEFM297E	1.00mm
	CEFM298A	CEFM298B	CEFM298C	CEFM298D	CEFM298E	2.00mm
Single-Use	CEFS296A	CEFS296B	CEFS296C	CEFS296D	CEFS296E	0.50mm
	CEFS297A	CEFS297B	CEFS297C	CEFS297D	CEFS297E	1.00mm
	CEFS298A	CEFS298B	CEFS298C	CEFS298D	CEFS298E	2.00mm



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM299A	CEFM299B	CEFM299C	CEFM299D	CEFM299E	0.50mm
	CEFM300A	CEFM300B	CEFM300C	CEFM300D	CEFM300E	1.00mm
	CEFM301A	CEFM301B	CEFM301C	CEFM301D	CEFM301E	2.00mm
Single-Use	CEFS299A	CEFS299B	CEFS299C	CEFS299D	CEFS299E	0.50mm
	CEFS300A	CEFS300B	CEFS300C	CEFS300D	CEFS300E	1.00mm
	CEFS301A	CEFS301B	CEFS301C	CEFS301D	CEFS301E	2.00mm



Regimintion Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM302C	CEFM302D	CEFM302E	0.50mm
	CEFM303C	CEFM303D	CEFM303E	1.00mm
	CEFM304C	CEFM304D	CEFM304E	2.00mm
Single-Use	CEFS302C	CEFS302D	CEFS302E	0.50mm
	CEFS303C	CEFS303D	CEFS303E	1.00mm
	CEFS304C	CEFS304D	CEFS304E	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
CEFM/S XXXX
Cosil Multi-Use/Single-Use Length



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM305W	CEFM305A	CEFM305B	CEFM305C	0.50mm
	CEFM306W	CEFM306A	CEFM306B	CEFM306C	1.00mm
	CEFM307W	CEFM307A	CEFM307B	CEFM307C	2.00mm
Single-Use	CEFS305W	CEFS305A	CEFS305B	CEFS305C	0.50mm
	CEFS306W	CEFS306A	CEFS306B	CEFS306C	1.00mm
	CEFS307W	CEFS307A	CEFS307B	CEFS307C	2.00mm



FOR ORDER Please mention tip code along with product code



Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM314A	CEFM314B	CEFM314C	CEFM314D	CEFM314E	0.50mm
	CEFM315A	CEFM315B	CEFM315C	CEFM315D	CEFM315E	1.00mm
	CEFM316A	CEFM316B	CEFM316C	CEFM316D	CEFM316E	2.00mm
Single-Use	CEFS314A	CEFS314B	CEFS314C	CEFS314D	CEFS314E	0.50mm
	CEFS315A	CEFS315B	CEFS315C	CEFS315D	CEFS315E	1.00mm
	CEFS316A	CEFS316B	CEFS316C	CEFS316D	CEFS316E	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM308A	CEFM308B	0.50mm
	CEFM309A	CEFM309B	1.00mm
	CEFM310A	CEFM310B	2.00mm
Single-Use	CEFS308A	CEFS308B	0.50mm
	CEFS309A	CEFS309B	1.00mm
	CEFS310A	CEFS310B	2.00mm



FOR ORDER Please mention tip code along with product code




Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM317A	CEFM317B	CEFM317C	0.50mm
	CEFM318A	CEFM318B	CEFM318C	1.00mm
	CEFM319A	CEFM319B	CEFM319C	2.00mm
Single-Use	CEFS317A	CEFS317B	CEFS317C	0.50mm
	CEFS318A	CEFS318B	CEFS318C	1.00mm
	CEFS319A	CEFS319B	CEFS319C	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM311A	CEFM311B	CEFM311C	0.50mm
	CEFM312A	CEFM312B	CEFM312C	1.00mm
	CEFM313A	CEFM313B	CEFM313C	2.00mm
Single-Use	CEFS311A	CEFS311B	CEFS311C	0.50mm
	CEFS312A	CEFS312B	CEFS312C	1.00mm
	CEFS313A	CEFS313B	CEFS313C	2.00mm



FOR ORDER Please mention tip code along with product code



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM320A	CEFM320B	CEFM320C	CEFM320D	CEFM320E	0.50mm
	CEFM321A	CEFM321B	CEFM321C	CEFM321D	CEFM321E	1.00mm
	CEFM322A	CEFM322B	CEFM322C	CEFM322D	CEFM322E	2.00mm
Single-Use	CEFS320A	CEFS320B	CEFS320C	CEFS320D	CEFS320E	0.50mm
	CEFS321A	CEFS321B	CEFS321C	CEFS321D	CEFS321E	1.00mm
	CEFS322A	CEFS322B	CEFS322C	CEFS322D	CEFS322E	2.00mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 20 Units UMAX: 0.5KVP

Also Available
Tip Sizes
0.25mm
0.75mm
1.25mm
2.50mm

Electrosurgical Forceps Serial
C E F M / S XXX X
Cosil Multi-Use/Single-Use Length

Electrosurgical Insulated Monopolar Forceps

4.8 Plug Connection



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM323A	CEFM323B	CEFM323C	CEFM323D	CEFM323E	0.50mm
	CEFM324A	CEFM324B	CEFM324C	CEFM324D	CEFM324E	1.00mm
	CEFM325A	CEFM325B	CEFM325C	CEFM325D	CEFM325E	2.00mm
Single-Use	CEFS323A	CEFS323B	CEFS323C	CEFS323D	CEFS323E	0.50mm
	CEFS324A	CEFS324B	CEFS324C	CEFS324D	CEFS324E	1.00mm
	CEFS325A	CEFS325B	CEFS325C	CEFS325D	CEFS325E	2.00mm



FOR ORDER
Please mention tip code along with product code



Regminton Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM326C	CEFM326D	CEFM326E	0.50mm
	CEFM327C	CEFM327D	CEFM327E	1.00mm
	CEFM328C	CEFM328D	CEFM328E	2.00mm
Single-Use	CEFS326C	CEFS326D	CEFS326E	0.50mm
	CEFS327C	CEFS327D	CEFS327E	1.00mm
	CEFS328C	CEFS328D	CEFS328E	2.00mm



FOR ORDER
Please mention tip code along with product code



Monopolar Grasping Holder 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM328A	CEFM328B	CEFM328C	CEFM328D	CEFM328E
	Angled CEFM329A	CEFM329B	CEFM329C	CEFM329D	CEFM329E
Single-Use	Straight CEFS328A	CEFS328B	CEFS328C	CEFS328D	CEFS328E
	Angled CEFS329A	CEFS329B	CEFS329C	CEFS329D	CEFS329E



Monopolar Dressing Scissors 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM330A	CEFM330B	CEFM330C	CEFM330D	CEFM330E
	Angled CEFM331A	CEFM331B	CEFM331C	CEFM331D	CEFM331E
Single-Use	Straight CEFS330A	CEFS330B	CEFS330C	CEFS330D	CEFS330E
	Angled CEFS331A	CEFS331B	CEFS331C	CEFS331D	CEFS331E



Monopolar Riche Forceps 48mm


	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM332A	CEFM332B	CEFM332C	CEFM332D	CEFM332E
	Angled CEFM333A	CEFM333B	CEFM333C	CEFM333D	CEFM333E
Single-Use	Straight CEFS332A	CEFS332B	CEFS332C	CEFS332D	CEFS332E
	Angled CEFS333A	CEFS333B	CEFS333C	CEFS333D	CEFS333E

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp Box of 20 Units
UMAX: 0.5KVP

Also Available
Tip Sizes: 0.25mm, 0.75mm, 1.25mm, 2.50mm

Electrosurgical Forceps Serial
CEFM/S XXXX
Cosil Multi-Use/Single-Use Length



Bipolar Vessel Sealer Stainless Steel

Curve

Angle



BOX OF 20 UNITS



BOX OF 20 UNITS



BOX OF 20 UNITS

	165mm	190mm	215mm	240mm
Multi-Use	Curve CPFM101K	CPFM101L	CPFM101M	CPFM101N
	Angle CPFM102K	CPFM102K	CPFM102M	CPFM102N
Single-Use	Curve CPFS101K	CPFS101L	CPFS101M	CPFS101N
	Angle CPFS102K	CPFS102L	CPFS102M	CPFS102N

	165mm	190mm	215mm	240mm
Multi-Use	Curve CPFM103K	CPFM103L	CPFM103M	CPFM103N
	Angle CPFM104K	CPFM104L	CPFM104M	CPFM104N
Single-Use	Curve CPFS103K	CPFS103L	CPFS103M	CPFS103N
	Angle CPFS104K	CPFS104L	CPFS104M	CPFS104N

	165mm	190mm	215mm	240mm
Multi-Use	Curve CPFM105K	CPFM105L	CPFM105M	CPFM105N
	Angle CPFM106K	CPFM106L	CPFM106M	CPFM106N
Single-Use	Curve CPFS105K	CPFS105L	CPFS105M	CPFS105N
	Angle CPFS106K	CPFS106L	CPFS106M	CPFS106N

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Premium Forceps Serial
C P F M / S XXX X
 Cosil Multi-Use/Single-Use Length



	165mm	190mm	215mm	240mm
Curve	CPFM107K	CPFM107L	CPFM107M	CPFM107N
Angle	CPFM108K	CPFM108L	CPFM108M	CPFM108N
Curve	CPFS107K	CPFS107L	CPFS107M	CPFS107N
Angle	CPFS108K	CPFS108L	CPFS108M	CPFS108N





	165mm	190mm	215mm	240mm
Curve	CPFM109K	CPFM109L	CPFM109M	CPFM109N
Angle	CPFM110K	CPFM110L	CPFM110M	CPFM110N
Curve	CPFS109K	CPFS109L	CPFS109M	CPFS109N
Angle	CPFS110K	CPFS110L	CPFS110M	CPFS110N



	165mm	190mm	215mm	240mm
Curve	CPFM111K	CPFM111L	CPFM111M	CPFM111N
Angle	CPFM112K	CPFM112L	CPFM112M	CPFM112N
Curve	CPFS111K	CPFS111L	CPFS111M	CPFS111N
Angle	CPFS112K	CPFS112L	CPFS112M	CPFS112N

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand





Bipolar Cut-Coag Scissor -Ceramics



BOX OF 20 UNITS



BOX OF 20 UNITS



BOX OF 20 UNITS

	180mm	200mm	220mm
Multi-Use	CPFM113P	CPFM113Q	CPFM113R
Single-Use	CPFS113P	CPFS113Q	CPFS113R

	180mm	200mm	220mm
Multi-Use	CPFM114P	CPFM114Q	CPFM114R
Single-Use	CPFS114P	CPFS114Q	CPFS114R

	180mm	200mm	220mm
Multi-Use	CPFM115P	CPFM115Q	CPFM115R
Single-Use	CPFS115P	CPFS115Q	CPFS115R

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5KVP

Premium Forceps Serial
C P F M / S XXX X
Cosil Multi-Use/Single-Use Length



	102mm	115mm	
Multi-Use	CPFM116T	CPFM116V	0.50mm
	CPFM117T	CPFM117V	1.00mm
	CPFM118T	CPFM118V	2.00mm
Single-Use	CPFS116T	CPFS116V	0.50mm
	CPFS117T	CPFS117V	1.00mm
	CPFS118T	CPFS118V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM119T	CPFM119V	0.50mm
	CPFM120T	CPFM120V	1.00mm
	CPFM121T	CPFM121V	2.00mm
Single-Use	CPFS119T	CPFS119V	0.50mm
	CPFS120T	CPFS120V	1.00mm
	CPFS121T	CPFS121V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM122X	CPFM122A	CPFM122B	CPFM122C	CPFM122D	CPFM122P	0.50mm
	CPFM123X	CPFM123A	CPFM123B	CPFM123C	CPFM123D	CPFM123P	1.00mm
	CPFM124X	CPFM124A	CPFM124B	CPFM124C	CPFM124D	CPFM124P	2.00mm
Single-Use	CPFS122X	CPFS122A	CPFS122B	CPFS122C	CPFS122D	CPFS122P	0.50mm
	CPFS123X	CPFS123A	CPFS123B	CPFS123C	CPFS123D	CPFS123P	1.00mm
	CPFS124X	CPFS124A	CPFS124B	CPFS124C	CPFS124D	CPFS124P	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM125T	CPFM125V	0.50mm
	CPFM126T	CPFM126V	1.00mm
	CPFM127T	CPFM127V	2.00mm
Single-Use	CPFS125T	CPFS125V	0.50mm
	CPFS126T	CPFS126V	1.00mm
	CPFS127T	CPFS127V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM128T	CPFM128V	0.50mm
	CPFM129T	CPFM129V	1.00mm
	CPFM130T	CPFM130V	2.00mm
Single-Use	CPFS128T	CPFS128V	0.50mm
	CPFS129T	CPFS129V	1.00mm
	CPFS130T	CPFS130V	2.00mm



Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM131X	CPFM131A	CPFM131B	CPFM131C	CPFM131D	CPFM131P	0.50mm
	CPFM132X	CPFM132A	CPFM132B	CPFM132C	CPFM132D	CPFM132P	1.00mm
	CPFM133X	CPFM133A	CPFM133B	CPFM133C	CPFM133D	CPFM133P	2.00mm
Single-Use	CPFS131X	CPFS131A	CPFS131B	CPFS131C	CPFS131D	CPFS131P	0.50mm
	CPFS132X	CPFS132A	CPFS132B	CPFS132C	CPFS132D	CPFS132P	1.00mm
	CPFS133X	CPFS133A	CPFS133B	CPFS133C	CPFS133D	CPFS133P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand





	102mm	115mm	
Multi-Use	CPFM134T	CPFM134V	0.50mm
	CPFM135T	CPFM135V	1.00mm
	CPFM136T	CPFM136V	2.00mm
Single-Use	CPFS134T	CPFS134V	0.50mm
	CPFS135T	CPFS135V	1.00mm
	CPFS136T	CPFS136V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM143T	CPFM143V	0.50mm
	CPFM144T	CPFM144V	1.00mm
	CPFM145T	CPFM145V	2.00mm
Single-Use	CPFS143T	CPFS143V	0.50mm
	CPFS144T	CPFS144V	1.00mm
	CPFS145T	CPFS145V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM137T	CPFM137V	0.50mm
	CPFM138T	CPFM138V	1.00mm
	CPFM139T	CPFM139V	2.00mm
Single-Use	CPFS137T	CPFS137V	0.50mm
	CPFS138T	CPFS138V	1.00mm
	CPFS139T	CPFS139V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM146T	CPFM146V	0.50mm
	CPFM147T	CPFM147V	1.00mm
	CPFM148T	CPFM148V	2.00mm
Single-Use	CPFS146T	CPFS146V	0.50mm
	CPFS147T	CPFS147V	1.00mm
	CPFS148T	CPFS148V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM140X	CPFM140A	CPFM140B	CPFM140C	CPFM140D	CPFM140P	0.50mm
	CPFM141X	CPFM141A	CPFM141B	CPFM141C	CPFM141D	CPFM141P	1.00mm
	CPFM142X	CPFM142A	CPFM142B	CPFM142C	CPFM142D	CPFM142P	2.00mm
Single-Use	CPFS140X	CPFS140A	CPFS140B	CPFS140C	CPFS140D	CPFS140P	0.50mm
	CPFS141X	CPFS141A	CPFS141B	CPFS141C	CPFS141D	CPFS141P	1.00mm
	CPFS142X	CPFS142A	CPFS142B	CPFS142C	CPFS142D	CPFS142P	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM149X	CPFM149A	CPFM149B	CPFM149C	CPFM149D	CPFM149P	0.50mm
	CPFM150X	CPFM150A	CPFM150B	CPFM150C	CPFM150D	CPFM150P	1.00mm
	CPFM151X	CPFM151A	CPFM151B	CPFM151C	CPFM151D	CPFM151P	2.00mm
Single-Use	CPFS149X	CPFS149A	CPFS149B	CPFS149C	CPFS149D	CPFS149P	0.50mm
	CPFS150X	CPFS150A	CPFS150B	CPFS150C	CPFS150D	CPFS150P	1.00mm
	CPFS151X	CPFS151A	CPFS151B	CPFS151C	CPFS151D	CPFS151P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

max Vp UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial

C P F M / S X X X X

Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM152A	CPFM152B	CPFM152C	CPFM152D	CPFM152P	0.50mm
	CPFM153A	CPFM153B	CPFM153C	CPFM153D	CPFM153P	1.00mm
	CPFM154A	CPFM154B	CPFM154C	CPFM154D	CPFM154P	2.00mm
Single-Use	CPFS152A	CPFS152B	CPFS152C	CPFS152D	CPFS152P	0.50mm
	CPFS153A	CPFS153B	CPFS153C	CPFS153D	CPFS153P	1.00mm
	CPFS154A	CPFS154B	CPFS154C	CPFS154D	CPFS154P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM155A	CPFM155B	CPFM155C	CPFM155D	CPFM155P	0.50mm
	CPFM156A	CPFM156B	CPFM156C	CPFM156D	CPFM156P	1.00mm
	CPFM157A	CPFM157B	CPFM157C	CPFM157D	CPFM157P	2.00mm
Single-Use	CPFS155A	CPFS155B	CPFS155C	CPFS155D	CPFS155P	0.50mm
	CPFS156A	CPFS156B	CPFS156C	CPFS156D	CPFS156P	1.00mm
	CPFS157A	CPFS157B	CPFS157C	CPFS157D	CPFS157P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM158A	CPFM158B	CPFM158C	CPFM158D	CPFM158P	0.50mm
	CPFM159A	CPFM159B	CPFM159C	CPFM159D	CPFM159P	1.00mm
	CPFM160A	CPFM160B	CPFM160C	CPFM160D	CPFM160P	2.00mm
Single-Use	CPFS158A	CPFS158B	CPFS158C	CPFS158D	CPFS158P	0.50mm
	CPFS159A	CPFS159B	CPFS159C	CPFS159D	CPFS159P	1.00mm
	CPFS160A	CPFS160B	CPFS160C	CPFS160D	CPFS160P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM161A	CPFM161B	CPFM161C	CPFM161D	CPFM161P	0.50mm
	CPFM162A	CPFM162B	CPFM162C	CPFM162D	CPFM162P	1.00mm
	CPFM163A	CPFM163B	CPFM163C	CPFM163D	CPFM163P	2.00mm
Single-Use	CPFS161A	CPFS161B	CPFS161C	CPFS161D	CPFS161P	0.50mm
	CPFS162A	CPFS162B	CPFS162C	CPFS162D	CPFS162P	1.00mm
	CPFS163A	CPFS163B	CPFS163C	CPFS163D	CPFS163P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM164A	CPFM164B	CPFM164C	CPFM164D	CPFM164P	0.50mm
	CPFM165A	CPFM165B	CPFM165C	CPFM165D	CPFM165P	1.00mm
	CPFM166A	CPFM166B	CPFM166C	CPFM166D	CPFM166P	2.00mm
Single-Use	CPFS164A	CPFS164B	CPFS164C	CPFS164D	CPFS164P	0.50mm
	CPFS165A	CPFS165B	CPFS165C	CPFS165D	CPFS165P	1.00mm
	CPFS166A	CPFS166B	CPFS166C	CPFS166D	CPFS166P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

max UMAX: 0.5KVP

Also Available
 0.25mm
 0.75mm
 1.25mm
 2.50mm

Premium Forceps Serial
C P F M / S X X X X
 Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM167A	CPFM167B	CPFM167C	CPFM167D	CPFM167P	0.50mm
	CPFM168A	CPFM168B	CPFM168C	CPFM168D	CPFM168P	1.00mm
	CPFM169A	CPFM169B	CPFM169C	CPFM169D	CPFM169P	2.00mm
Single-Use	CPFS167A	CPFS167B	CPFS167C	CPFS167D	CPFS167P	0.50mm
	CPFS168A	CPFS168B	CPFS168C	CPFS168D	CPFS168P	1.00mm
	CPFS169A	CPFS169B	CPFS169C	CPFS169D	CPFS169P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM170A	CPFM170B	CPFM170C	CPFM170D	CPFM170P	0.50mm
	CPFM171A	CPFM171B	CPFM171C	CPFM171D	CPFM171P	1.00mm
	CPFM172A	CPFM172B	CPFM172C	CPFM172D	CPFM172P	2.00mm
Single-Use	CPFS170A	CPFS170B	CPFS170C	CPFS170D	CPFS170P	0.50mm
	CPFS171A	CPFS171B	CPFS171C	CPFS171D	CPFS171P	1.00mm
	CPFS172A	CPFS172B	CPFS172C	CPFS172D	CPFS172P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM173A	CPFM173B	CPFM173C	CPFM173D	CPFM173P	0.50mm
	CPFM174A	CPFM174B	CPFM174C	CPFM174D	CPFM174P	1.00mm
	CPFM175A	CPFM175B	CPFM175C	CPFM175D	CPFM175P	2.00mm
Single-Use	CPFS173A	CPFS173B	CPFS173C	CPFS173D	CPFS173P	0.50mm
	CPFS174A	CPFS174B	CPFS174C	CPFS174D	CPFS174P	1.00mm
	CPFS175A	CPFS175B	CPFS175C	CPFS175D	CPFS175P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM176A	CPFM176B	CPFM176C	CPFM176D	CPFM176P	0.50mm
	CPFM177A	CPFM177B	CPFM177C	CPFM177D	CPFM177P	1.00mm
	CPFM178A	CPFM178B	CPFM178C	CPFM178D	CPFM178P	2.00mm
Single-Use	CPFS176A	CPFS176B	CPFS176C	CPFS176D	CPFS176P	0.50mm
	CPFS177A	CPFS177B	CPFS177C	CPFS177D	CPFS177P	1.00mm
	CPFS178A	CPFS178B	CPFS178C	CPFS178D	CPFS178P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM179A	CPFM179B	CPFM179C	CPFM179D	CPFM179P	0.50mm
	CPFM180A	CPFM180B	CPFM180C	CPFM180D	CPFM180P	1.00mm
	CPFM181A	CPFM181B	CPFM181C	CPFM181D	CPFM181P	2.00mm
Single-Use	CPFS179A	CPFS179B	CPFS179C	CPFS179D	CPFS179P	0.50mm
	CPFS180A	CPFS180B	CPFS180C	CPFS180D	CPFS180P	1.00mm
	CPFS181A	CPFS181B	CPFS181C	CPFS181D	CPFS181P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available


- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial
C P F M / S XXX X
Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM182A	CPFM182B	CPFM182C	CPFM182D	CPFM182P	0.50mm
	CPFM183A	CPFM183B	CPFM183C	CPFM183D	CPFM183P	1.00mm
	CPFM184A	CPFM184B	CPFM184C	CPFM184D	CPFM184P	2.00mm
Single-Use	CPFS182A	CPFS182B	CPFS182C	CPFS182D	CPFS182P	0.50mm
	CPFS183A	CPFS183B	CPFS183C	CPFS183D	CPFS183P	1.00mm
	CPFS184A	CPFS184B	CPFS184C	CPFS184D	CPFS184P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM188A	CPFM188B	CPFM188C	CPFM188D	CPFM188P	0.50mm
	CPFM189A	CPFM189B	CPFM189C	CPFM189D	CPFM189P	1.00mm
	CPFM190A	CPFM190B	CPFM190C	CPFM190D	CPFM190P	2.00mm
Single-Use	CPFS188A	CPFS188B	CPFS188C	CPFS188D	CPFS188P	0.50mm
	CPFS189A	CPFS189B	CPFS189C	CPFS189D	CPFS189P	1.00mm
	CPFS190A	CPFS190B	CPFS190C	CPFS190D	CPFS190P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM185A	CPFM185B	CPFM185C	CPFM185D	CPFM185P	0.50mm
	CPFM186A	CPFM186B	CPFM186C	CPFM186D	CPFM186P	1.00mm
	CPFM187A	CPFM187B	CPFM187C	CPFM187D	CPFM187P	2.00mm
Single-Use	CPFS185A	CPFS185B	CPFS185C	CPFS185D	CPFS185P	0.50mm
	CPFS186A	CPFS186B	CPFS186C	CPFS186D	CPFS186P	1.00mm
	CPFS187A	CPFS187B	CPFS187C	CPFS187D	CPFS187P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM191A	CPFM191B	CPFM191C	CPFM191D	CPFM191P	0.50mm
	CPFM192A	CPFM192B	CPFM192C	CPFM192D	CPFM192P	1.00mm
	CPFM193A	CPFM193B	CPFM193C	CPFM193D	CPFM193P	2.00mm
Single-Use	CPFS191A	CPFS191B	CPFS191C	CPFS191D	CPFS191P	0.50mm
	CPFS192A	CPFS192B	CPFS192C	CPFS192D	CPFS192P	1.00mm
	CPFS193A	CPFS193B	CPFS193C	CPFS193D	CPFS193P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial
C P F M / S X X X X
 Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM194A	CPFM194B	CPFM194C	CPFM194D	CPFM194P	0.50mm
	CPFM195A	CPFM195B	CPFM195C	CPFM195D	CPFM195P	1.00mm
	CPFM196A	CPFM196B	CPFM196C	CPFM196D	CPFM196P	2.00mm
Single-Use	CPFS194A	CPFS194B	CPFS194C	CPFS194D	CPFS194P	0.50mm
	CPFS195A	CPFS195B	CPFS195C	CPFS195D	CPFS195P	1.00mm
	CPFS196A	CPFS196B	CPFS196C	CPFS196D	CPFS196P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM203A	CPFM203B	CPFM203C	CPFM203D	CPFM203P	0.50mm
	CPFM204A	CPFM204B	CPFM204C	CPFM204D	CPFM204P	1.00mm
	CPFM205A	CPFM205B	CPFM205C	CPFM205D	CPFM205P	2.00mm
Single-Use	CPFS203A	CPFS203B	CPFS203C	CPFS203D	CPFS203P	0.50mm
	CPFS204A	CPFS204B	CPFS204C	CPFS204D	CPFS204P	1.00mm
	CPFS205A	CPFS205B	CPFS205C	CPFS205D	CPFS205P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM197A	CPFM197B	CPFM197C	CPFM197D	CPFM197P	0.50mm
	CPFM198A	CPFM198B	CPFM198C	CPFM198D	CPFM198P	1.00mm
	CPFM199A	CPFM199B	CPFM199C	CPFM199D	CPFM199P	2.00mm
Single-Use	CPFS197A	CPFS197B	CPFS197C	CPFS197D	CPFS197P	0.50mm
	CPFS198A	CPFS198B	CPFS198C	CPFS198D	CPFS198P	1.00mm
	CPFS199A	CPFS199B	CPFS199C	CPFS199D	CPFS199P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM206A	CPFM206B	CPFM206C	CPFM206D	CPFM206P	0.50mm
	CPFM207A	CPFM207B	CPFM207C	CPFM207D	CPFM207P	1.00mm
	CPFM208A	CPFM208B	CPFM208C	CPFM208D	CPFM208P	2.00mm
Single-Use	CPFS206A	CPFS206B	CPFS206C	CPFS206D	CPFS206P	0.50mm
	CPFS207A	CPFS207B	CPFS207C	CPFS207D	CPFS207P	1.00mm
	CPFS208A	CPFS208B	CPFS208C	CPFS208D	CPFS208P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM200A	CPFM200B	CPFM200C	CPFM200D	CPFM200P	0.50mm
	CPFM201A	CPFM201B	CPFM201C	CPFM201D	CPFM201P	1.00mm
	CPFM202A	CPFM202B	CPFM202C	CPFM202D	CPFM202P	2.00mm
Single-Use	CPFS200A	CPFS200B	CPFS200C	CPFS200D	CPFS200P	0.50mm
	CPFS201A	CPFS201B	CPFS201C	CPFS201D	CPFS201P	1.00mm
	CPFS202A	CPFS202B	CPFS202C	CPFS202D	CPFS202P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM209A	CPFM209B	CPFM209C	CPFM209D	CPFM209P	0.50mm
	CPFM210A	CPFM210B	CPFM210C	CPFM210D	CPFM210P	1.00mm
	CPFM211A	CPFM211B	CPFM211C	CPFM211D	CPFM211P	2.00mm
Single-Use	CPFS209A	CPFS209B	CPFS209C	CPFS209D	CPFS209P	0.50mm
	CPFS210A	CPFS210B	CPFS210C	CPFS210D	CPFS210P	1.00mm
	CPFS211A	CPFS211B	CPFS211C	CPFS211D	CPFS211P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

max UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial
C P F M / S X X X X
 Cosil Multi-Use/Single-Use Length



	102mm	115mm	
Multi-Use	CPFM212T	CPFM212V	0.50mm
	CPFM213T	CPFM213V	1.00mm
	CPFM214T	CPFM214V	2.00mm
Single-Use	CPFS212T	CPFS212V	0.50mm
	CPFS213T	CPFS213V	1.00mm
	CPFS214T	CPFS214V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM218X	CPFM218A	CPFM218B	CPFM218C	CPFM218D	CPFM218P	CPFM218F	0.50mm
	CPFM219X	CPFM219A	CPFM219B	CPFM219C	CPFM219D	CPFM219P	CPFM219F	1.00mm
	CPFM220X	CPFM220A	CPFM220B	CPFM220C	CPFM220D	CPFM220P	CPFM220F	2.00mm
Single-Use	CPFS218X	CPFS218A	CPFS218B	CPFS218C	CPFS218D	CPFS218P	CPFS218F	0.50mm
	CPFS219X	CPFS219A	CPFS219B	CPFS219C	CPFS219D	CPFS219P	CPFS219F	1.00mm
	CPFS220X	CPFS220A	CPFS220B	CPFS220C	CPFS220D	CPFS220P	CPFS220F	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM215T	CPFM215V	0.50mm
	CPFM216T	CPFM216V	1.00mm
	CPFM217T	CPFM217V	2.00mm
Single-Use	CPFS215T	CPFS215V	0.50mm
	CPFS216T	CPFS216V	1.00mm
	CPFS217T	CPFS217V	2.00mm


Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM221X	CPFM221A	CPFM221B	CPFM221C	CPFM221D	CPFM221P	CPFM221F	0.50mm
	CPFM222X	CPFM222A	CPFM222B	CPFM222C	CPFM222D	CPFM222P	CPFM222F	1.00mm
	CPFM223X	CPFM223A	CPFM223B	CPFM223C	CPFM223D	CPFM223P	CPFM223F	2.00mm
Single-Use	CPFS221X	CPFS221A	CPFS221B	CPFS221C	CPFS221D	CPFS221P	CPFS221F	0.50mm
	CPFS222X	CPFS222A	CPFS222B	CPFS222C	CPFS222D	CPFS222P	CPFS222F	1.00mm
	CPFS223X	CPFS223A	CPFS223B	CPFS223C	CPFS223D	CPFS223P	CPFS223F	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Premium Forceps Serial
C P F M / S XXX X
 Cosil Multi-Use/Single-Use Length



	102mm	115mm	
Multi-Use	CPFM224T	CPFM224V	0.50mm
	CPFM225T	CPFM225V	1.00mm
	CPFM226T	CPFM226V	2.00mm
Single-Use	CPFS224T	CPFS224V	0.50mm
	CPFS225T	CPFS225V	1.00mm
	CPFS226T	CPFS226V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM230X	CPFM230A	CPFM230B	CPFM230C	CPFM230D	CPFM230P	CPFM230F	0.50mm
	CPFM231X	CPFM231A	CPFM231B	CPFM231C	CPFM231D	CPFM231P	CPFM231F	1.00mm
	CPFM232X	CPFM232A	CPFM232B	CPFM232C	CPFM232D	CPFM232P	CPFM232F	2.00mm
Single-Use	CPFS230X	CPFS230A	CPFS230B	CPFS230C	CPFS230D	CPFS230P	CPFS230F	0.50mm
	CPFS231X	CPFS231A	CPFS231B	CPFS231C	CPFS231D	CPFS231P	CPFS231F	1.00mm
	CPFS232X	CPFS232A	CPFS232B	CPFS232C	CPFS232D	CPFS232P	CPFS232F	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM227T	CPFM227V	0.50mm
	CPFM228T	CPFM228V	1.00mm
	CPFM229T	CPFM229V	2.00mm
Single-Use	CPFS227T	CPFS227V	0.50mm
	CPFS228T	CPFS228V	1.00mm
	CPFS229T	CPFS229V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM233X	CPFM233A	CPFM233B	CPFM233C	CPFM233D	CPFM233P	CPFM233F	0.50mm
	CPFM234X	CPFM234A	CPFM234B	CPFM234C	CPFM234D	CPFM234P	CPFM234F	1.00mm
	CPFM235X	CPFM235A	CPFM235B	CPFM235C	CPFM235D	CPFM235P	CPFM235F	2.00mm
Single-Use	CPFS233X	CPFS233A	CPFS233B	CPFS233C	CPFS233D	CPFS233P	CPFS233F	0.50mm
	CPFS234X	CPFS234A	CPFS234B	CPFS234C	CPFS234D	CPFS234P	CPFS234F	1.00mm
	CPFS235X	CPFS235A	CPFS235B	CPFS235C	CPFS235D	CPFS235P	CPFS235F	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5KVP

Also Available
0.25mm
0.75mm
1.25mm
2.50mm

Premium Forceps Serial
C P F M / S XXX X
Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM236A	CPFM236B	CPFM236C	CPFM236D	CPFM236P	0.50mm
	CPFM237A	CPFM237B	CPFM237C	CPFM237D	CPFM237P	1.00mm
	CPFM238A	CPFM238B	CPFM238C	CPFM238D	CPFM238P	2.00mm
Single-Use	CPFS236A	CPFS236B	CPFS236C	CPFS236D	CPFS236P	0.50mm
	CPFS237A	CPFS237B	CPFS237C	CPFS237D	CPFS237P	1.00mm
	CPFS238A	CPFS238B	CPFS238C	CPFS238D	CPFS238P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM239A	CPFM239B	CPFM239C	CPFM239D	CPFM239P	0.50mm
	CPFM240A	CPFM240B	CPFM240C	CPFM240D	CPFM240P	1.00mm
	CPFM241A	CPFM241B	CPFM241C	CPFM241D	CPFM241P	2.00mm
Single-Use	CPFS239A	CPFS239B	CPFS239C	CPFS239D	CPFS239P	0.50mm
	CPFS240A	CPFS240B	CPFS240C	CPFS240D	CPFS240P	1.00mm
	CPFS241A	CPFS241B	CPFS241C	CPFS241D	CPFS241P	2.00mm

Box of 20 Units




	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM242A	CPFM242B	CPFM242C	CPFM242D	CPFM242P	0.50mm
	CPFM243A	CPFM243B	CPFM243C	CPFM243D	CPFM243P	1.00mm
	CPFM245A	CPFM245B	CPFM245C	CPFM245D	CPFM245P	2.00mm
Single-Use	CPFS242A	CPFS242B	CPFS242C	CPFS242D	CPFS242P	0.50mm
	CPFS243A	CPFS243B	CPFS243C	CPFS243D	CPFS243P	1.00mm
	CPFS245A	CPFS245B	CPFS245C	CPFS245D	CPFS245P	2.00mm

Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5KVP



Premium Forceps Serial
C P F M / S XXX X
 Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM246A	CPFM246B	CPFM246C	CPFM246D	CPFM246P	0.50mm
	CPFM247A	CPFM247B	CPFM247C	CPFM247D	CPFM247P	1.00mm
	CPFM248A	CPFM248B	CPFM248C	CPFM248D	CPFM248P	2.00mm
Single-Use	CPFS246A	CPFS246B	CPFS246C	CPFS246D	CPFS246P	0.50mm
	CPFS247A	CPFS247B	CPFS247C	CPFS247D	CPFS247P	1.00mm
	CPFS248A	CPFS248B	CPFS248C	CPFS248D	CPFS248P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM249A	CPFM249B	CPFM249C	CPFM249D	CPFM249P	0.50mm
	CPFM250A	CPFM250B	CPFM250C	CPFM250D	CPFM250P	1.00mm
	CPFM251A	CPFM251B	CPFM251C	CPFM251D	CPFM251P	2.00mm
Single-Use	CPFS249A	CPFS249B	CPFS249C	CPFS249D	CPFS249P	0.50mm
	CPFS250A	CPFS250B	CPFS250C	CPFS250D	CPFS250P	1.00mm
	CPFS251A	CPFS251B	CPFS251C	CPFS251D	CPFS251P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM252A	CPFM252B	CPFM252C	CPFM252D	CPFM252P	0.50mm
	CPFM253A	CPFM253B	CPFM253C	CPFM253D	CPFM253P	1.00mm
	CPFM254A	CPFM254B	CPFM254C	CPFM254D	CPFM254P	2.00mm
Single-Use	CPFS252A	CPFS252B	CPFS252C	CPFS252D	CPFS252P	0.50mm
	CPFS253A	CPFS253B	CPFS253C	CPFS253D	CPFS253P	1.00mm
	CPFS254A	CPFS254B	CPFS254C	CPFS254D	CPFS254P	2.00mm

Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand

UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial
C P F M / S XXX X
 Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM255A	CPFM255B	CPFM255C	CPFM255D	CPFM255P	0.50mm
	CPFM256A	CPFM256B	CPFM256C	CPFM256D	CPFM256P	1.00mm
	CPFM257A	CPFM257B	CPFM257C	CPFM257D	CPFM257P	2.00mm
Single-Use	CPFS255A	CPFS255B	CPFS255C	CPFS255D	CPFS255P	0.50mm
	CPFS256A	CPFS256B	CPFS256C	CPFS256D	CPFS256P	1.00mm
	CPFS257A	CPFS257B	CPFS257C	CPFS257D	CPFS257P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM261A	CPFM261B	CPFM261C	CPFM261D	CPFM261P	0.50mm
	CPFM262A	CPFM262B	CPFM262C	CPFM262D	CPFM262P	1.00mm
	CPFM263A	CPFM263B	CPFM263C	CPFM263D	CPFM263P	2.00mm
Single-Use	CPFS261A	CPFS261B	CPFS261C	CPFS261D	CPFS261P	0.50mm
	CPFS262A	CPFS262B	CPFS262C	CPFS262D	CPFS262P	1.00mm
	CPFS263A	CPFS263B	CPFS263C	CPFS263D	CPFS263P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM258A	CPFM258B	CPFM258C	CPFM258D	CPFM258P	0.50mm
	CPFM259A	CPFM259B	CPFM259C	CPFM259D	CPFM259P	1.00mm
	CPFM260A	CPFM260B	CPFM260C	CPFM260D	CPFM260P	2.00mm
Single-Use	CPFS258A	CPFS258B	CPFS258C	CPFS258D	CPFS258P	0.50mm
	CPFS259A	CPFS259B	CPFS259C	CPFS259D	CPFS259P	1.00mm
	CPFS260A	CPFS260B	CPFS260C	CPFS260D	CPFS260P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM264A	CPFM264B	CPFM264C	CPFM264D	CPFM264P	0.50mm
	CPFM265A	CPFM265B	CPFM265C	CPFM265D	CPFM265P	1.00mm
	CPFM266A	CPFM266B	CPFM266C	CPFM266D	CPFM266P	2.00mm
Single-Use	CPFS264A	CPFS264B	CPFS264C	CPFS264D	CPFS264P	0.50mm
	CPFS265A	CPFS265B	CPFS265C	CPFS265D	CPFS265P	1.00mm
	CPFS266A	CPFS266B	CPFS266C	CPFS266D	CPFS266P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
 0.25mm
 0.75mm
 1.25mm
 2.50mm

Premium Forceps Serial
C P F M / S XXX X
 Cosil Multi-Use/Single-Use Length



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM267A	CPFM267B	CPFM267C	CPFM267D	CPFM267P	0.50mm
	CPFM268A	CPFM268B	CPFM268C	CPFM268D	CPFM268P	1.00mm
	CPFM269A	CPFM269B	CPFM269C	CPFM269D	CPFM269P	2.00mm
Single-Use	CPFS267A	CPFS267B	CPFS267C	CPFS267D	CPFS267P	0.50mm
	CPFS268A	CPFS268B	CPFS268C	CPFS268D	CPFS268P	1.00mm
	CPFS269A	CPFS269B	CPFS269C	CPFS269D	CPFS269P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM273A	CPFM273B	CPFM273C	CPFM273D	CPFM273P	0.50mm
	CPFM274A	CPFM274B	CPFM274C	CPFM274D	CPFM274P	1.00mm
	CPFM275A	CPFM275B	CPFM275C	CPFM275D	CPFM275P	2.00mm
Single-Use	CPFS273A	CPFS273B	CPFS273C	CPFS273D	CPFS273P	0.50mm
	CPFS274A	CPFS274B	CPFS274C	CPFS274D	CPFS274P	1.00mm
	CPFS275A	CPFS275B	CPFS275C	CPFS275D	CPFS275P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM270A	CPFM270B	CPFM270C	CPFM270D	CPFM270P	0.50mm
	CPFM271A	CPFM271B	CPFM271C	CPFM271D	CPFM271P	1.00mm
	CPFM272A	CPFM272B	CPFM272C	CPFM272D	CPFM272P	2.00mm
Single-Use	CPFS270A	CPFS270B	CPFS270C	CPFS270D	CPFS270P	0.50mm
	CPFS271A	CPFS271B	CPFS271C	CPFS271D	CPFS271P	1.00mm
	CPFS272A	CPFS272B	CPFS272C	CPFS272D	CPFS272P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM276A	CPFM276B	CPFM276C	CPFM276D	CPFM276P	0.50mm
	CPFM277A	CPFM277B	CPFM277C	CPFM277D	CPFM277P	1.00mm
	CPFM278A	CPFM278B	CPFM278C	CPFM278D	CPFM278P	2.00mm
Single-Use	CPFS276A	CPFS276B	CPFS276C	CPFS276D	CPFS276P	0.50mm
	CPFS277A	CPFS277B	CPFS277C	CPFS277D	CPFS277P	1.00mm
	CPFS278A	CPFS278B	CPFS278C	CPFS278D	CPFS278P	2.00mm

Box of 20 Units

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Premium Forceps Serial
C P F M / S XXX X
Cosil Multi-Use/Single-Use Length



DIATHERMY ELECTROSURGICAL CABLES & PLUGS




DIATHERMY

ELECTROSURGICAL CABLES & PLUGS

Electrosurgery connecting cables serve as essential accessories that seamlessly link our precision instruments to diathermy machines, ensuring a stable and efficient transfer of electrical energy. Designed with high-quality materials for durability and reliability, these cables provide a secure connection that supports uninterrupted functionality and enhanced safety throughout electrosurgical procedures.



The cables are often made with high-grade insulation materials such as silicone or thermoplastic elastomers, which withstand high temperatures and reduce the risk of accidental shocks or burns to surrounding tissues.



Designed to be flexible yet durable, these cables allow smooth maneuverability without compromising signal integrity or reliability, even in demanding surgical settings.



The conductors inside the cables are usually made from high-quality copper or similar conductive materials, minimizing resistance and ensuring accurate energy transfer.

They play a critical role in maintaining safety and effectiveness during procedures by securely transmitting high-frequency electrical currents with minimal loss or interference.

Instruments Side ▶ ▼▼ Generator Side	Flat Plug	Flat Plug	Flat Plug	Two Pool	Two Pin	Two Pin Ang.	Dual Plug	Flat Ang.
	 <p>SIMENS: SUTTER BM50, ERBE T-SERIES: ERBE ACC/ICC 4MM</p>	 <p>CDCM101 - 3 Meter CDCM102 - 5 Meter</p>	 <p>CDCM103 - 3 Meter CDCM104 - 5 Meter</p>	 <p>CDCM105 - 3 Meter CDCM106 - 5 Meter</p>	 <p>CDCM107 - 3 Meter CDCM108 - 5 Meter</p>	 <p>CDCM109 - 3 Meter CDCM110 - 5 Meter</p>	 <p>CDCM111 - 3 Meter CDCM112 - 5 Meter</p>	 <p>CDCM113 - 3 Meter CDCM114 - 5 Meter</p>
 <p>ASSCULAP: WOLF: MARTIN BERCHTOLD: SUTTER BM780</p>	 <p>CDCM117 - 3 Meter CDCM118 - 5 Meter</p>	 <p>CDCM119 - 3 Meter CDCM120 - 5 Meter</p>	 <p>CDCM121 - 3 Meter CDCM122 - 5 Meter</p>	 <p>CDCM123 - 3 Meter CDCM124 - 5 Meter</p>	 <p>CDCM125 - 3 Meter CDCM126 - 5 Meter</p>	 <p>CDCM127 - 3 Meter CDCM128 - 5 Meter</p>	 <p>CDCM129 - 3 Meter CDCM130 - 5 Meter</p>	 <p>CDCM131 - 3 Meter CDCM132 - 5 Meter</p>
 <p>STANDARD INTERNATIONAL 4MM BANANA PLUG FITTING</p>	 <p>CDCM133 - 3 Meter CDCM134 - 5 Meter</p>	 <p>CDCM135 - 3 Meter CDCM136 - 5 Meter</p>	 <p>CDCM137 - 3 Meter CDCM138 - 5 Meter</p>	 <p>CDCM139 - 3 Meter CDCM140 - 5 Meter</p>	 <p>CDCM141 - 3 Meter CDCM142 - 5 Meter</p>	 <p>CDCM143 - 3 Meter CDCM144 - 5 Meter</p>	 <p>CDCM145 - 3 Meter CDCM146 - 5 Meter</p>	 <p>CDCM147 - 3 Meter CDCM148 - 5 Meter</p>
 <p>VALLEYLAB: LAMIDEY: EMC BOWA, ETHICON</p>	 <p>CDCM149 - 3 Meter CDCM150 - 5 Meter</p>	 <p>CDCM151 - 3 Meter CDCM152 - 5 Meter</p>	 <p>CDCM153 - 3 Meter CDCM154 - 5 Meter</p>	 <p>CDCM155 - 3 Meter CDCM156 - 5 Meter</p>	 <p>CDCM157 - 3 Meter CDCM158 - 5 Meter</p>	 <p>CDCM159 - 3 Meter CDCM160 - 5 Meter</p>	 <p>CDCM161 - 3 Meter CDCM162 - 5 Meter</p>	 <p>CDCM163 - 3 Meter CDCM164 - 5 Meter</p>
 <p>EMC: DOLLEY, ELLMAN</p>	 <p>CDCM165 - 3 Meter CDCM166 - 5 Meter</p>	 <p>CDCM167 - 3 Meter CDCM168 - 5 Meter</p>	 <p>CDCM169 - 3 Meter CDCM170 - 5 Meter</p>	 <p>CDCM171 - 3 Meter CDCM172 - 5 Meter</p>	 <p>CDCM173 - 3 Meter CDCM174 - 5 Meter</p>	 <p>CDCM175 - 3 Meter CDCM176 - 5 Meter</p>	 <p>CDCM177 - 3 Meter CDCM178 - 5 Meter</p>	 <p>CDCM179 - 3 Meter CDCM180 - 5 Meter</p>
 <p>VALLEYLAB (AUTOBIPOLAR MODELS ONLY)</p>	 <p>CDCM181 - 3 Meter CDCM182 - 5 Meter</p>	 <p>CDCM183 - 3 Meter CDCM184 - 5 Meter</p>	 <p>CDCM185 - 3 Meter CDCM186 - 5 Meter</p>	 <p>CDCM187 - 3 Meter CDCM188 - 5 Meter</p>	 <p>CDCM189 - 3 Meter CDCM190 - 5 Meter</p>	 <p>CDCM191 - 3 Meter CDCM192 - 5 Meter</p>	 <p>CDCM193 - 3 Meter CDCM194 - 5 Meter</p>	 <p>CDCM195 - 3 Meter CDCM196 - 5 Meter</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp Box of 25 Units

Diathermy Cables Serial
CDCM/S XXX
Cosil Multi-Use/Single-Use

Instruments Side ▶▶				
▼▼ Generator Side	 4mm Female Plugs	 4mm Male Plugs Fitting	 2.4mm Female Ethicon	 Block Connector Fitting
 SIMENS: SUTTER BM50, ERBE T-SERIES: ERBE ACC/CC 4MM	Multi-Use CDCM236 - 3 Meter CDCM237 - 5 Meter	Multi-Use CDCM238 - 3 Meter CDCM239 - 5 Meter	Multi-Use CDCM240 - 3 Meter CDCM241 - 5 Meter	Multi-Use CDCM242 - 3 Meter CDCM243 - 5 Meter
	Single-Use CDCS236 - 3 Meter CDCS237 - 5 Meter	Single-Use CDCS238 - 3 Meter CDCS239 - 5 Meter	Single-Use CDCS240 - 3 Meter CDCS241 - 5 Meter	Single-Use CDCS242 - 3 Meter CDCS243 - 5 Meter
 ASSCULAP: WOLF: MARTIN BERCHTOLD: SUTTER BM780	Multi-Use CDCM244 - 3 Meter CDCM245 - 5 Meter	Multi-Use CDCM246 - 3 Meter CDCM247 - 5 Meter	Multi-Use CDCM248 - 3 Meter CDCM249 - 5 Meter	Multi-Use CDCM250 - 3 Meter CDCM251 - 5 Meter
	Single-Use CDCS244 - 3 Meter CDCS245 - 5 Meter	Single-Use CDCS246 - 3 Meter CDCS247 - 5 Meter	Single-Use CDCS248 - 3 Meter CDCS249 - 5 Meter	Single-Use CDCS250 - 3 Meter CDCS251 - 5 Meter
 STANDARD INTERNATIONAL 4MM BANANA PLUG FITTING	Multi-Use CDCM252 - 3 Meter CDCM253 - 5 Meter	Multi-Use CDCM254 - 3 Meter CDCM255 - 5 Meter	Multi-Use CDCM256 - 3 Meter CDCM257 - 5 Meter	Multi-Use CDCM258 - 3 Meter CDCM259 - 5 Meter
	Single-Use CDCS252 - 3 Meter CDCS253 - 5 Meter	Single-Use CDCS254 - 3 Meter CDCS255 - 5 Meter	Single-Use CDCS256 - 3 Meter CDCS257 - 5 Meter	Single-Use CDCS258 - 3 Meter CDCS259 - 5 Meter
 VALLEYLAB: LAMIDEY: EMC BOWA, ETHICON	Multi-Use CDCM260 - 3 Meter CDCM261 - 5 Meter	Multi-Use CDCM262 - 3 Meter CDCM263 - 5 Meter	Multi-Use CDCM264 - 3 Meter CDCM265 - 5 Meter	Multi-Use CDCM266 - 3 Meter CDCM267 - 5 Meter
	Single-Use CDCS260 - 3 Meter CDCS261 - 5 Meter	Single-Use CDCS262 - 3 Meter CDCS263 - 5 Meter	Single-Use CDCS264 - 3 Meter CDCS265 - 5 Meter	Single-Use CDCS266 - 3 Meter CDCS267 - 5 Meter
 EMC: DOLLEY, ELLMAN	Multi-Use CDCM268 - 3 Meter CDCM269 - 5 Meter	Multi-Use CDCM270 - 3 Meter CDCM271 - 5 Meter	Multi-Use CDCM272 - 3 Meter CDCM273 - 5 Meter	Multi-Use CDCM274 - 3 Meter CDCM275 - 5 Meter
	Single-Use CDCS268 - 3 Meter CDCS269 - 5 Meter	Single-Use CDCS270 - 3 Meter CDCS271 - 5 Meter	Single-Use CDCS272 - 3 Meter CDCS273 - 5 Meter	Single-Use CDCS274 - 3 Meter CDCS275 - 5 Meter
 VALLEYLAB (AUTOBIPOLAR MODELS ONLY)	Multi-Use CDCM276 - 3 Meter CDCM277 - 5 Meter	Multi-Use CDCM278 - 3 Meter CDCM279 - 5 Meter	Multi-Use CDCM280 - 3 Meter CDCM281 - 5 Meter	Multi-Use CDCM282 - 3 Meter CDCM283 - 5 Meter
	Single-Use CDCS276 - 3 Meter CDCS277 - 5 Meter	Single-Use CDCS278 - 3 Meter CDCS279 - 5 Meter	Single-Use CDCS280 - 3 Meter CDCS281 - 5 Meter	Single-Use CDCS282 - 3 Meter CDCS283 - 5 Meter

▼▼ Generator Side

▼▼ Instruments Side

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM284 - 3 Meter
CDCM285 - 5 Meter

Single-Use

CDCS284 - 3 Meter
CDCS285 - 5 Meter

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM286 - 3 Meter
CDCM287 - 5 Meter

Single-Use

CDCS286 - 3 Meter
CDCS287 - 5 Meter

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM288 - 3 Meter
CDCM289 - 5 Meter

Single-Use

CDCS288 - 3 Meter
CDCS289 - 5 Meter

4mm banana to 6mm Alsa Excel connector adapter for Alsa Excelgenerator



Multi-Use

CDCM290 - 3 Meter
CDCM291 - 5 Meter










Single-Use

CDCS290 - 3 Meter
CDCS291 - 5 Meter


Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp Box of 25 Units

Diathermy Cables Serial
C D C M / S XXX
Cosil Multi-Use/Single-Use

Instruments Side ▶▶ ▼▼ Generator Side	 4mm Plug		 Special Edition		 3mm Plug		 4.8mm Plug	
	Multi-Use	Single-Use	Multi-Use	Single-Use	Multi-Use	Single-Use	Multi-Use	Single-Use
 MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CDCM197 - 3 Meter	CDCM198 - 5 Meter	CDCM199 - 3 Meter	CDCM200 - 5 Meter	CDCM201 - 3 Meter	CDCM202 - 5 Meter	CDCM203 - 3 Meter	CDCM204 - 5 Meter
	CDCS197 - 3 Meter	CDCS198 - 5 Meter	CDCS199 - 3 Meter	CDCS200 - 5 Meter	CDCS201 - 3 Meter	CDCS202 - 5 Meter	CDCS203 - 3 Meter	CDCS204 - 5 Meter
 ERBE ACC/ICC	CDCM205 - 3 Meter	CDCM206 - 5 Meter	CDCM207 - 3 Meter	CDCM208 - 5 Meter	CDCM209 - 3 Meter	CDCM210 - 5 Meter	CDCM211 - 3 Meter	CDCM212 - 5 Meter
	CDCS205 - 3 Meter	CDCS206 - 5 Meter	CDCS207 - 3 Meter	CDCS208 - 5 Meter	CDCS209 - 3 Meter	CDCS210 - 5 Meter	CDCS-211 - 3 Meter	CDCS-212 - 5 Meter
 VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CDCM213 - 3 Meter	CDCM214 - 5 Meter	CDCM215 - 3 Meter	CDCM216 - 5 Meter	CDCM217 - 3 Meter	CDCM218 - 5 Meter	CDCM219 - 3 Meter	CDCM220 - 5 Meter
	CDCS213 - 3 Meter	CDCS214 - 5 Meter	CDCS215 - 3 Meter	CDCS216 - 5 Meter	CDCS217 - 3 Meter	CDCS218 - 5 Meter	CDCS219 - 3 Meter	CDCS220 - 5 Meter
 MARTIN: BERCHTOLD:	CDCM221 - 3 Meter	CDCM222 - 5 Meter	CDCM223 - 3 Meter	CDCM224 - 5 Meter	CDCM225 - 3 Meter	CDCM226 - 5 Meter	CDCM227 - 3 Meter	CDCM228 - 5 Meter
	CDCS221 - 3 Meter	CDCS222 - 5 Meter	CDCS223 - 3 Meter	CDCS224 - 5 Meter	CDCS225 - 3 Meter	CDCS226 - 5 Meter	CDCS227 - 3 Meter	CDCS228 - 5 Meter
 INTERNATIONAL 3 PIN PLUG DIATHERMY	CDCM229 - 3 Meter	CDCM230 - 5 Meter	N/A	N/A	N/A	N/A	CDCM231 - 3 Meter	CDCM232 - 5 Meter
	CDCS229 - 3 Meter	CDCS230 - 5 Meter	N/A	N/A	N/A	N/A	CDCS231 - 3 Meter	CDCS232 - 5 Meter

Connecting Finger switch handles with 3 pin plugs, for ERBE T-Series generators.



CDC233

Connecting Finger switch handles with 3 pin plugs, for ERBE ACC / ICC, ERBE VIO generators.



CDC234

Connecting Finger switch handles with 3 pin plugs, for Martin, Integra generators.




CDC235

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max 3vp Box of 25 Units

Diathermy Cables Serial
CDCM/S XXX
 Cosil Multi-Use/Single-Use




NON CONDUCTIVE
**GYNÆCOLOGY
INSTRUMENTS**

NON CONDUCTIVE
**GYNÆCOLOGY
INSTRUMENTS**


Insulated Gynecology Instruments: Our premium quality of instruments, including speculums, retractors, and hooks, are designed to enhance safety and precision in gynecology procedures. With advanced insulation technology, these instruments minimize the risk of electrical burns by preventing unintended electrical contact with surrounding tissues. The ergonomic design ensures ease of handling, while the insulated tips offer reliable performance for delicate procedures, providing optimal control and minimizing patient discomfort.



The locking mechanism holds the blades open at the desired angle, allowing hands-free operation during examination.



The screw enables slight adjustments, ensuring the speculum remains securely in position.



Often curved or ergonomically shaped for comfortable grip and control.





Pederson Vaginal Speculum with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM101	SMALL	CGIS101
CGIM102	MEDIUM	CGIS102
CGIM103	LARGE	CGIS103

BOX OF 1 UNITS



Graves Vaginal Speculum with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM104	SMALL	CGIS104
CGIM105	MEDIUM	CGIS105
CGIM106	LARGE	CGIS106

BOX OF 1 UNITS



Cusco Speculum with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM107	SMALL	CGIS107
CGIM108	MEDIUM	CGIS108
CGIM109	LARGE	CGIS109


BOX OF 1 UNITS



Collin Speculum with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM110	SMALL	CGIS110
CGIM111	MEDIUM	CGIS111
CGIM112	LARGE	CGIS112

BOX OF 1 UNITS



Lateral Vaginal Retractor

Multi-Use		Single-Use
CGIM113	SMALL	CGIS113
CGIM114	MEDIUM	CGIS114
CGIM115	LARGE	CGIS115

BOX OF 1 UNITS



Kogan Endocervical

Multi-Use		Single-Use
CGIM116	SMALL	CGIS116
CGIM117	MEDIUM	CGIS117
CGIM118	LARGE	CGIS118

BOX OF 1 UNITS

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP


Gyne Instruments Serial
CGI M/S XXX
Cosil Multi-Use/Single-Use



Rigby Speculum

Multi-Use		Single-Use
CGIM119	SMALL	CGIS119
CGIM120	MEDIUM	CGIS120
CGIM121	LARGE	CGIS121

BOX OF 1 UNITS



Lateral Retractor

Multi-Use		Single-Use
CGIM122	SMALL	CGIS122
CGIM123	MEDIUM	CGIS123
CGIM124	LARGE	CGIS124

BOX OF 1 UNITS



Champion Forceps
Mat Finishing

Multi-Use		Single-Use
CGIM125	SMALL	CGIS125
CGIM126	MEDIUM	CGIS126
CGIM127	LARGE	CGIS127

BOX OF 1 UNITS



Ring Forceps
Mat Finishing

Multi-Use		Single-Use
CGIM128	SMALL	CGIS128
CGIM129	MEDIUM	CGIS129
CGIM130	LARGE	CGIS130


BOX OF 1 UNITS



Tooth Forceps 1/1
Mat Finishing

Multi-Use		Single-Use
CGIM131	SMALL	CGIS131
CGIM132	MEDIUM	CGIS132
CGIM133	Large	CGIS133

BOX OF 1 UNITS



Gelpi Retractors
Smooth Ratchet Integrated

Multi-Use		Single-Use
CGIM134	SMALL	CGIS134
CGIM135	MEDIUM	CGIS135
CGIM136	Large	CGIS136

BOX OF 1 UNITS

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5KVP

Gyne Instruments Serial
CGI M/S XXX
Cosil Multi-Use/Single-Use



Dressing Forcep

Multi-Use		Single-Use
CGIM137	155mm	CGIS137
CGIM138	180mm	CGIS138
CGIM139	210mm	CGIS139
CGIM140	250mm	CGIS140

BOX OF 20 UNITS



Tissue Forcep

Multi-Use		Single-Use
CGIM141	155mm	CGIS141
CGIM142	180mm	CGIS142
CGIM143	210mm	CGIS143
CGIM144	250mm	CGIS144


BOX OF 20 UNITS



Iris Hook

Multi-Use		Single-Use
CGIM145	155mm	CGIS145
CGIM146	180mm	CGIS146
CGIM147	210mm	CGIS147
CGIM148	250mm	CGIS148


BOX OF 20 UNITS



Straight Hook

Multi-Use		Single-Use
CGIM149	155mm	CGIS149
CGIM150	180mm	CGIS150
CGIM151	210mm	CGIS151
CGIM152	250mm	CGIS152


BOX OF 20 UNITS



Straight Hook 2

Multi-Use		Single-Use
CGIM153	155mm	CGIS153
CGIM154	180mm	CGIS154
CGIM155	210mm	CGIS155
CGIM156	250mm	CGIS156

BOX OF 20 UNITS




Straight Hook 3

Multi-Use		Single-Use
CGIM157	155mm	CGIS157
CGIM158	180mm	CGIS158
CGIM159	210mm	CGIS159
CGIM160	250mm	CGIS160

BOX OF 20 UNITS

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand





DIATHERMY CAUTERY
**HAND PIECE &
ELECTRODES**


DIATHERMY CAUTERY

**HAND PIECE &
ELECTRODES**

The electrosurgical active handpiece delivers high-frequency energy for cutting and coagulating tissue. Ergonomically designed, it supports various electrodes and offers finger or foot control. Available in single-use and reusable options, it ensures precision and efficiency in surgeries.



Some designs incorporate a button or switch near the upper grip for fingertip control, allowing quick activation of various functions (e.g., cauterization, suction). The button is usually recessed slightly to prevent accidental activation.



The upper section often includes a textured or grooved grip surface that enhances handling, ensuring the user's fingers won't slip even in delicate procedures. Some models might use a rubberized or silicone coating for added comfort.



The bottom part often includes a secure locking or threaded mechanism that connects the shaft to the instrument handle, power source, or control console. This may be a twist-lock or snap-fit system, ensuring stability during use and preventing accidental detachment.





BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use	Single-Use
CDPM101 3M	CDPS101
CDPM102 5M	CDPS102



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM103 3M	CDPS103
CDPM104 5M	CDPS104



Valleylab, Bowa, ERBE ICC ERBE Vio, Eschmann, Conmed

Multi-Use	Single-Use
CDPM105 3M	CDPS105
CDPM106 5M	CDPS106



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM107 3M	CDPS107
CDPM108 5M	CDPS108



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use	Single-Use
CDPM109 3M	CDPS109
CDPM110 5M	CDPS110



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM111 3M	CDPS111
CDPM112 5M	CDPS112



Valleylab, Bowa, ERBE ICC ERBE Vio, Eschmann, Conmed

Multi-Use	Single-Use
CDPM113 3M	CDPS113
CDPM114 5M	CDPS114



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM115 3M	CDPS115
CDPM116 5M	CDPS116



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use	Single-Use
CDPM117 3M	CDPS117
CDPM118 5M	CDPS118



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM119 3M	CDPS119
CDPM120 5M	CDPS120



Valleylab, Bowa, ERBE ICC ERBE Vio, Eschmann, Conmed

Multi-Use	Single-Use
CDPM121 3M	CDPS121
CDPM122 5M	CDPS122



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM123 3M	CDPS123
CDPM124 5M	CDPS124



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use	Single-Use
CDPM125 3M	CDPS125
CDPM126 5M	CDPS126



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM127 3M	CDPS127
CDPM128 5M	CDPS128



Valleylab, Bowa, ERBE ICC ERBE Vio, Eschmann, Conmed

Multi-Use	Single-Use
CDPM129 3M	CDPS129
CDPM130 5M	CDPS130



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM131 3M	CDPS131
CDPM132 5M	CDPS132

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Pencils
Serial
CDPM/S XXX
Cosil Multi-Use/Single-Use

Foot Control Cautery Pencil

BOX OF 20 UNITS



4mm International, Bowa, ICC ERBE Vio , Eschmann, Conmed

Multi-Use	Single-Use
CDPM133 3M	CDPS133
CDPM134 5M	CDPS134



8mm ERBE T-Series Valleylab, ERBE Int'l

Multi-Use	Single-Use
CDPM135 3M	CDPS135
CDPM136 5M	CDPS136



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM137 3M	CDPS137
CDPM138 5M	CDPS138



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM139 3M	CDPS139
CDPM140 5M	CDPS140

BOX OF 20 UNITS



4mm International, Bowa, ICC ERBE Vio , Eschmann, Conmed

Multi-Use	Single-Use
CDPM141 3M	CDPS141
CDPM142 5M	CDPS142



8mm ERBE T-Series Valleylab, ERBE Int'l

Multi-Use	Single-Use
CDPM143 3M	CDPS143
CDPM144 5M	CDPS144



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM145 3M	CDPS145
CDPM146 5M	CDPS146



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM147 3M	CDPS147
CDPM148 5M	CDPS148

BOX OF 20 UNITS



4mm International, Bowa, ICC ERBE Vio , Eschmann, Conmed

Multi-Use	Single-Use
CDPM149 3M	CDPS149
CDPM150 5M	CDPS150



8mm ERBE T-Series Valleylab, ERBE Int'l

Multi-Use	Single-Use
CDPM151 3M	CDPS151
CDPM152 5M	CDPS152



5mm ERBE ACC Series Standard ICC Series Vio, Karl Storz

Multi-Use	Single-Use
CDPM153 3M	CDPS153
CDPM154 5M	CDPS154



Martin, ERBE T-Series, Integra Berchtold, Aesculap

Multi-Use	Single-Use
CDPM155 3M	CDPS155
CDPM156 5M	CDPS156

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Pencils Serial
C D P M / S XXX
 Cosil Multi-Use/Single-Use

Foot Control Cautery Pencil



BOX OF 20 UNITS

BOX OF 20 UNITS

MARTIN: BERCHTOLD: ERBE T-SERIES, Aesculap, Sutter		VALLEYLAB: CONMED: BOVIE 8mm plug, Sowa, Sard	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM157	3M	CDPS157	
CDPM158	5M	CDPS158	

MARTIN: BERCHTOLD: ERBE T-SERIES, Aesculap, Sutter		VALLEYLAB: CONMED: BOVIE 8mm plug, Sowa, Sard	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM167	3M	CDPS167	
CDPM168	5M	CDPS168	



INTERNATIONAL 3 PIN PLUG DIATHERMY

INTERNATIONAL 3 PIN PLUG DIATHERMY

ERBE ACC/ICC		MARTIN: BERCHTOLD:	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM161	3M	CDPS161	
CDPM162	5M	CDPS162	

ERBE ACC/ICC		MARTIN: BERCHTOLD:	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM171	3M	CDPS171	
CDPM172	5M	CDPS172	

Multi-Use	Single-Use	
CDPM165	3M	CDPS165
CDPM166	5M	CDPS166

Multi-Use	Single-Use	
CDPM175	3M	CDPS175
CDPM176	5M	CDPS176



BOX OF 20 UNITS

MARTIN: BERCHTOLD: ERBE T-SERIES, Aesculap, Sutter		VALLEYLAB: CONMED: BOVIE 8mm plug, Sowa, Sard	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM177	3M	CDPS177	
CDPM178	5M	CDPS178	

ERBE ACC/ICC		MARTIN: BERCHTOLD:	
Multi-Use	Single-Use	Multi-Use	Single-Use
CDPM181	3M	CDPS181	
CDPM182	5M	CDPS182	

















INTERNATIONAL 3 PIN PLUG DIATHERMY

Multi-Use	Single-Use	
CDPM185	3M	CDPS185
CDPM186	5M	CDPS186

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

















Diathermy Pencils Serial
C D P M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM190</p> <p>Single-Use CDES190</p>	 <p>0.5x20mm</p> <p>Needle electrodes, Straight</p>	<p>4mm Shaft</p>	 <p>3.0x25mm</p> <p>Knife Electrodes, Straight</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM191</p> <p>Single-Use CDES191</p>
<p>Multi-Use CDEM191</p> <p>Single-Use CDES192</p>	 <p>2.5x25mm</p> <p>Lancet electrodes, Straight</p>	<p>4mm Shaft</p>	 <p>2.5x20mm</p> <p>Blade Electrodes, Straight</p>	<p>Multi-Use CDEM193</p> <p>Single-Use CDES193</p>
<p>Multi-Use CDEM194</p> <p>Single-Use CNES194</p>	 <p>5 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>4mm Shaft</p>	 <p>10 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>Multi-Use CDEM195</p> <p>Single-Use CDES195</p>
<p>Multi-Use CDEM196</p> <p>Single-Use CDES196</p>	 <p>15 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>4mm Shaft</p>	 <p>2 mm</p> <p>Ball Electrodes, Straight</p>	<p>Multi-Use CDEM197</p> <p>Single-Use CDES197</p>
<p>Multi-Use CDEM198</p> <p>Single-Use CDES198</p>	 <p>4 mm</p> <p>Ball Electrodes, Straight</p>	<p>4mm Shaft</p>	 <p>6 mm</p> <p>Ball Electrodes, Straight</p>	<p>Multi-Use CDEM199</p> <p>Single-Use CDES199</p>
<p>Multi-Use CDEM200</p> <p>Single-Use CDES200</p>	 <p>0.5x20mm</p> <p>Needle electrodes, Straight</p>	<p>2.4mm Shaft</p>	 <p>3.0x25mm</p> <p>Knife Electrodes, Straight</p>	<p>Multi-Use CDEM201</p> <p>Single-Use CDES201</p>
<p>Multi-Use CDEM202</p> <p>Single-Use CDES202</p>	 <p>2.5x25mm</p> <p>Lancet Electrodes, Straight</p>	<p>2.4mm Shaft</p>	 <p>2.5x20mm</p> <p>Blade Electrodes, Straight</p>	<p>Multi-Use CDEM203</p> <p>Single-Use CDES203</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
CDEM/S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM204</p> <p>Single-Use CDES204</p>	 <p>5 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>2.4mm Shaft</p>	 <p>10 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM205</p> <p>Single-Use CDES205</p>
<p>Multi-Use CDEM206</p> <p>Single-Use CDES206</p>	 <p>15 mm</p> <p>Wire Loop electrodes, Straight</p>	<p>2.4mm Shaft</p>	 <p>2 mm</p> <p>Ball Electrodes, Straight</p>	<p>Multi-Use CDEM207</p> <p>Single-Use CDES207</p>
<p>Multi-Use CDEM208</p> <p>Single-Use CNES208</p>	 <p>4 mm</p> <p>Ball Electrodes, Straight</p>	<p>2.4mm Shaft</p>	 <p>6 mm</p> <p>Ball Electrodes, Straight</p>	<p>Multi-Use CDEM209</p> <p>Single-Use CDES1209</p>
<p>Multi-Use CDEM210</p> <p>Single-Use CDES210</p>	 <p>0.5 mm</p> <p>Needle electrodes, Straight</p>	<p>2.38mm Shaft</p>	 <p>2.5x20mm</p> <p>Blade Electrodes, Straight</p>	<p>Multi-Use CDEM211</p> <p>Single-Use CDES211</p>
<p>Multi-Use CDEM212</p> <p>Single-Use CDES212</p>	 <p>0.3mm</p> <p>Needle electrodes, Straight</p>	<p>2.38mm Shaft</p>	 <p>0.5 mm</p> <p>Needle Electrodes, Straight</p>	<p>Multi-Use CDEM213</p> <p>Single-Use CDES213</p>
<p>Multi-Use CDEM214</p> <p>Single-Use CDES214</p>	 <p>0.3 mm</p> <p>Needle electrodes, Straight</p>	<p>2.38mm Shaft</p>	 <p>0.5 mm</p> <p>Needle electrodes, Straight</p>	<p>Multi-Use CDEM215</p> <p>Single-Use CDES215</p>
<p>Multi-Use CDEM216</p> <p>Single-Use CDES216</p>	 <p>0.5 mm</p> <p>Needle electrodes, Straight</p>	<p>2.38mm Shaft</p>	 <p>2 mm</p> <p>Ball Electrodes, Straight</p>	<p>Multi-Use CDEM217</p> <p>Single-Use CDES217</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 4/1 KVP















Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM218</p> <p>Single-Use CDES218</p>	 <p>Needle Electrodes, Insulated 45° Angled</p>	<p>1.6mm Shaft</p>	 <p>Needle Electrodes, Insulated 45° Angled</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM219</p> <p>Single-Use CDES219</p>
<p>Multi-Use CDEM220</p> <p>Single-Use CDES220</p>	 <p>Needle Electrodes, Insulated 45° Angled</p>	<p>1.6mm Shaft</p>	 <p>Needle Electrodes, Insulated 45° Angled</p>	<p>Multi-Use CDEM221</p> <p>Single-Use CDES221</p>
<p>Multi-Use CDEM222</p> <p>Single-Use CNES222</p>	 <p>Needle Electrodes, Insulated 90° Angled</p>	<p>1.6mm Shaft</p>	 <p>Wire Loop electrodes, Straight Triangle</p>	<p>Multi-Use CDEM223</p> <p>Single-Use CDES223</p>
<p>Multi-Use CDEM224</p> <p>Single-Use CDES224</p>	 <p>Wire Loop electrodes</p>	<p>1.6mm Shaft</p>	 <p>Wire Loop electrodes</p>	<p>Multi-Use CDEM225</p> <p>Single-Use CDES225</p>
<p>Multi-Use CDEM226</p> <p>Single-Use CDES226</p>	 <p>Wire Loop electrodes</p>	<p>1.6mm Shaft</p>		

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 1KVP















Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM232</p> <p>Single-Use CDES232</p>	 <p>Knife Electrodes</p>	<p>4mm Shaft</p>	 <p>Lancet-Cutting Electrodes</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM233</p> <p>Single-Use CDES233</p>
<p>Multi-Use CDEM234</p> <p>Single-Use CDES234</p>	 <p>Lancet-Cutting Electrodes</p>	<p>4mm Shaft</p>	 <p>Lancet-Cutting Electrodes</p>	<p>Multi-Use CDEM235</p> <p>Single-Use CDES235</p>
<p>Multi-Use CDEM236</p> <p>Single-Use CDES236</p>	 <p>Spatula Electrodes, Oval Shape</p>	<p>4mm Shaft</p>	 <p>Knife Electrodes, Classic</p>	<p>Multi-Use CDEM237</p> <p>Single-Use CDES237</p>
<p>Multi-Use CDEM238</p> <p>Single-Use CDES238</p>	 <p>Knife Electrodes</p>	<p>4mm Shaft</p>	 <p>Knife Electrodes, Angular</p>	<p>Multi-Use CDEM239</p> <p>Single-Use CDES239</p>
<p>Multi-Use CDEM240</p> <p>Single-Use CDES240</p>	 <p>Knife Electrode, Lancet</p>	<p>4mm Shaft</p>	 <p>Blade Electrodes</p>	<p>Multi-Use CDEM241</p> <p>Single-Use CDES241</p>
<p>Multi-Use CDEM242</p> <p>Single-Use CDES242</p>	 <p>Needle Electrodes, Insulated 10MM WL</p>	<p>4mm Shaft</p>	 <p>Needle Electrodes</p>	<p>Multi-Use CDEM243</p> <p>Single-Use CDES243</p>
<p>Multi-Use CDEM244</p> <p>Single-Use CDES244</p>	 <p>Needle Electrode</p>	<p>4mm Shaft</p>	 <p>Needle Electrode</p>	<p>Multi-Use CDEM245</p> <p>Single-Use CDES245</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4 KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM246</p> <p>Single-Use CDES246</p>	 <p>0.8X20mm</p> <p>Needle Electrode</p>	<p>4mm Shaft</p>	 <p>0.5X10mm</p> <p>Needle Electrode</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM247</p> <p>Single-Use CDES247</p>
<p>Multi-Use CDEM248</p> <p>Single-Use CDES248</p>	 <p>0.5X20mm</p> <p>Needle Electrode</p>	<p>4mm Shaft</p>	 <p>5mm</p> <p>Wire Electrode</p>	<p>Multi-Use CDEM249</p> <p>Single-Use CDES249</p>
<p>Multi-Use CDEM250</p> <p>Single-Use CDES250</p>	 <p>10mm</p> <p>Wire Electrode</p>	<p>4mm Shaft</p>	 <p>10mm</p> <p>Ribbon Loop Electrode</p>	<p>Multi-Use CDEM251</p> <p>Single-Use CDES251</p>
<p>Multi-Use CDEM252</p> <p>Single-Use CDES252</p>	 <p>15mm</p> <p>Ribbon Loop Electrode</p>	<p>4mm Shaft</p>	 <p>2mm</p> <p>Ball Electrode</p>	<p>Multi-Use CDEM253</p> <p>Single-Use CDES253</p>
<p>Multi-Use CDEM254</p> <p>Single-Use CDES254</p>	 <p>4mm</p> <p>Ball Electrode</p>	<p>4mm Shaft</p>	 <p>6mm</p> <p>Ball Electrode</p>	<p>Multi-Use CDEM255</p> <p>Single-Use CDES255</p>
<p>Multi-Use CDEM256</p> <p>Single-Use CDES256</p>	 <p>2mm</p> <p>Ball Electrode, Angular</p>	<p>4mm Shaft</p>	 <p>4mm</p> <p>Ball Electrode, Angular</p>	<p>Multi-Use CDEM257</p> <p>Single-Use CDES257</p>
<p>Multi-Use CDEM258</p> <p>Single-Use CDES258</p>	 <p>0.5X5mm</p> <p>Tungsten Electrode, Straight</p>	<p>4mm Shaft</p>	 <p>0.5X5mm</p> <p>Tungsten Electrode, Angle</p>	<p>Multi-Use CDEM259</p> <p>Single-Use CDES259</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP















Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM260</p> <p>Single-Use CDES260</p>	<p>0.2mm, 7mm</p> <p>Wire Loop Electrode, Triangular</p>	<p>4mm Shaft</p>	<p>0.2mm, 9mm</p> <p>Wire Loop Electrode, Triangular</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM261</p> <p>Single-Use CDES261</p>
<p>Multi-Use CDEM262</p> <p>Single-Use CDES262</p>	<p>0.2mm</p> <p>Wire Loop Electrode, Rhombic</p>	<p>4mm Shaft</p>	<p>5mm</p> <p>Wire Loop Electrode</p>	<p>Multi-Use CDEM263</p> <p>Single-Use CDES263</p>
<p>Multi-Use CDEM264</p> <p>Single-Use CNES264</p>	<p>10 mm</p> <p>Wire Loop Electrode</p>	<p>4mm Shaft</p>	<p>15 mm</p> <p>Wire Loop Electrode</p>	<p>Multi-Use CDEM265</p> <p>Single-Use CDES265</p>
<p>Multi-Use CDEM266</p> <p>Single-Use CDES266</p>	<p>0.25mm</p> <p>Needle Electrode</p>	<p>4mm Shaft</p>	<p>1mm</p> <p>Ball Electrode</p>	<p>Multi-Use CDEM267</p> <p>Single-Use CDES267</p>
<p>Multi-Use CDEM268</p> <p>Single-Use CDES268</p>	<p>2mm</p> <p>Ball Electrode</p>	<p>4mm Shaft</p>	<p>4mm</p> <p>Ball Electrode</p>	<p>Multi-Use CDEM269</p> <p>Single-Use CDES269</p>
<p>Multi-Use CDEM270</p> <p>Single-Use CDES270</p>	<p>6mm</p> <p>Ball Electrode</p>	<p>4mm Shaft</p>	<p>2X8mm</p> <p>Spatula Electrode</p>	<p>Multi-Use CDEM271</p> <p>Single-Use CDES271</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand





95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM272</p> <p>Single-Use CDES272</p>	 <p>2mm</p> <p>Ball electrode</p>	<p>4mm Shaft</p>	 <p>4mm</p> <p>Ball electrode</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM273</p> <p>Single-Use CDES273</p>
<p>Multi-Use CDEM274</p> <p>Single-Use CDES274</p>	 <p>6mm</p> <p>Ball electrode</p>	<p>4mm Shaft</p>	 <p>5mm</p> <p>Wire loop electrode</p>	<p>Multi-Use CDEM275</p> <p>Single-Use CDES275</p>
<p>Multi-Use CDEM276</p> <p>Single-Use CDES276</p>	 <p>8mm</p> <p>Wire loop electrode</p>	<p>4mm Shaft</p>	 <p>10mm</p> <p>Wire loop electrode</p>	<p>Multi-Use CDEM277</p> <p>Single-Use CDES277</p>
<p>Multi-Use CDEM278</p> <p>Single-Use CDES278</p>	 <p>5mm</p> <p>Ribbon loop electrode</p>	<p>4mm Shaft</p>	 <p>10mm 15mm</p> <p>Ribbon loop electrode</p>	<p>Multi-Use CDEM279</p> <p>Single-Use CDES279</p>
<p>Multi-Use CDEM280</p> <p>Single-Use CDES280</p>	 <p>15mm</p> <p>Ribbon loop electrode</p>	<p>4mm Shaft</p>	 <p>1.6mm, 19mm</p> <p>Lancet electrode, Stright</p>	<p>Multi-Use CDEM281</p> <p>Single-Use CDES281</p>
<p>Multi-Use CDEM282</p> <p>Single-Use CDES282</p>	 <p>1.6x19mm</p> <p>LANCET ELECTRODE ANGULAR</p>	<p>4mm Shaft</p>	 <p>3.3x25mm</p> <p>Knife electrode</p>	<p>Multi-Use CDEM283</p> <p>Single-Use CDES283</p>
<p>Multi-Use CDEM284</p> <p>Single-Use CDES284</p>	 <p>3.3x24mm</p> <p>Knife electrode</p>	<p>4mm Shaft</p>	 <p>22mm, 30mm</p> <p>Needle electrode</p>	<p>Multi-Use CDEM285</p> <p>Single-Use CDES285</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand








<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM286 Single-Use CDES286</p>	 <p>4mm</p> <p>Ball electrode</p>	<p>2.4mm Shaft</p>	 <p>2mm</p> <p>Ball electrode</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM287 Single-Use CDES287</p>
<p>Multi-Use CDEM288 Single-Use CDES288</p>	 <p>4mm</p> <p>Ball electrode</p>	<p>2.4mm Shaft</p>	 <p>6mm</p> <p>Ball electrode</p>	<p>Multi-Use CDEM289 Single-Use CDES289</p>
<p>Multi-Use CDEM290 Single-Use CNES290</p>	 <p>2mm</p> <p>Ball electrode</p>	<p>2.4mm Shaft</p>	 <p>4mm</p>	<p>Multi-Use CDEM291 Single-Use CDES291</p>
<p>Multi-Use CDEM292 Single-Use CDES292</p>	 <p>6mm</p>	<p>2.4mm Shaft</p>	 <p>2mm</p> <p>Ball electrode</p>	<p>Multi-Use CDEM293 Single-Use CDES293</p>
<p>Multi-Use CDEM294 Single-Use CDES294</p>	 <p>2.5X20mm</p> <p>Blade Electrode</p>	<p>2.4mm Shaft</p>	 <p>2.5X20mm</p> <p>Spatula electrode</p>	<p>Multi-Use CDEM295 Single-Use CDES295</p>
<p>Multi-Use CDEM296 Single-Use CDES296</p>	 <p>0.8X20mm</p> <p>Needle electrode</p>	<p>2.4mm Shaft</p>	 <p>0.5X20mm</p> <p>Needle electrode</p>	<p>Multi-Use CDEM297 Single-Use CDES297</p>
<p>Multi-Use CDEM298 Single-Use CDES298</p>	 <p>0.5X20mm</p> <p>Tungsten electrode, Angled</p>	<p>2.4mm Shaft</p>	 <p>0.5X20mm</p> <p>Tungsten electrode</p>	<p>Multi-Use CDEM299 Single-Use CDES299</p>


Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 2KVP


Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM300</p> <p>Single-Use CDES300</p>	 <p>0.5mm 5mm</p> <p>Needle electrode with ceramic isolating & Tungsten needle</p>	<p>2.4mm Shaft</p>	 <p>0.5X20mm</p> <p>Needle electrode</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM301</p> <p>Single-Use CDES301</p>
<p>Multi-Use CDEM302</p> <p>Single-Use CDES302</p>	 <p>0.2mm 7mm</p> <p>Wire loop electrode, Triangular</p>	<p>2.4mm Shaft</p>	 <p>0.2mm</p> <p>Wire loop electrode, Rhombic</p>	<p>Multi-Use CDEM303</p> <p>Single-Use CDES303</p>
<p>Multi-Use CDEM304</p> <p>Single-Use CDES304</p>	 <p>2X8mm</p> <p>Spatula electrode</p>	<p>2.4mm Shaft</p>	 <p>0.2mm 9mm</p> <p>Wire loop electrode, Triangular</p>	<p>Multi-Use CDEM305</p> <p>Single-Use CDES305</p>
<p>Multi-Use CDEM306</p> <p>Single-Use CDES306</p>	 <p>10mm</p> <p>Wire loop electrode</p>	<p>2.4mm Shaft</p>	 <p>5mm</p> <p>Wire loop electrode</p>	<p>Multi-Use CDEM307</p> <p>Single-Use CDES307</p>
<p>Multi-Use CDEM308</p> <p>Single-Use CDES308</p>	 <p>5mm</p> <p>Ribbon loop electrode</p>	<p>2.4mm Shaft</p>	 <p>15mm</p> <p>Wire loop electrode</p>	<p>Multi-Use CDEM309</p> <p>Single-Use CDES309</p>
<p>Multi-Use CDEM310</p> <p>Single-Use CDES310</p>	 <p>15mm</p> <p>Ribbon loop electrode</p>	<p>2.4mm Shaft</p>	 <p>10mm</p> <p>Ribbon loop electrode</p>	<p>Multi-Use CDEM311</p> <p>Single-Use CDES311</p>
<p>Multi-Use CDEM312</p> <p>Single-Use CDES312</p>	 <p>2mm</p> <p>Ball electrode</p>	<p>2.4mm Shaft</p>	 <p>4mm</p> <p>Ball electrode</p>	<p>Multi-Use CDEM313</p> <p>Single-Use CDES313</p>


Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand











Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM314</p> <p>Single-Use CDES314</p>	 <p>Electrode coagulating and fluguration</p> <p>1.6mm Shaft</p>	 <p>Wire loop electrode 3mm</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM315</p> <p>Single-Use CDES315</p>
<p>Multi-Use CDEM316</p> <p>Single-Use CDES316</p>	 <p>Wire loop electrode</p> <p>1.6mm Shaft</p>	 <p>Wire loop electrode</p>	<p>Multi-Use CDEM317</p> <p>Single-Use CDES317</p>
<p>Multi-Use CDEM318</p> <p>Single-Use CNES318</p>	 <p>Needle electrode wired WL 5mm</p> <p>1.6mm Shaft</p>	 <p>Needle electrode wired, Angled 45° WL 5mm</p>	<p>Multi-Use CDEM319</p> <p>Single-Use CDES319</p>
<p>Multi-Use CDEM320</p> <p>Single-Use CDES320</p>	 <p>Coagulation Needle electrode, Angled 45° WL 5mm</p> <p>1.6mm Shaft</p>	 <p>Ball electrode, Angled 45°</p>	<p>Multi-Use CDEM321</p> <p>Single-Use CDES321</p>
<p>Multi-Use CDEM322</p> <p>Single-Use CDES322</p>	 <p>Ball electrode</p> <p>1.6mm Shaft</p>	 <p>Ball electrode</p>	<p>Multi-Use CDEM323</p> <p>Single-Use CDES323</p>
<p>Multi-Use CDEM324</p> <p>Single-Use CDES324</p>	 <p>Ball electrode</p> <p>1.6mm Shaft</p>	 <p>Ball electrode</p>	<p>Multi-Use CDEM325</p> <p>Single-Use CDES325</p>
<p>Multi-Use CDEM326</p> <p>Single-Use CDES326</p>	 <p>Ball electrode, Angled 45°</p> <p>2.38mm Shaft</p>	 <p>Wire loop electrode, Triangular</p>	<p>Multi-Use CDEM327</p> <p>Single-Use CDES327</p>


Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

Multi-Use CDEM328	Single-Use CNES328	 <p>0.2mm 9mm</p> <p>Wire loop electrode, Triangular</p>	2.38mm Shaft	 <p>5mm</p> <p>Wire loop electrode</p>	Multi-Use CDEM329	Single-Use CNES329
Multi-Use CDEM330	Single-Use CNES330	 <p>10mm</p> <p>Wire loop electrode, Rohimbic</p>	2.38mm Shaft	 <p>15mm</p> <p>Wire loop electrode</p>	Multi-Use CDEM331	Single-Use CNES331
Multi-Use CDEM332	Single-Use CNES332	 <p>0.2mm</p> <p>Wire loop electrode, Rohimbic</p>	2.38mm Shaft	 <p>1mm</p> <p>Ball electrode</p>	Multi-Use CDEM333	Single-Use CNES333
Multi-Use CDEM334	Single-Use CNES334	 <p>2X8mm</p> <p>Spatula electrode</p>	2.38mm Shaft	 <p>0.25mm</p> <p>Needle electrodes, Straight</p>	Multi-Use CDEM335	Single-Use CNES335

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM336</p> <p>Single-Use CDES336</p>	 <p>160mm</p> <p>5mm</p> <p>Wire loop electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM337</p> <p>Single-Use CDES337</p>	 <p>160mm</p> <p>10mm</p> <p>Wire loop electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM338</p> <p>Single-Use CDES338</p>	 <p>160mm</p> <p>15mm</p> <p>Wire loop electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM339</p> <p>Single-Use CDES339</p>	 <p>140mm</p> <p>2mm</p> <p>Ball electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM340</p> <p>Single-Use CDES340</p>	 <p>140mm</p> <p>4mm</p> <p>Ball electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM341</p> <p>Single-Use CDES341</p>	 <p>140mm</p> <p>8mm</p> <p>Wire loop electrode</p> <p>4mm Shaft</p>
<p>Multi-Use CDEM342</p> <p>Single-Use CNES342</p>	 <p>140mm</p> <p>0.8X12mm</p> <p>Needle electrodes</p> <p>4mm Shaft</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

BOX OF 50 UNITS		140mm		
Multi-Use CDEM343	Single-Use CDES343		4mm Shaft	Needle electrode, Angled 0.8X12mm
		140mm		
Multi-Use CDEM344	Single-Use CDES344		4mm Shaft	Knife electrode 2.4X10mm
		140mm		
Multi-Use CDEM345	Single-Use CDES345		4mm Shaft	Knife electrode, Angled 2.4X10mm
		140mm		
Multi-Use CDEM346	Single-Use CDES346		4mm Shaft	Knife electrode 3.3X25mm
		140mm		
Multi-Use CDEM347	Single-Use CDES347		4mm Shaft	Knife Electrode 3.3X20mm
		140mm		
Multi-Use CDEM348	Single-Use CDES348		4mm Shaft	Knife Electrode 2.4X20mm
		140mm		
Multi-Use CDEM349	Single-Use CNES349		4mm Shaft	Lancet Electrode 16X20mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

Multi-Use CDEM350	Single-Use CDES350	140mm	4mm Shaft	2.2X20mm	Knife Electrode
Multi-Use CDEM351	Single-Use CDES351	100mm	4mm Shaft	2X16mm	Spatula Electrode
Multi-Use CDEM352	Single-Use CDES352	160mm	4mm Shaft	10X10mm	Loop Electrode
Multi-Use CDEM353	Single-Use CDES353	160mm	4mm Shaft	15X15mm	Loop Electrode
Multi-Use CDEM354	Single-Use CDES354	160mm	4mm Shaft	15X20mm	Loop Electrode
Multi-Use CDEM355	Single-Use CDES355	160mm	4mm Shaft	20X20mm	Loop Electrode
Multi-Use CDEM356	Single-Use CDES356	160mm	4mm Shaft	20X25mm	Loop Electrode

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM357</p> <p>Single-Use CDES357</p>	 <p>160mm</p> <p>10mm</p> <p>4mm Shaft</p> <p>Loop Electrode</p>	<p>BOX OF 50 UNITS</p>
<p>Multi-Use CDEM358</p> <p>Single-Use CDES358</p>	 <p>160mm</p> <p>15mm</p> <p>4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM39</p> <p>Single-Use CDES359</p>	 <p>160mm</p> <p>20mm</p> <p>4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM3560</p> <p>Single-Use CDES360</p>	 <p>160mm</p> <p>25mm</p> <p>4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM361</p> <p>Single-Use CDES3561</p>	 <p>160mm</p> <p>30mm</p> <p>4mm Shaft</p> <p>Loop Electrode</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM362</p> <p>Single-Use CDES362</p>	<p>160mm</p> <p>4mm Shaft</p> <p>10X25mm</p> <p>Loop Electrode</p>	<p>BOX OF 50 UNITS</p>
<p>Multi-Use CDEM363</p> <p>Single-Use CDES363</p>	<p>160mm</p> <p>4mm Shaft</p> <p>15X25mm</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM364</p> <p>Single-Use CDES364</p>	<p>160mm</p> <p>4mm Shaft</p> <p>20X25mm</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM365</p> <p>Single-Use CDES365</p>	<p>140mm</p> <p>4mm Shaft</p> <p>1X10mm</p> <p>Coagulation Needle electrodes</p>	
<p>Multi-Use CDEM366</p> <p>Single-Use CDES366</p>	<p>140mm</p> <p>4mm Shaft</p> <p>4mm</p> <p>Coagulation Ball electrodes</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5/4KVP

Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM370</p> <p>Single-Use CDES370</p>	<p>100mm, 150mm</p> <p>4mm Extension To Insert 2.4mm Electrode</p> <p>Shaft With Cage</p> <p>4mm Shaft</p>		
<p>Multi-Use CDEM371</p> <p>Single-Use CDES371</p>	<p>4mm 2.4mm</p> <p>4mm Extension To Insert 1.6mm Electrode</p> <p>4mm Shaft</p>	<p>4mm 1.6mm</p> <p>4mm Extension To Insert 2.4mm Electrode</p>	<p>Multi-Use CDEM372</p> <p>Single-Use CDES372</p>
<p>Multi-Use CDEM373</p> <p>Single-Use CNES373</p>	<p>100mm, 150mm</p> <p>2.4mm Extension To Insert 4mm Electrode</p> <p>Shaft With Cage</p> <p>4mm Shaft</p>		
<p>Multi-Use CDEM374</p> <p>Single-Use CDES374</p>	<p>100mm, 150mm</p> <p>2.4mm Extension To Insert 1.6mm Electrode</p> <p>Shaft With Cage</p>		
<p>Multi-Use CDEM375</p> <p>Single-Use CDES375</p>	<p>140mm</p> <p>Spatula Electrode</p> <p>2.4mm Shaft</p> <p>2x1.6mm 2.4x20mm</p>		
<p>Multi-Use CDEM376</p> <p>Single-Use CDES376</p>	<p>140mm</p> <p>Knife Electrode</p> <p>2.4mm Shaft</p> <p>2.4x10mm 3.3x25mm</p>		
<p>Multi-Use CDEM377</p> <p>Single-Use CDES377</p>	<p>140mm</p> <p>Needle Electrode</p> <p>2.4mm Shaft</p> <p>0.8x12mm 0.5x20mm</p>		

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 4KVP

Diathermy Electrodes Serial
CDEM/S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM378</p> <p>Single-Use CDES378</p>	<p>140mm</p> <p>2mm</p> <p>2.4mm Shaft</p> <p>Ball electrode</p>	
<p>Multi-Use CDEM379</p> <p>Single-Use CDES379</p>	<p>140mm</p> <p>4mm</p> <p>2.4mm Shaft</p> <p>Ball electrode</p>	
<p>Multi-Use CDEM380</p> <p>Single-Use CNES380</p>	<p>140mm</p> <p>6mm</p> <p>2.4mm Shaft</p> <p>Ball electrode</p>	
<p>Multi-Use CDEM381</p> <p>Single-Use CDES381</p>	<p>140mm</p> <p>8mm</p> <p>2.4mm Shaft</p> <p>Ball electrode</p>	
<p>Multi-Use CDEM382</p> <p>Single-Use CDES382</p>	<p>160mm</p> <p>10x25mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM383</p> <p>Single-Use CDES383</p>	<p>160mm</p> <p>15x25mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM384</p> <p>Single-Use CDES384</p>	<p>160mm</p> <p>20x25mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 4/0.5KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM385 Single-Use CDES385</p>	<p>160mm</p> <p>10x10mm</p> <p>2.4mm Shaft Loop Electrode</p>	
<p>Multi-Use CDEM386 Single-Use CDES386</p>	<p>160mm</p> <p>15x15mm</p> <p>2.4mm Shaft Loop Electrode</p>	
<p>Multi-Use CDEM387 Single-Use CNES387</p>	<p>160mm</p> <p>15x20mm</p> <p>2.4mm Shaft Loop Electrode</p>	
<p>Multi-Use CDEM388 Single-Use CDES388</p>	<p>160mm</p> <p>20x20mm</p> <p>2.4mm Shaft Loop Electrode</p>	
<p>Multi-Use CDEM389 Single-Use CDES389</p>	<p>160mm</p> <p>20x25mm</p> <p>2.4mm Shaft Loop Electrode</p>	
<p>Multi-Use CDEM390 Single-Use CDES390</p>	<p>160mm</p> <p>20x25mm</p> <p>2.4mm Shaft Loop Electrode</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP















Diathermy Electrodes Serial
C D E M / S X X X
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS Multi-Use CDEM391 Single-Use CDES391</p>	<p>160mm</p>  <p>10mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM392 Single-Use CDES392</p>	<p>160mm</p>  <p>15mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM393 Single-Use CNES393</p>	<p>160mm</p>  <p>20mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	
<p>Multi-Use CDEM394 Single-Use CDES394</p>	<p>160mm</p>  <p>25mm</p> <p>2.4mm Shaft</p> <p>Loop Electrode</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 0.5KVP







Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM398</p> <p>Single-Use CDES398</p>	 <p>Microdissection Needle Electrode, Flexible, 45°, Tungsten</p> <p>2.4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM399</p> <p>Single-Use CDES399</p>
<p>Multi-Use CDEM400</p> <p>Single-Use CDES400</p>	 <p>Microdissection Needle Electrode, Flexible, 45°, Tungsten</p> <p>2.4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	<p>Multi-Use CDEM401</p> <p>Single-Use CDES401</p>
<p>Multi-Use CDEM402</p> <p>Single-Use CDES402</p>	 <p>Microdissection Needle Electrode, Flexible, 45°, Tungsten</p> <p>2.4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	
<p>Multi-Use CDEM403</p> <p>Single-Use CDES403</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>2.4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	
<p>Multi-Use CDEM404</p> <p>Single-Use CDES404</p>	 <p>Microdissection Needle Electrode, Flexible, 45°, Tungsten</p> <p>2.4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	<p>Multi-Use CDEM405</p> <p>Single-Use CDES405</p>
<p>Multi-Use CDEM406</p> <p>Single-Use CDES406</p>	 <p>Microdissection Needle Electrode, Flexible, 45°, Tungsten</p> <p>4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	
<p>Multi-Use CDEM407</p> <p>Single-Use CDES407</p>	 <p>Microdissection Needle Electrode, Needle Flexible, Tungsten</p> <p>4mm Shaft</p> <p>0.25x4mm</p>	 <p>Microdissection Needle Electrode, Needle Flexible, Tungsten</p> <p>100mm</p> <p>0.25x4mm</p>	

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP



Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

BOX OF 50 UNITS			
Multi-Use CDEM408	Single-Use CDES408		0.25x4mm
		Microdissection Needle Electrode, Flexible, Tungsten	
		100mm	
Multi-Use CDEM409	Single-Use CDES409		0.25x4mm
		Microdissection Needle Electrode Flexible, Tungsten	4mm Shaft
		100mm	
Multi-Use CDEM410	Single-Use CDES410		0.25x4mm
		Microdissection Needle Electrode Flexible, Tungsten	4mm Shaft
		100mm	
Multi-Use CDEM411	Single-Use CDES411		0.25x4mm
		Microdissection Needle Electrode, Flexible, 45°, Tungsten	4mm Shaft
Multi-Use CDEM412	Single-Use CDES412		0.25x4mm
		Microdissection Needle Electrode, Flexible, 45°, Tungsten	4mm Shaft
		100mm	
Multi-Use CDEM413	Single-Use CDES413		0.25x4mm
		Microdissection Needle Electrode, Flexible, 45°, Tungsten	4mm Shaft

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Electrodes Serial
C D E M / S XXX
 Cosil Multi-Use/Single-Use

<p>BOX OF 50 UNITS</p> <p>Multi-Use CDEM414 Single-Use CDES414</p>	 <p>100mm</p> <p>4mm Shaft</p> <p>0.5x4mm</p> <p>Tungsten Needle Electrode, With Ceramic, Angular</p>
<p>Multi-Use CDEM415 Single-Use CDES415</p>	 <p>100mm</p> <p>4mm Shaft</p> <p>0.5x4mm</p> <p>Tungsten Needle Electrode, With Ceramic</p>

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 2KVP

Diathermy Electrodes Serial
C DE M/S XXX
Cosil Multi-Use/Single-Use



Needle Electrode

BOX OF 50 UNITS

Multi-Use	Single-Use
CDEM416	CDES416



Multi-Use	Single-Use
CDEM417	CDES417



Multi-Use	Single-Use
CDEM418	CDES418



Multi-Use	Single-Use
CDEM419	CDES419



Multi-Use	Single-Use
CDEM420	CNES420



Multi-Use	Single-Use
CDEM421	CDES421



Multi-Use	Single-Use
CDEM422	CDES422



Multi-Use	Single-Use
CDEM423	CDES423



Multi-Use	Single-Use
CDEM424	CDES424



Needle Electrode

BOX OF 50 UNITS

Multi-Use	Single-Use
CDEM425	CDES425



Multi-Use	Single-Use
CDEM426	CDES426



Multi-Use	Single-Use
CDEM427	CDES427



Multi-Use	Single-Use
CDEM428	CNES428



Multi-Use	Single-Use
CDEM429	CDES429



Multi-Use	Single-Use
CDEM430	CDES430



Multi-Use	Single-Use
CDEM431	CDES431








Multi-Use	Single-Use
CDEM432	CDES432



Multi-Use	Single-Use
CDEM433	CDES433

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



BOX OF 50 UNITS		
Multi-Use CDEM434	Single-Use CDES434	Needle Electrode 1 90° Angular With Handle 2.4mm Shaft
		
Multi-Use CDEM435	Single-Use CDES435	Needle Electrode 90° Angular With Handle 2.4mm Shaft
		
Multi-Use CDEM436	Single-Use CNES436	Knife Electrode 45° Angled Length 115 Mm 2.4mm Shaft
		
Multi-Use CDEM413	Single-Use CDES413	Button Electrode With Handle 90° Angled Length 115 Mm 2.4mm Shaft
		
Multi-Use CDEM437	Single-Use CDES437	Ball Electrode With Handle 2.4mm Shaft

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

 UMAX: 4KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

BOX OF 50 UNITS		← 115mm →	4X GB/M
Multi-Use CDEM438	Single-Use CDES438	Needle Electrode Angular with handle	4mm Shaft
Multi-Use CDEM439		← 115mm →	4MM
Multi-Use CDEM439	Single-Use CDES439	Knife Electrode 45° Angled	4mm Shaft
Multi-Use CDEM440		← 115mm →	4MM 13MM
Multi-Use CDEM440	Single-Use CNES440	Button Electrode with handle 90° Angled	4mm Shaft
Multi-Use CDEM441		← 115mm →	3MM
Multi-Use CDEM441	Single-Use CDES441	Ball Electrode with handle	4mm Shaft
Multi-Use CDEM442		← 115mm →	15X 0.9MM
Multi-Use CDEM442	Single-Use CDES442	Needle Electrode 90° Angular, with handle	4mm Shaft


Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 2KVP

Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use

Bipolar ENT Electrodes Steel

BOX OF 50 UNITS



200mm
0.50x10mm
Straight Bipolar Electrodes Hand Pieces
Multi-Use CDEM443 Single-Use CDES443

Curved Bipolar Electrodes Hand Pieces Angle Tip

0.50x10mm
Multi-Use CDEM444 Single-Use CDES444

Short Angled Bipolar Electrodes Hand Pieces Angle Tip

0.50x10mm
Multi-Use CDEM445 Single-Use CDES445

Angled Bipolar Electrodes Hand Pieces Straight Tip


0.50x10mm
Multi-Use CDEM446 Single-Use CDES446

Long Angled Bipolar Electrodes Hand Pieces Angle Tip

0.50x10mm
Multi-Use CDEM447 Single-Use CDES447

Bipolar ENT Electrodes Teflon

BOX OF 50 UNITS



200mm
0.60x1.5mm
Arthroscopic Bipolar Cutting Electrode Hook Shaft
Multi-Use CDEM448 Single-Use CDES448

Bipolar Ent Electrode 40° Horizontal 16.5mm From Handle

0.50x15mm
Multi-Use CDEM449 Single-Use CDES449

Bipolar Ent Electrode 20° Vertical 16.5mm From Handle

0.50x20mm
Multi-Use CDEM450 Single-Use CDES450

Bipolar Ent Electrode 30° Turbiantomie 16.5mm From Handle

0.50x30mm
Multi-Use CDEM451 Single-Use CDES451

Concentric Bipolar Electrode 50° For Therapy Of Root Of Tongue

0.8x10mm
Multi-Use CDEM452 Single-Use CDES452

Concentric Bipolar Electrode 30° For Therapy Of Polyps

0.6x5mm
Multi-Use CDEM453 Single-Use CDES453

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Laproscopic Electrodes without Suction

BOX OF 50 UNITS

330mm

Lab Electrode Fully Insulated With 4 Mm Plug

Multi-Use CDEM454 Single-Use CDES454



5X330MM Ball Electrode

Multi-Use CDEM455 Single-Use CDES455



5X330MM Needle Electrode

Multi-Use CDEM456 Single-Use CNES456



5X330MM Spatula Electrode

Multi-Use CDEM457 Single-Use CDES457



5X330MM J-hook Electrode

Multi-Use CDEM458 Single-Use CDES458

Laproscopic Electrodes with Suction

BOX OF 50 UNITS

330mm

Electrodes With Ceramic, Fully Insulated With Suction, Elect. Connector

Multi-Use CDEM459 Single-Use CDES459



5X330MM Ball Electrode

Multi-Use CDEM460 Single-Use CDES460



5X330MM Needle Electrode

Multi-Use CDEM461 Single-Use CNES461



5X330MM Spatula Electrode

Multi-Use CDEM462 Single-Use CDES462



5X330MM J-hook Electrode

Multi-Use CDEM463 Single-Use CDES463

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max Vp UMAX: 2/4KVP


Diathermy Electrodes Serial
C D E M / S XXX
Cosil Multi-Use/Single-Use




ELECTROSURGICAL
**DIATHERMY
ELECTRODES PAD**

ELECTROSURGICAL
**DIATHERMY
ELECTRODES PAD**


Electrosurgical diathermy pad electrodes, or return electrodes, safely disperse electrical energy in monopolar electrosurgery, preventing burns. Made from conductive materials with adhesive hydrogel, they ensure secure skin contact. Available in single-use or reusable options, some include contact quality monitoring (CQM) for added safety. Typically placed on areas with sufficient muscle mass, they are vital for patient safety during procedures.



Beneath the electrode surface lies a conductive layer, usually composed of a flexible metal mesh or foil. This layer distributes the electrical current evenly across the pad's surface, ensuring that heat is generated uniformly. The design of the conductive layer is critical to the pad's



The electromagnetic coil is responsible for generating the oscillating high-frequency currents that produce the therapeutic heat. Located within the pad, the coil creates a controlled electromagnetic field that transfers energy deep into tissues, stimulating circulation and promoting healing.



Diathermy pads contain an insulation layer that surrounds the conductive elements. This layer is essential for patient safety, as it helps to prevent electrical leakage or accidental burns. High-quality, heat-resistant materials, such as silicone or specialized foams, are used in this layer to maximize safety and comfort.



Classic Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Adult Plate Dimension 260x150x3.3mm




Martin, ERBE, Berchtold, Tekno Series
Generators
CDP101



Valleylab, Conmed, BOWA Series
Generators
CDP102

Classic Silicon Rubber Patient Return Plate – Without Cable



Adult Plate Dimension
260x150x3.3mm
CDP103

Child Child Plate Dimension
200x80x3.3mm
CDP104

Upgraded – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Martin, ERBE, Berchtold, Tekno Series
Generators
CDP105



Valleylab, Conmed, BOWA Series
Generators
CDP106

Upgraded – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors




Martin, ERBE, Berchtold, Tekno Series
Generators
CDP107



Valleylab, Conmed, BOWA Series
Generators
CDP108

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Diathermy Pads
C DP XXX
Cosil Serial

Premium – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Bipolar Adult Plate Dimension 200x150x3.3mm




ERBE, Berchtold, Tekno Series
Generators
CDP109



Valleylab, Conmed, BOWA Series
Generators
CDP110

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



ERBE, Berchtold, Tekno Series
Generators

CDP113

Single-Use Monopolar
Patient Return Plates - Child

CDP114

Single-Use Monopolar
Patient Return Plates - Adult



Valleylab, Conmed, BOWA Series
Generators

CDP115

CDP116

Premium – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Monopolar Adult Plate Dimension 200x150x3.3mm




ERBE, Berchtold, Tekno Series
Generators
CDP111



Valleylab, Conmed, BOWA Series
Generators
CDP112

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



ERBE, Berchtold, Tekno Series
Generators

CDP117

Single-Use Monopolar
Patient Return Plates - Child

CDP118

Single-Use Monopolar
Patient Return Plates - Adult




Valleylab, Conmed, BOWA Series
Generators

CDP119


CDP120

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Premium Pad Serial
PP/WC- XXX
Without connector

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors




Single-Use Monopolar Patient Return Plates - Child

CDP121

Single-Use Monopolar Patient Return Plates - Adult

CDP122

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



Single-Use Monopolar Patient Return Plates - Child

CDP123

Single-Use Monopolar Patient Return Plates - Adult

CDP124



ERBE, Berchtold, Tekno Series Generators



Valleylab, Conmed, BOWA Series Generators

Multi-Use CDP125
Single-Use CDP126

High Frequency Jack for Split/Round Bipolar Plates

CDP127
CDP128



Multi-Use CDP129
Single-Use CDP130

High Frequency Jack for Split/Round Bipolar Plates

CDP131
CDP132




Tip Cleaner Pads – Available in 3 Different Sizes



Single-Use Tip Cleaner Pad

CDP133 45x45mm
CDP134 48x48mm
CDP135 50x50mm

Note for above products:
Special lengths+Tips and/or colours can be custom ordered as per the demand



Cable Pad Reuse/Dispose
WOC/CP/TP R/D - XXX
Without Connector: Tip Cleaner: Serial



SINGLE USE
**FORCEPS
INTEGRATED
CABLE**



SINGLE USE

**FORCEPS
INTEGRATED CABLE**

Electrosurgical forceps with an integrated cable are precise surgical tools that combine tissue manipulation with electrical energy for cutting and coagulation. Their ergonomic design, durable materials, and seamless integration enhance safety, efficiency, and ease of use across various medical procedures.



The middle area of non-stick forceps often features a textured or grooved surface to provide a secure and controlled grip for the user. This section is usually reinforced for added stability and balance, allowing for precise manipulation even during delicate procedures.



Forceps tips are precision-engineered for targeted energy delivery, featuring durable, conductive materials and insulation to ensure accuracy and minimize thermal damage.



Forceps Plug is an essential accessory for medical and laboratory equipment, designed to enhance the functionality of standard forceps by adding capabilities.



Electrosurgical forceps with an integrated cable design combine precision and functionality, eliminating the need for detachable cables. This streamlined design enhances safety, reduces clutter, and ensures reliable energy transmission for cutting and coagulation during surgical procedures.



07

SU Bipolar Forceps Integrated Cable

Regular Insulated Forceps-Blue

Regular Single Use Mcherson Bipolar Forceps with Cable-3m



	102mm	108mm	115mm	
International 4mm Banana Plug Fitting	CSUB101T	CSUB101U	CSUB101V	0.50mm
	CSUB102T	CSUB102U	CSUB102V	1.00mm
	CSUB103T	CSUB103U	CSUB103V	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB104T	CSUB104U	CSUB104V	0.50mm
	CSUB105T	CSUB105U	CSUB105V	1.00mm
	CSUB106T	CSUB106U	CSUB106V	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUB107T	CSUB107U	CSUB107V	0.50mm
	CSUB108T	CSUB108U	CSUB108V	1.00mm
	CSUB109T	CSUB109U	CSUB109V	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB110T	CSUB110U	CSUB110V	0.50mm
	CSUB111T	CSUB111U	CSUB111V	1.00mm
	CSUB112T	CSUB112U	CSUB112V	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB113T	CSUB113U	CSUB113V	0.50mm
	CSUB114T	CSUB114U	CSUB114V	1.00mm
	CSUB115T	CSUB115U	CSUB115V	2.00mm

Box of 20 Units

Regular Single Use Jeweller Bipolar Forceps with Cable-3m



	102mm	115mm	125mm	
International 4mm Banana Plug Fitting	CSUB131T	CSUB131V	CSUB131Y	0.50mm
	CSUB132T	CSUB132V	CSUB132Y	1.00mm
	CSUB133T	CSUB133V	CSUB133Y	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB134T	CSUB134V	CSUB134Y	0.50mm
	CSUB135T	CSUB135V	CSUB135Y	1.00mm
	CSUB136T	CSUB136V	CSUB136Y	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUB137T	CSUB137V	CSUB137Y	0.50mm
	CSUB138T	CSUB138V	CSUB138Y	1.00mm
	CSUB139T	CSUB139V	CSUB139Y	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB140T	CSUB140V	CSUB140Y	0.50mm
	CSUB141T	CSUB141V	CSUB141Y	1.00mm
	CSUB142T	CSUB142V	CSUB142Y	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB143T	CSUB143V	CSUB143Y	0.50mm
	CSUB144T	CSUB144V	CSUB144Y	1.00mm
	CSUB145T	CSUB145V	CSUB145Y	2.00mm

Box of 20 Units

Regular Single Use Mcherson Bipolar Forceps Angled with Cable-3m



	102mm	108mm	115mm	
International 4mm Banana Plug Fitting	CSUB116T	CSUB116U	CSUB116V	0.50mm
	CSUB117T	CSUB117U	CSUB117V	1.00mm
	CSUB118T	CSUB118U	CSUB118V	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB119T	CSUB119U	CSUB119V	0.50mm
	CSUB120T	CSUB120U	CSUB120V	1.00mm
	CSUB121T	CSUB121U	CSUB121V	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUB122T	CSUB122U	CSUB122V	0.50mm
	CSUB123T	CSUB123U	CSUB123V	1.00mm
	CSUB124T	CSUB124U	CSUB124V	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB125T	CSUB125U	CSUB125V	0.50mm
	CSUB126T	CSUB126U	CSUB126V	1.00mm
	CSUB127T	CSUB127U	CSUB127V	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB128T	CSUB128U	CSUB128V	0.50mm
	CSUB129T	CSUB129U	CSUB129V	1.00mm
	CSUB130T	CSUB130U	CSUB130V	2.00mm

Box of 20 Units

Regular Single Use Jeweller Bipolar Forceps Angled with Cable-3m



	102mm	115mm	125mm	
International 4mm Banana Plug Fitting	CSUB146T	CSUB146V	CSUB146Y	0.50mm
	CSUB147T	CSUB147V	CSUB147Y	1.00mm
	CSUB148T	CSUB148V	CSUB148Y	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB149T	CSUB149V	CSUB149Y	0.50mm
	CSUB150T	CSUB150V	CSUB150Y	1.00mm
	CSUB151T	CSUB151V	CSUB151Y	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUB152T	CSUB152V	CSUB152Y	0.50mm
	CSUB153T	CSUB153V	CSUB153Y	1.00mm
	CSUB154T	CSUB154V	CSUB154Y	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB155T	CSUB155V	CSUB155Y	0.50mm
	CSUB156T	CSUB156V	CSUB156Y	1.00mm
	CSUB157T	CSUB157V	CSUB157Y	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB158T	CSUB158V	CSUB158Y	0.50mm
	CSUB159T	CSUB159V	CSUB159Y	1.00mm
	CSUB160T	CSUB160V	CSUB160Y	2.00mm

Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available

- 0.25mm
- 0.75mm
- 1.25mm
- 2.50mm

Also Available

- 3mm Regular
- 5mm
- On Demand

Single Use Serial
C S U B XXX X
 Cosil Blue Length



07

SU Bipolar Forceps Integrated Cable

Regular Insulated Forceps-Blue

Regular Single Use Adson Bipolar Forceps with Cable-3m



	121mm	155mm	180mm	200mm	
International 4mm Banana Plug Fitting	CSUB161W	CSUB161A	CSUB161B	CSUB161C	0.50mm
	CSUB162W	CSUB162A	CSUB162B	CSUB162C	1.00mm
	CSUB163W	CSUB163A	CSUB163B	CSUB163C	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB164W	CSUB164A	CSUB164B	CSUB164C	0.50mm
	CSUB165W	CSUB165A	CSUB165B	CSUB165C	1.00mm
	CSUB166W	CSUB166A	CSUB166B	CSUB166C	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB167W	CSUB167A	CSUB167B	CSUB167C	0.50mm
	CSUB168W	CSUB168A	CSUB168B	CSUB168C	1.00mm
	CSUB169W	CSUB169A	CSUB169B	CSUB169C	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB170W	CSUB170A	CSUB170B	CSUB170C	0.50mm
	CSUB171W	CSUB171A	CSUB171B	CSUB171C	1.00mm
	CSUB172W	CSUB172A	CSUB172B	CSUB172C	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB173W	CSUB173A	CSUB173B	CSUB173C	0.50mm
	CSUB174W	CSUB174A	CSUB174B	CSUB174C	1.00mm
	CSUB175W	CSUB175A	CSUB175B	CSUB175C	2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB191X	CSUB191Z	0.50mm
	CSUB192X	CSUB192Z	1.00mm
	CSUB193X	CSUB193Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB194X	CSUB194Z	0.50mm
	CSUB195X	CSUB195Z	1.00mm
	CSUB196X	CSUB196Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB197X	CSUB197Z	0.50mm
	CSUB198X	CSUB198Z	1.00mm
	CSUB199X	CSUB199Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB200X	CSUB200Z	0.50mm
	CSUB201X	CSUB201Z	1.00mm
	CSUB202X	CSUB202Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB203X	CSUB203Z	0.50mm
	CSUB204X	CSUB204Z	1.00mm
	CSUB205X	CSUB205Z	2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB176X	CSUB176Z	0.50mm
	CSUB177X	CSUB177Z	1.00mm
	CSUB178X	CSUB178Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB179X	CSUB179Z	0.50mm
	CSUB180X	CSUB180Z	1.00mm
	CSUB181X	CSUB181Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB182X	CSUB182Z	0.50mm
	CSUB183X	CSUB183Z	1.00mm
	CSUB184X	CSUB184Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB185X	CSUB185Z	0.50mm
	CSUB186X	CSUB186Z	1.00mm
	CSUB187X	CSUB187Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB188X	CSUB188Z	0.50mm
	CSUB189X	CSUB189Z	1.00mm
	CSUB190X	CSUB190Z	2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps Angled with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB206X	CSUB206Z	0.50mm
	CSUB207X	CSUB207Z	1.00mm
	CSUB208X	CSUB208Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB209X	CSUB209Z	0.50mm
	CSUB210X	CSUB210Z	1.00mm
	CSUB211X	CSUB211Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB212X	CSUB212Z	0.50mm
	CSUB213X	CSUB213Z	1.00mm
	CSUB214X	CSUB214Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB215X	CSUB215Z	0.50mm
	CSUB216X	CSUB216Z	1.00mm
	CSUB217A	CSUB217Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB218X	CSUB218Z	0.50mm
	CSUB219X	CSUB219Z	1.00mm
	CSUB220X	CSUB220Z	2.00mm

Box of 20 Units

Note for above products: Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
Tip Sizes: 0.25mm, 0.75mm, 1.25mm, 2.50mm

Also Available
Wire Sizes: 3mm Regular, 5mm, On Demand

Single Use Serial
C S U B XXX X
Cosil Blue Length



07

SU Bipolar Forceps Integrated Cable

Regular Insulated Forceps-Blue

Regular Single Use Straight Bipolar Forceps with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB221A	CSUB221B	CSUB221C	CSUB221D	CSUB221E	0.50mm
	CSUB222A	CSUB222B	CSUB222C	CSUB222D	CSUB222E	1.00mm
	CSUB223A	CSUB223B	CSUB223C	CSUB223D	CSUB223E	2.00mm
BOWA, Limidey, EMC Valleylab, Ethicon	CSUB224A	CSUB224B	CSUB224C	CSUB224D	CSUB224E	0.50mm
	CSUB225A	CSUB225B	CSUB225C	CSUB225D	CSUB225E	1.00mm
	CSUB226A	CSUB226B	CSUB226C	CSUB226D	CSUB226E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB227A	CSUB227B	CSUB227C	CSUB227D	CSUB227E	0.50mm
	CSUB228A	CSUB228B	CSUB228C	CSUB228D	CSUB228E	1.00mm
	CSUB229A	CSUB229B	CSUB229C	CSUB229D	CSUB229E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB230A	CSUB230B	CSUB230C	CSUB230D	CSUB230E	0.50mm
	CSUB231A	CSUB231B	CSUB231C	CSUB231D	CSUB231E	1.00mm
	CSUB232A	CSUB232B	CSUB232C	CSUB232D	CSUB132E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB233A	CSUB233B	CSUB233C	CSUB233D	CSUB233E	0.50mm
	CSUB234A	CSUB234B	CSUB234C	CSUB234D	CSUB234E	1.00mm
	CSUB235A	CSUB235B	CSUB235C	CSUB235D	CSUB235E	2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps Angled with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB236A	CSUB236B	CSUB236C	CSUB236D	CSUB236E	0.50mm
	CSUB237A	CSUB237B	CSUB237C	CSUB237D	CSUB237E	1.00mm
	CSUB238A	CSUB238B	CSUB238C	CSUB238D	CSUB238E	2.00mm
BOWA, Limidey, EMC Valleylab, Ethicon	CSUB239A	CSUB239B	CSUB239C	CSUB239D	CSUB239E	0.50mm
	CSUB240A	CSUB240B	CSUB240C	CSUB240D	CSUB240E	1.00mm
	CSUB241A	CSUB241B	CSUB241C	CSUB241D	CSUB241E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB242A	CSUB242B	CSUB242C	CSUB242D	CSUB242E	0.50mm
	CSUB243A	CSUB243B	CSUB243C	CSUB243D	CSUB243E	1.00mm
	CSUB244A	CSUB244B	CSUB244C	CSUB244D	CSUB244E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB245A	CSUB245B	CSUB245C	CSUB245D	CSUB245E	0.50mm
	CSUB246A	CSUB246B	CSUB246C	CSUB246D	CSUB246E	1.00mm
	CSUB247A	CSUB247B	CSUB247C	CSUB247D	CSUB247E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB248A	CSUB248B	CSUB248C	CSUB248D	CSUB248E	0.50mm
	CSUB249A	CSUB249B	CSUB249C	CSUB249D	CSUB249E	1.00mm
	CSUB250A	CSUB250B	CSUB250C	CSUB250D	CSUB250E	2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB251A	CSUB251B	CSUB251C	CSUB251D	CSUB251E	0.50mm
	CSUB252A	CSUB252B	CSUB252C	CSUB252D	CSUB252E	1.00mm
	CSUB253A	CSUB253B	CSUB253C	CSUB253D	CSUB253E	2.00mm
BOWA, Limidey, EMC Valleylab, Ethicon	CSUB254A	CSUB254B	CSUB254C	CSUB254D	CSUB254E	0.50mm
	CSUB255A	CSUB255B	CSUB255C	CSUB255D	CSUB255E	1.00mm
	CSUB256A	CSUB256B	CSUB256C	CSUB256D	CSUB256E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB257A	CSUB257B	CSUB257C	CSUB257D	CSUB257E	0.50mm
	CSUB258A	CSUB258B	CSUB258C	CSUB258D	CSUB258E	1.00mm
	CSUB259A	CSUB259B	CSUB259C	CSUB259D	CSUB259E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB260A	CSUB260B	CSUB260C	CSUB260D	CSUB260E	0.50mm
	CSUB261A	CSUB261B	CSUB261C	CSUB261D	CSUB261E	1.00mm
	CSUB262A	CSUB262B	CSUB262C	CSUB262D	CSUB132E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB263A	CSUB263B	CSUB263C	CSUB263D	CSUB263E	0.50mm
	CSUB264A	CSUB264B	CSUB264C	CSUB264D	CSUB264E	1.00mm
	CSUB265A	CSUB265B	CSUB265C	CSUB265D	CSUB265E	2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps Angled with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB266A	CSUB266B	CSUB266C	CSUB266D	CSUB266E	0.50mm
	CSUB267A	CSUB267B	CSUB267C	CSUB267D	CSUB267E	1.00mm
	CSUB268A	CSUB268B	CSUB268C	CSUB268D	CSUB268E	2.00mm
BOWA, Limidey, EMC Valleylab, Ethicon	CSUB269A	CSUB269B	CSUB269C	CSUB269D	CSUB269E	0.50mm
	CSUB270A	CSUB270B	CSUB270C	CSUB270D	CSUB270E	1.00mm
	CSUB271A	CSUB271B	CSUB271C	CSUB271D	CSUB271E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUB272A	CSUB272B	CSUB272C	CSUB272D	CSUB272E	0.50mm
	CSUB273A	CSUB273B	CSUB273C	CSUB273D	CSUB273E	1.00mm
	CSUB274A	CSUB274B	CSUB274C	CSUB274D	CSUB274E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB275A	CSUB275B	CSUB275C	CSUB275D	CSUB275E	0.50mm
	CSUB276A	CSUB276B	CSUB276C	CSUB276D	CSUB276E	1.00mm
	CSUB277A	CSUB277B	CSUB277C	CSUB277D	CSUB277E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB278A	CSUB278B	CSUB278C	CSUB278D	CSUB278E	0.50mm
	CSUB279A	CSUB279B	CSUB279C	CSUB279D	CSUB279E	1.00mm
	CSUB280A	CSUB280B	CSUB280C	CSUB280D	CSUB280E	2.00mm


Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand



UMAX: 0.5KVP



Single Use Serial
C S U B XXX X
Cosil Blue Length





07

SU Bipolar Forceps Integrated Cable

Non-Stick Insulated Forceps-Pink

Non-Stick Single Use Mcepherson Bipolar Forceps with Cable-3m



	102mm	108mm	115mm	
International 4mm Banana Plug Fitting	CSUP281T	CSUP281U	CSUP281V	0.50mm
	CSUP282T	CSUP282U	CSUP282V	1.00mm
	CSUP283T	CSUP283U	CSUP283V	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP284T	CSUP284U	CSUP284V	0.50mm
	CSUP285T	CSUP285U	CSUP285V	1.00mm
	CSUP286T	CSUP286U	CSUP286V	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUP287T	CSUP287U	CSUP287V	0.50mm
	CSUP288T	CSUP288U	CSUP288V	1.00mm
	CSUP289T	CSUP289U	CSUP289V	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP290T	CSUP290U	CSUP290V	0.50mm
	CSUP291T	CSUP291U	CSUP291V	1.00mm
	CSUP292T	CSUP292U	CSUP292V	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP293T	CSUP293U	CSUP293V	0.50mm
	CSUP294T	CSUP294U	CSUP294V	1.00mm
	CSUP295T	CSUP295U	CSUP295V	2.00mm

Box of 20 Units

Non-Stick Single Use Jeweller Bipolar Forceps with Cable-3m



	102mm	115mm	125mm	
International 4mm Banana Plug Fitting	CSUB311T	CSUB311V	CSUB311Y	0.50mm
	CSUB312T	CSUB312V	CSUB312Y	1.00mm
	CSUB313T	CSUB313V	CSUB313Y	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB314T	CSUB314V	CSUB314Y	0.50mm
	CSUB315T	CSUB315V	CSUB315Y	1.00mm
	CSUB316T	CSUB316V	CSUB316Y	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUB317T	CSUB317V	CSUB317Y	0.50mm
	CSUB318T	CSUB318V	CSUB318Y	1.00mm
	CSUB319T	CSUB319V	CSUB319Y	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB320T	CSUB320V	CSUB320Y	0.50mm
	CSUB321T	CSUB321V	CSUB321Y	1.00mm
	CSUB322T	CSUB322V	CSUB322Y	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB323T	CSUB323V	CSUB323Y	0.50mm
	CSUB324T	CSUB324V	CSUB324Y	1.00mm
	CSUB325T	CSUB325V	CSUB325Y	2.00mm

Box of 20 Units

Non-Stick Single Use Mcepherson Bipolar Forceps Angled with Cable-3m



	102mm	108mm	115mm	
International 4mm Banana Plug Fitting	CSUP296T	CSUP296U	CSUP296V	0.50mm
	CSUP297T	CSUP297U	CSUP297V	1.00mm
	CSUP298T	CSUP298U	CSUP298V	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP299T	CSUP299U	CSUP299V	0.50mm
	CSUP300T	CSUP300U	CSUP300V	1.00mm
	CSUP301T	CSUP301U	CSUP301V	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUP302T	CSUP302U	CSUP302V	0.50mm
	CSUP303T	CSUP303U	CSUP303V	1.00mm
	CSUP304T	CSUP304U	CSUP304V	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP305T	CSUP305U	CSUP305V	0.50mm
	CSUP306T	CSUP306U	CSUP306V	1.00mm
	CSUP307T	CSUP307U	CSUP307V	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP308T	CSUP308U	CSUP308V	0.50mm
	CSUP309T	CSUP309U	CSUP309V	1.00mm
	CSUP310T	CSUP310U	CSUP310V	2.00mm

Box of 20 Units

Non-Stick Single Use Jeweller Bipolar Forceps Angled with Cable-3m




	102mm	115mm	125mm	
International 4mm Banana Plug Fitting	CSUP326T	CSUP326V	CSUP326Y	0.50mm
	CSUP327T	CSUP327V	CSUP327Y	1.00mm
	CSUP328T	CSUP328V	CSUP328Y	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP329T	CSUP329V	CSUP329Y	0.50mm
	CSUP330T	CSUP330V	CSUP330Y	1.00mm
	CSUP331T	CSUP331V	CSUP331Y	2.00mm
ERBE 4mm, Siemens Shutter BMSQ	CSUP332T	CSUP332V	CSUP332Y	0.50mm
	CSUP333T	CSUP333V	CSUP333Y	1.00mm
	CSUP334T	CSUP334V	CSUP334Y	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP335T	CSUP335V	CSUP335Y	0.50mm
	CSUP336T	CSUP336V	CSUP336Y	1.00mm
	CSUP337T	CSUP337V	CSUP337Y	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP338T	CSUP338V	CSUP338Y	0.50mm
	CSUP339T	CSUP339V	CSUP339Y	1.00mm
	CSUP340T	CSUP340V	CSUP340Y	2.00mm

Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP



Single Use Serial
C S U P XXX X
Cosil Pink Length



07

SU Bipolar Forceps Integrated Cable

Non-Stick Insulated Forceps-Pink

Non-Stick Single Use Adson Bipolar Forceps with Cable-3m



	121mm	155mm	180mm	200mm	
International 4mm Banana Plug Fitting	CUSP341W	CUSP341A	CUSP341B	CUSP341C	0.50mm
	CUSP342W	CUSP342A	CUSP342B	CUSP342C	1.00mm
	CUSP343W	CUSP343A	CUSP343B	CUSP343C	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CUSP344W	CUSP344A	CUSP344B	CUSP344C	0.50mm
	CUSP345W	CUSP345A	CUSP345B	CUSP345C	1.00mm
	CUSP346W	CUSP346A	CUSP346B	CUSP346C	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CUSP347W	CUSP347A	CUSP347B	CUSP347C	0.50mm
	CUSP348W	CUSP348A	CUSP348B	CUSP348C	1.00mm
	CUSP349W	CUSP349A	CUSP349B	CUSP349C	2.00mm
Aesculap, Berchtold, Wolf, Martin	CUSP350W	CUSP350A	CUSP350B	CUSP350C	0.50mm
	CUSP351W	CUSP351A	CUSP351B	CUSP351C	1.00mm
	CUSP352W	CUSP352A	CUSP352B	CUSP352C	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CUSP353W	CUSP353A	CUSP353B	CUSP353C	0.50mm
	CUSP354W	CUSP354A	CUSP354B	CUSP354C	1.00mm
	CUSP355W	CUSP355A	CUSP355B	CUSP355C	2.00mm

Box of 20 Units

Non-Stick Single Use Bayonet Bipolar Forceps with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUP371X	CSUP371Z	0.50mm
	CSUP372X	CSUP372Z	1.00mm
	CSUP373X	CSUP373Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP374X	CSUP374Z	0.50mm
	CSUP375X	CSUP375Z	1.00mm
	CSUP376X	CSUP376Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP377X	CSUP377Z	0.50mm
	CSUP378X	CSUP378Z	1.00mm
	CSUP379X	CSUP379Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP380X	CSUP380Z	0.50mm
	CSUP381X	CSUP381Z	1.00mm
	CSUP382X	CSUP382Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP383X	CSUP383Z	0.50mm
	CSUP384X	CSUP384Z	1.00mm
	CSUP385X	CSUP385Z	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUP356X	CSUP356Z	0.50mm
	CSUP357X	CSUP357Z	1.00mm
	CSUP358X	CSUP358Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP359X	CSUP359Z	0.50mm
	CSUP360X	CSUP360Z	1.00mm
	CSUP361X	CSUP361Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP362X	CSUP362Z	0.50mm
	CSUP363X	CSUP363Z	1.00mm
	CSUP364X	CSUP364Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP365X	CSUP365Z	0.50mm
	CSUP366X	CSUP366Z	1.00mm
	CSUP367X	CSUP367Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP368X	CSUP368Z	0.50mm
	CSUP369X	CSUP369Z	1.00mm
	CSUP370X	CSUP370Z	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps Angled with Cable-3m



	140mm	165mm	
International 4mm Banana Plug Fitting	CSUP386X	CSUP386Z	0.50mm
	CSUP387X	CSUP387Z	1.00mm
	CSUP388X	CSUP388Z	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP389X	CSUP389Z	0.50mm
	CSUP390X	CSUP390Z	1.00mm
	CSUP391X	CSUP391Z	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP392X	CSUP392Z	0.50mm
	CSUP393X	CSUP393Z	1.00mm
	CSUP394X	CSUP394Z	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP395X	CSUP395Z	0.50mm
	CSUP396X	CSUP396Z	1.00mm
	CSUP397A	CSUP397Z	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP398X	CSUP398Z	0.50mm
	CSUP399X	CSUP399Z	1.00mm
	CSUP400X	CSUP400Z	2.00mm

Box of 20 Units

Note for above products: Special lengths+Tips and/or colours can be custom ordered as per the demand

95 137 max vp UMAX: 0.5KVP

Also Available
Tip Sizes: 0.25mm, 0.75mm, 1.25mm, 2.50mm

Also Available
Wire Sizes: 3mm, 5mm, On Demand

Single Use Serial
C S U P XXX X
Cosil Pink Length



07

SU Bipolar Forceps Integrated Cable

Non-Stick Insulated Forceps-Pink

Non-Stick Single Use Straight Bipolar Forceps with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP401A	CSUP401B	CSUP401C	CSUP401D	CSUP401E	0.50mm
	CSUP402A	CSUP402B	CSUP402C	CSUP402D	CSUP402E	1.00mm
	CSUP403A	CSUP403B	CSUP403C	CSUP403D	CSUP403E	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP404A	CSUP404B	CSUP404C	CSUP404D	CSUP404E	0.50mm
	CSUP405A	CSUP405B	CSUP405C	CSUP405D	CSUP405E	1.00mm
	CSUP406A	CSUP406B	CSUP406C	CSUP406D	CSUP406E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP407A	CSUP407B	CSUP407C	CSUP407D	CSUP407E	0.50mm
	CSUP408A	CSUP408B	CSUP408C	CSUP408D	CSUP408E	1.00mm
	CSUP409A	CSUP409B	CSUP409C	CSUP409D	CSUP409E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP410A	CSUP410B	CSUP410C	CSUP410D	CSUP410E	0.50mm
	CSUP411A	CSUP411B	CSUP411C	CSUP411D	CSUP411E	1.00mm
	CSUP412A	CSUP412B	CSUP412C	CSUP412D	CSUP412E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP413A	CSUP413B	CSUP413C	CSUP413D	CSUP413E	0.50mm
	CSUP414A	CSUP414B	CSUP414C	CSUP414D	CSUP414E	1.00mm
	CSUP415A	CSUP415B	CSUP415C	CSUP415D	CSUP415E	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps Angled with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP416A	CSUP416B	CSUP416C	CSUP416D	CSUP416E	0.50mm
	CSUP417A	CSUP417B	CSUP417C	CSUP417D	CSUP417E	1.00mm
	CSUP418A	CSUP418B	CSUP418C	CSUP418D	CSUP418E	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP419A	CSUP419B	CSUP419C	CSUP419D	CSUP419E	0.50mm
	CSUP420A	CSUP420B	CSUP420C	CSUP420D	CSUP420E	1.00mm
	CSUP421A	CSUP421B	CSUP421C	CSUP421D	CSUP421E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP422A	CSUP422B	CSUP422C	CSUP422D	CSUP422E	0.50mm
	CSUP423A	CSUP423B	CSUP423C	CSUP423D	CSUP423E	1.00mm
	CSUP424A	CSUP424B	CSUP424C	CSUP424D	CSUP424E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP425A	CSUP425B	CSUP425C	CSUP425D	CSUP425E	0.50mm
	CSUP426A	CSUP426B	CSUP426C	CSUP426D	CSUP426E	1.00mm
	CSUP427A	CSUP427B	CSUP427C	CSUP427D	CSUP427E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP428A	CSUP428B	CSUP428C	CSUP428D	CSUP428E	0.50mm
	CSUP429A	CSUP429B	CSUP429C	CSUP429D	CSUP429E	1.00mm
	CSUP430A	CSUP430B	CSUP430C	CSUP430D	CSUP430E	2.00mm

Box of 20 Units

Non-Stick Single Use Bayonet Bipolar Forceps with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP431A	CSUP431B	CSUP431C	CSUP431D	CSUP431E	0.50mm
	CSUP432A	CSUP432B	CSUP432C	CSUP432D	CSUP432E	1.00mm
	CSUP433A	CSUP433B	CSUP433C	CSUP433D	CSUP433E	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP434A	CSUP434B	CSUP434C	CSUP434D	CSUP434E	0.50mm
	CSUP435A	CSUP435B	CSUP435C	CSUP435D	CSUP435E	1.00mm
	CSUP436A	CSUP436B	CSUP436C	CSUP436D	CSUP436E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP437A	CSUP437B	CSUP437C	CSUP437D	CSUP437E	0.50mm
	CSUP438A	CSUP438B	CSUP438C	CSUP438D	CSUP438E	1.00mm
	CSUP439A	CSUP439B	CSUP439C	CSUP439D	CSUP439E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP440A	CSUP440B	CSUP440C	CSUP440D	CSUP440E	0.50mm
	CSUP441A	CSUP441B	CSUP441C	CSUP441D	CSUP441E	1.00mm
	CSUP442A	CSUP442B	CSUP442C	CSUP442D	CSUP442E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP443A	CSUP443B	CSUP443C	CSUP443D	CSUP443E	0.50mm
	CSUP444A	CSUP444B	CSUP444C	CSUP444D	CSUP444E	1.00mm
	CSUP445A	CSUP445B	CSUP445C	CSUP445D	CSUP445E	2.00mm

Box of 20 Units

Non-Stick Single Use Bayonet Bipolar Forceps Angled with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP446A	CSUP446B	CSUP446C	CSUP446D	CSUP446E	0.50mm
	CSUP447A	CSUP447B	CSUP447C	CSUP447D	CSUP447E	1.00mm
	CSUP448A	CSUP448B	CSUP448C	CSUP448D	CSUP448E	2.00mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUP449A	CSUP449B	CSUP449C	CSUP449D	CSUP449E	0.50mm
	CSUP450A	CSUP450B	CSUP450C	CSUP450D	CSUP450E	1.00mm
	CSUP451A	CSUP451B	CSUP451C	CSUP451D	CSUP451E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP452A	CSUP452B	CSUP452C	CSUP452D	CSUP452E	0.50mm
	CSUP453A	CSUP453B	CSUP453C	CSUP453D	CSUP453E	1.00mm
	CSUP454A	CSUP454B	CSUP454C	CSUP454D	CSUP454E	2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUP455A	CSUP455B	CSUP455C	CSUP455D	CSUP455E	0.50mm
	CSUP456A	CSUP456B	CSUP456C	CSUP456D	CSUP456E	1.00mm
	CSUP457A	CSUP457B	CSUP457C	CSUP457D	CSUP457E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP458A	CSUP458B	CSUP458C	CSUP458D	CSUP458E	0.50mm
	CSUP459A	CSUP459B	CSUP459C	CSUP459D	CSUP459E	1.00mm
	CSUP460A	CSUP460B	CSUP460C	CSUP460D	CSUP460E	2.00mm


Box of 20 Units

Note for above products:

Special lengths+Tips and/or colours can be custom ordered as per the demand



UMAX: 0.5KVP



Single Use Serial

C S U P X X X X

Cosil Pink Length






ELECTROSURGICAL
GENERATORS



ELECTROSURGICAL
GENERATORS


Our Diathermy Electrocautery Machines are designed to provide reliable and precise energy delivery for a wide range of surgical procedures. Equipped with advanced technology, these machines offer adjustable power settings, ensuring optimal performance for both cutting and coagulation. With user-friendly interfaces and safety features, they deliver consistent results while maintaining patient safety. Ideal for various medical specialties, our machines are engineered for durability, ease of use, and enhanced surgical outcomes.



Convenient Operation
by Finger-suit Buttons



Simple Error Code for
Fast Trouble Shooting



REM (Return Electrode Monitoring)
System to avoid any burnt or damage
to the surgeon and patient.



It can be connected to
endoscopes for more
clinical applications

CEG101

Smart Electrosurgical Unit

Power Supply:

110V, 50/60Hz or 220V, 50/60Hz

Dimension:

Product: 330(L) x 280(W) x 100(H) mm

Package: 440(L) x 390(W) x 260(H) mm

Electrosurgical Unit

150 WATTS

Weight:

Net Weight: 4.7KG

Gross Weight: 8.9KG

Compact design for both monopolar & bipolar.
 Max 4 working modes: Pure/Blend, Standard/Bipolar.
 Fast trouble-shooting by default Error alarm system
 REM (Return Electrode Monitoring) system for highest safety; Excellent operation performance with minimal sticking & carbonization. Different & adjustable operation sounds in case of misoperation; Wide clinical application, Suitable for different surgery



Easy-to-Use

Convenient Operation by Finger-suit Buttons



Error Code

Simple Error Code for Fast Trouble-shooting



High Safety

Rem (Return Electrode Monitoring) system to avoid any burnt or damage to the surgeon and patient.



Wide Application


It can be connected to endoscopes for more clinical applications.

Standard Accessory



Optional

Note for above products:
 Special lengths+Tips and/or colours can be custom ordered as per the demand



Electrosurgical Generator
C E G XXX
 Cosil Serial

CEG102

Power Supply:

110V, 50/60Hz or 220V, 50/60Hz

Dimension:

Product: 421(L) x 357(W) x 162(H) mm
 Package: 520(L) x 445(W) x 305(H) mm


Electrosurgical Unit 400WATTS

Weight:

Net Weight: 8.7KG
 Gross Weight: 14.5KG

High-frequency Electrosurgical Unit

Easy-to-use buttons for simple operation
 High power output 400W for more operation selection
 Up to 7 work modes: Pure/Blend1/Blend2, Standard/Spray, Bipo1/Bipo2
 Increment of 5W micro power control system to ensure optimal effect
 Minimal sticking and carbonization, reduce the surgery damage
 Fast trouble-shooting by default error alarm system
 Auto-terminate of power output after 5 seconds to avoid misoperation
 Support endoscopes connection and underwater cutting (TURP etc) OEM projects



Multiple Functions
Clear Color Indicators for Work



Easy-to-Use
Convenient Operation by Finger-suit Buttons



Error Code
Simple Error Code for Fast Trouble-shooting

Mode & Alarm



High Safety
Rem(Return Electrode Monitoring) system to avoid any burnt or damage to the surgeon and patient.




Wide Application
It can be connected to endoscopes for more clinical applications.


Standard Accessory



Optional



Note for above products:
 Special lengths+Tips and/or colours can be custom ordered as per the demand



Electrosurgical Generator
CEG XXX
 Cosil Serial



We Provide
24/7 Customer Support
cosilinstruments@gmail.com

sales@cosilinstruments.com

www.cosilinstruments.com