



PRODUCT
CATALOGUE
UPDATED EDITION

Instrumente für die Elektrochirurgie
Instruments for electrosurgery
Instruments pour électrochirurgie
Instrumentos para la electrocirugía
Strumenti per elettrochirurgia



VISION

To be a globally trusted manufacturer of electrosurgical instruments and surgical solutions, recognized for innovation, precision, and uncompromising quality. We strive to empower healthcare professionals worldwide with reliable, high-performance products that enhance surgical outcomes and patient safety. Through continuous innovation, advanced manufacturing, and strict adherence to international quality standards, we aim to build lasting partnerships with customers and establish ourselves as a leading name in the global medical device industry.

MISSION

Our mission is to design, develop, and manufacture high-quality electrosurgical instruments that deliver exceptional performance, precision, and patient safety. We are committed to meeting the evolving needs of healthcare professionals through advanced manufacturing technologies, continuous innovation, and rigorous quality control. By maintaining international standards, providing responsive customer support, and building long-term partnerships, we strive to deliver reliable surgical solutions that improve clinical outcomes and create lasting value for customers worldwide.

OUR QUALITY APPROACH

Quality is the foundation of everything we do. Every electrosurgical instrument is manufactured under strict quality-controlled processes using premium materials, advanced manufacturing technologies, and rigorous inspection procedures. Our commitment to continuous improvement, regulatory compliance, and customer satisfaction ensures that every product delivers exceptional precision, reliability, safety, and consistent performance while meeting international quality standards.



CONTENT

cosil instrument company

04

ELECTROSURGICAL
FORCEPS

20

C-TECH
INSTRUMENTS

36

DIATHERMY
CABLES AND PLUGS

40

GYNAECOLOGY
INSTRUMENTS

44

CAUTERY
AND ELECTRODES

76

PAD
ELECTRODE

80

SINGLE USE
BIPOLAR FORCEPS

90

ELECTROSURGICAL
GENERATORS



01

ELECTROSURGICAL INSULATED FORCEPS

PRECISE. INSULATED. RELIABLE.

INSULATED BIPOLAR FORCEPS

FLAT PLUG CONNECTION



UMAX: 0.5KVP



Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM101A	CEFM101B	CEFM101C	CEFM101D	CEFM101E	0.50mm
	CEFM102A	CEFM102B	CEFM102C	CEFM102D	CEFM102E	1.00mm
	CEFM103A	CEFM103B	CEFM103C	CEFM103D	CEFM103E	2.00mm
Single-Use	CEFS101A	CEFS101B	CEFS101C	CEFS101D	CEFS101E	0.50mm
	CEFS102A	CEFS102B	CEFS102C	CEFS102D	CEFS102E	1.00mm
	CEFS103A	CEFS103B	CEFS103C	CEFS103D	CEFS103E	2.00mm

Box of 20 Units



Insulated Bayonet Tip Down Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM110A	CEFM110B	CEFM110C	CEFM110D	CEFM110E	0.50mm
	CEFM111A	CEFM111B	CEFM111C	CEFM111D	CEFM111E	1.00mm
	CEFM112A	CEFM112B	CEFM112C	CEFM112D	CEFM112E	2.00mm
Single-Use	CEFS110A	CEFS110B	CEFS110C	CEFS110D	CEFS110E	0.50mm
	CEFS111A	CEFS111B	CEFS111C	CEFS111D	CEFS111E	1.00mm
	CEFS112A	CEFS112B	CEFS112C	CEFS112D	CEFS112E	2.00mm

Box of 20 Units



Insulated Straight Angle Tip Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM104A	CEFM104B	CEFM104C	CEFM104D	CEFM104E	0.50mm
	CEFM105A	CEFM105B	CEFM105C	CEFM105D	CEFM105E	1.00mm
	CEFM106A	CEFM106B	CEFM106C	CEFM106D	CEFM106E	2.00mm
Single-Use	CEFS104A	CEFS104B	CEFS104C	CEFS104D	CEFS104E	0.50mm
	CEFS105A	CEFS105B	CEFS105C	CEFS105D	CEFS105E	1.00mm
	CEFS106A	CEFS106B	CEFS106C	CEFS106D	CEFS106E	2.00mm

Box of 20 Units



Insulated Bayonet Tip Up Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM113A	CEFM113B	CEFM113C	CEFM113D	CEFM113E	0.50mm
	CEFM114A	CEFM114B	CEFM114C	CEFM114D	CEFM114E	1.00mm
	CEFM115A	CEFM115B	CEFM115C	CEFM115D	CEFM115E	2.00mm
Single-Use	CEFS113A	CEFS113B	CEFS113C	CEFS113D	CEFS113E	0.50mm
	CEFS114A	CEFS114B	CEFS114C	CEFS114D	CEFS114E	1.00mm
	CEFS115A	CEFS115B	CEFS115C	CEFS115D	CEFS115E	2.00mm

Box of 20 Units



Insulated Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM107A	CEFM107B	CEFM107C	CEFM107D	CEFM107E	0.50mm
	CEFM108A	CEFM108B	CEFM108C	CEFM108D	CEFM108E	1.00mm
	CEFM109A	CEFM109B	CEFM109C	CEFM109D	CEFM109E	2.00mm
Single-Use	CEFS107A	CEFS107B	CEFS107C	CEFS107D	CEFS107E	0.50mm
	CEFS108A	CEFS108B	CEFS108C	CEFS108D	CEFS108E	1.00mm
	CEFS109A	CEFS109B	CEFS109C	CEFS109D	CEFS109E	2.00mm

Box of 20 Units



01 ELECTROSURGICAL FORCEPS

INSULATED BIPOLAR FORCEPS

FLAT PLUG CONNECTION



UMAX: 0.5KVP



Mchpherson Bipolar Forceps

Multi-Use	102mm	
	CEFM116T	0.50mm
	CEFM117T	1.00mm
Single-Use	CEFM118T	2.00mm
	CEFS116T	0.50mm
	CEFS117T	1.00mm
	CEFS118T	2.00mm

Box of 20 Units



Jeweler#7 Bipolar Forceps

Multi-Use	115mm	
	CEFM119U	0.50mm
	CEFM120U	1.00mm
Single-Use	CEFM121U	2.00mm
	CEFS119U	0.50mm
	CEFS120U	1.00mm
	CEFS121U	2.00mm

Box of 20 Units



Jeweler#4 Ang. Tip Bipolar Forceps

Multi-Use	108mm	
	CEFM122U	0.50mm
	CEFM123U	1.00mm
Single-Use	CEFM124U	2.00mm
	CEFS122U	0.50mm
	CEFS123U	1.00mm
	CEFS124U	2.00mm

Box of 20 Units



Jeweler#4 Str. Tip Bipolar Forceps

Multi-Use	115mm	
	CEFM125V	0.50mm
	CEFM126V	1.00mm
Single-Use	CEFM127V	2.00mm
	CEFS125V	0.50mm
	CEFS126V	1.00mm
	CEFS127V	2.00mm

Box of 20 Units



Adson Bipolar Forceps

Multi-Use	121mm	155mm	180mm	200mm	
	CEFM128W	CEFM128A	CEFM128B	CEFM128C	0.50mm
	CEFM129W	CEFM129A	CEFM129B	CEFM129C	1.00mm
Single-Use	CEFM130W	CEFM130A	CEFM130B	CEFM130C	2.00mm
	CEFS128W	CEFS128A	CEFS128B	CEFS128C	0.50mm
	CEFS129W	CEFS129A	CEFS129B	CEFS129C	1.00mm
	CEFS130W	CEFS130A	CEFS130B	CEFS130C	2.00mm

05



INSULATED BIPOLAR FORCEPS

TWO PIN CONNECTION



UMAX: 0.5KVP



Iris Str. Bipolar Forceps

108mm	
Multi-Use	CEFM167U 0.50mm
	CEFM168U 1.00mm
	CEFM169U 2.00mm
Box of 20 Units	
Single-Use	CEFS167U 0.50mm
	CEFS168U 1.00mm
	CEFS169U 2.00mm



Iris Ang. Bipolar Forceps

108mm	
Multi-Use	CEFM170U 0.50mm
	CEFM171U 1.00mm
	CEFM172U 2.00mm
Box of 20 Units	
Single-Use	CEFS170U 0.50mm
	CEFS171U 1.00mm
	CEFS172U 2.00mm



Jeweler#5 Str. Bipolar Forceps

115mm	
Multi-Use	CEFM173V 0.50mm
	CEFM174V 1.00mm
	CEFM175V 2.00mm
Box of 20 Units	
Single-Use	CEFS173V 0.50mm
	CEFS174V 1.00mm
	CEFS175V 2.00mm



Jeweler#4 Str. Bipolar Forceps

115mm	
Multi-Use	CEFM176V 0.50mm
	CEFM177V 1.00mm
	CEFM178V 2.00mm
Box of 20 Units	
Single-Use	CEFS176V 0.50mm
	CEFS177V 1.00mm
	CEFS178V 2.00mm



Jeweler#4 Pointed Bipolar Forceps

115mm	
Multi-Use	CEFM179V 0.50mm
	CEFM180V 1.00mm
	CEFM181V 2.00mm
Box of 20 Units	
Single-Use	CEFS179V 0.50mm
	CEFS180V 1.00mm
	CEFS181V 2.00mm



Jeweler#7 Bipolar Forceps

115mm	
Multi-Use	CEFM182V 0.50mm
	CEFM183V 1.00mm
	CEFM184V 2.00mm
Box of 20 Units	
Single-Use	CEFS182V 0.50mm
	CEFS183V 1.00mm
	CEFS184V 2.00mm



INSULATED BIPOLAR FORCEPS

TWO PIN CONNECTION



UMAX: 0.5KVP



Coaptation Special Bipolar Forceps

115mm	
Multi-Use	CEFM185V 0.50mm
	CEFM186V 1.00mm
	CEFM187V 2.00mm
Box of 20 Units	
Single-Use	CEFS185V 0.50mm
	CEFS186V 1.00mm
	CEFS187V 2.00mm



Jeweler#9 Bipolar Forceps

108mm	
Multi-Use	CEFM188U 0.50mm
	CEFM189U 1.00mm
	CEFM190U 2.00mm
Box of 20 Units	
Single-Use	CEFS188U 0.50mm
	CEFS189U 1.00mm
	CEFS190U 2.00mm



Jeweler#4 Monopolar Forceps

115mm	
Multi-Use	CEFM191V 0.50mm
	CEFM192V 1.00mm
	CEFM193V 2.00mm
Box of 20 Units	
Single-Use	CEFS191V 0.50mm
	CEFS192V 1.00mm
	CEFS193V 2.00mm



Jeweler#7 Monopolar Forceps

108mm	
Multi-Use	CEFM194U 0.50mm
	CEFM195U 1.00mm
	CEFM196U 2.00mm
Box of 20 Units	
Single-Use	CEFS194U 0.50mm
	CEFS195U 1.00mm
	CEFS196U 2.00mm



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM197W	CEFM197A	CEFM197B	CEFM197C	 0.50mm
	CEFM198W	CEFM198A	CEFM198B	CEFM198C	 1.00mm
	CEFM199W	CEFM199A	CEFM199B	CEFM199C	 2.00mm
Box of 20 Units					
Single-Use	CEFS197W	CEFS197A	CEFS197B	CEFS197C	 0.50mm
	CEFS198W	CEFS198A	CEFS198B	CEFS198C	 1.00mm
	CEFS199W	CEFS199A	CEFS199B	CEFS199C	 2.00mm



INSULATED BIPOLAR FORCEPS

TWO PIN CONNECTION



UMAX: 0.5KVP



Insulated Malis Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM131A	CEFM131B	CEFM131C	CEFM131D	CEFM131E	0.50mm	
	CEFM132A	CEFM132B	CEFM132C	CEFM132D	CEFM132E	1.00mm	
	CEFM133A	CEFM133B	CEFM133C	CEFM133D	CEFM133E	2.00mm	
Single-Use	CEFS131A	CEFS131B	CEFS131C	CEFS131D	CEFS131E	0.50mm	
	CEFS132A	CEFS132B	CEFS132C	CEFS132D	CEFS132E	1.00mm	
	CEFS133A	CEFS133B	CEFS133C	CEFS133D	CEFS133E	2.00mm	

Box of 20 Units



Insulated Hardy Angled Shape Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM140A	CEFM140B	CEFM140C	CEFM140D	CEFM140E	0.50mm	
	CEFM141A	CEFM141B	CEFM141C	CEFM141D	CEFM141E	1.00mm	
	CEFM142A	CEFM142B	CEFM142C	CEFM142D	CEFM142E	2.00mm	
Single-Use	CEFS140A	CEFS140B	CEFS140C	CEFS140D	CEFS140E	0.50mm	
	CEFS141A	CEFS141B	CEFS141C	CEFS141D	CEFS141E	1.00mm	
	CEFS142A	CEFS142B	CEFS142C	CEFS142D	CEFS142E	2.00mm	

Box of 20 Units



Insulated Jensen Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM134A	CEFM134B	CEFM134C	CEFM134D	CEFM134E	0.50mm	
	CEFM135A	CEFM135B	CEFM135C	CEFM135D	CEFM135E	1.00mm	
	CEFM136A	CEFM136B	CEFM136C	CEFM136D	CEFM136E	2.00mm	
Single-Use	CEFS134A	CEFS134B	CEFS134C	CEFS134D	CEFS134E	0.50mm	
	CEFS135A	CEFS135B	CEFS135C	CEFS135D	CEFS135E	1.00mm	
	CEFS136A	CEFS136B	CEFS136C	CEFS136D	CEFS136E	2.00mm	

Box of 20 Units



Insulated Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM143A	CEFM143B	CEFM143C	CEFM143D	CEFM143E	0.50mm	
	CEFM144A	CEFM144B	CEFM144C	CEFM144D	CEFM144E	1.00mm	
	CEFM145A	CEFM145B	CEFM145C	CEFM145D	CEFM145E	2.00mm	
Single-Use	CEFS143A	CEFS143B	CEFS143C	CEFS143D	CEFS143E	0.50mm	
	CEFS144A	CEFS144B	CEFS144C	CEFS144D	CEFS144E	1.00mm	
	CEFS145A	CEFS145B	CEFS145C	CEFS145D	CEFS145E	2.00mm	

Box of 20 Units



Insulated Hardy Bayonet Straight Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM137A	CEFM137B	CEFM137C	CEFM137D	CEFM137E	0.50mm	
	CEFM138A	CEFM138B	CEFM138C	CEFM138D	CEFM138E	1.00mm	
	CEFM139A	CEFM139B	CEFM139C	CEFM139D	CEFM139E	2.00mm	
Single-Use	CEFS137A	CEFS137B	CEFS137C	CEFS137D	CEFS137E	0.50mm	
	CEFS138A	CEFS138B	CEFS138C	CEFS138D	CEFS138E	1.00mm	
	CEFS139A	CEFS139B	CEFS139C	CEFS139D	CEFS139E	2.00mm	

Box of 20 Units



Insulated Cushing Grush Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm		
Multi-Use	CEFM146A	CEFM146B	CEFM146C	CEFM146D	CEFM146E	0.50mm	
	CEFM147A	CEFM147B	CEFM147C	CEFM147D	CEFM147E	1.00mm	
	CEFM148A	CEFM148B	CEFM148C	CEFM148D	CEFM148E	2.00mm	
Single-Use	CEFS146A	CEFS146B	CEFS146C	CEFS146D	CEFS146E	0.50mm	
	CEFS147A	CEFS147B	CEFS147C	CEFS147D	CEFS147E	1.00mm	
	CEFS148A	CEFS148B	CEFS148C	CEFS148D	CEFS148E	2.00mm	

Box of 20 Units



INSULATED BIPOLAR FORCEPS

TWO PIN CONNECTION



UMAX: 0.5KVp



Insulated Gerald Straight Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM149A	CEFM149B	CEFM149C	CEFM149D	CEFM149E	0.50mm
	CEFM150A	CEFM150B	CEFM150C	CEFM150D	CEFM150E	1.00mm
	CEFM151A	CEFM151B	CEFM151C	CEFM151D	CEFM151E	2.00mm
Single-Use	CEFS149A	CEFS149B	CEFS149C	CEFS149D	CEFS149E	0.50mm
	CEFS150A	CEFS150B	CEFS150C	CEFS150D	CEFS150E	1.00mm
	CEFS151A	CEFS151B	CEFS151C	CEFS151D	CEFS151E	2.00mm

Box of 20 Units



Adson Bipolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM158W	CEFM158A	CEFM158B	CEFM158C	0.50mm
	CEFM159W	CEFM159A	CEFM159B	CEFM159C	1.00mm
	CEFM160W	CEFM160A	CEFM160B	CEFM160C	2.00mm
Single-Use	CEFS158W	CEFS158A	CEFS158B	CEFS158C	0.50mm
	CEFS159W	CEFS159A	CEFS159B	CEFS159C	1.00mm
	CEFS160W	CEFS160A	CEFS160B	CEFS160C	2.00mm

Box of 20 Units



Insulated Gerald Straight Angled Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM152A	CEFM152B	CEFM152C	CEFM152D	CEFM152E	0.50mm
	CEFM153A	CEFM153B	CEFM153C	CEFM153D	CEFM153E	1.00mm
	CEFM154A	CEFM154B	CEFM154C	CEFM154D	CEFM154E	2.00mm
Single-Use	CEFS152A	CEFS152B	CEFS152C	CEFS152D	CEFS152E	0.50mm
	CEFS153A	CEFS153B	CEFS153C	CEFS153D	CEFS153E	1.00mm
	CEFS154A	CEFS154B	CEFS154C	CEFS154D	CEFS154E	2.00mm

Box of 20 Units



Smart Iris Str. Bipolar Forceps

	102mm	
Multi-Use	CEFM161T	0.50mm
	CEFM162T	1.00mm
	CEFM163T	2.00mm
Single-Use	CEFS161T	0.50mm
	CEFS162T	1.00mm
	CEFS163T	2.00mm

Box of 20 Units



Insulated Gerald Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM155A	CEFM155B	CEFM155C	CEFM155D	CEFM155E	0.50mm
	CEFM156A	CEFM156B	CEFM156C	CEFM156D	CEFM156E	1.00mm
	CEFM157A	CEFM157B	CEFM157C	CEFM157D	CEFM157E	2.00mm
Single-Use	CEFS155A	CEFS155B	CEFS155C	CEFS155D	CEFS155E	0.50mm
	CEFS156A	CEFS156B	CEFS156C	CEFS156D	CEFS156E	1.00mm
	CEFS157A	CEFS157B	CEFS157C	CEFS157D	CEFS157E	2.00mm

Box of 20 Units



Smart Iris Ang. Bipolar Forceps

	102mm	
Multi-Use	CEFM164T	0.50mm
	CEFM165T	1.00mm
	CEFM166T	2.00mm
Single-Use	CEFS164T	0.50mm
	CEFS165T	1.00mm
	CEFS166T	2.00mm

Box of 20 Units



01 ELECTROSURGICAL FORCEPS

INSULATED BIPOLAR FORCEPS

TWO & FLAT FLAT CONNECTION



UMAX: 0.5KVP



Irrigation Bayonet Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM200A	CEFM200B	CEFM200C	CEFM200D	CEFM200E	0.50mm
	CEFM201A	CEFM201B	CEFM201C	CEFM201D	CEFM201E	1.00mm
	CEFM202A	CEFM202B	CEFM202C	CEFM202D	CEFM202E	2.00mm
Single-Use	CEFS200A	CEFS200B	CEFS200C	CEFS200D	CEFS200E	0.50mm
	CEFS201A	CEFS201B	CEFS201C	CEFS201D	CEFS201E	1.00mm
	CEFS202A	CEFS202B	CEFS202C	CEFS202D	CEFS202E	2.00mm

Box of 20 Units



Irrigation Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM206A	CEFM206B	CEFM206C	CEFM206D	CEFM206E	0.50mm
	CEFM207A	CEFM207B	CEFM207C	CEFM207D	CEFM207E	1.00mm
	CEFM208A	CEFM208B	CEFM208C	CEFM208D	CEFM208E	2.00mm
Single-Use	CEFS206A	CEFS206B	CEFS206C	CEFS206D	CEFS206E	0.50mm
	CEFS207A	CEFS207B	CEFS207C	CEFS207D	CEFS207E	1.00mm
	CEFS208A	CEFS208B	CEFS208C	CEFS208D	CEFS208E	2.00mm

Box of 20 Units



Irrigation Straight Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM203A	CEFM203B	CEFM203C	CEFM203D	CEFM203E	0.50mm
	CEFM204A	CEFM204B	CEFM204C	CEFM204D	CEFM204E	1.00mm
	CEFM205A	CEFM205B	CEFM205C	CEFM205D	CEFM205E	2.00mm
Single-Use	CEFS203A	CEFS203B	CEFS203C	CEFS203D	CEFS203E	0.50mm
	CEFS204A	CEFS204B	CEFS204C	CEFS204D	CEFS204E	1.00mm
	CEFS205A	CEFS205B	CEFS205C	CEFS205D	CEFS205E	2.00mm

Box of 20 Units



Irrigation Jensen Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM209A	CEFM209B	CEFM209C	CEFM209D	CEFM209E	0.50mm
	CEFM210A	CEFM210B	CEFM210C	CEFM210D	CEFM210E	1.00mm
	CEFM211A	CEFM211B	CEFM211C	CEFM211D	CEFM211E	2.00mm
Single-Use	CEFS209A	CEFS209B	CEFS209C	CEFS209D	CEFS209E	0.50mm
	CEFS210A	CEFS210B	CEFS210C	CEFS210D	CEFS210E	1.00mm
	CEFS211A	CEFS211B	CEFS211C	CEFS211D	CEFS211E	2.00mm

Box of 20 Units



Irrigation Mals Bipolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM212A	CEFM212B	CEFM212C	CEFM212D	CEFM212E	0.50mm
	CEFM213A	CEFM213B	CEFM213C	CEFM213D	CEFM213E	1.00mm
	CEFM214A	CEFM214B	CEFM214C	CEFM214D	CEFM214E	2.00mm
Single-Use	CEFS212A	CEFS212B	CEFS212C	CEFS212D	CEFS212E	0.50mm
	CEFS213A	CEFS213B	CEFS213C	CEFS213D	CEFS213E	1.00mm
	CEFS214A	CEFS214B	CEFS214C	CEFS214D	CEFS214E	2.00mm

Box of 20 Units



INSULATED BIPOLAR FORCEPS

FEMALE BLOCK CONNECTION




UMAX: 0.5KVp




Bayonet Premium Tips Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM215A	CEFM215B	CEFM215C	CEFM215D	CEFM215E	0.50mm
	CEFM217A	CEFM217B	CEFM217C	CEFM217D	CEFM217E	1.00mm
	CEFM218A	CEFM218B	CEFM218C	CEFM218D	CEFM218E	2.00mm
Single-Use	CEFS215A	CEFS215B	CEFS215C	CEFS215D	CEFS215E	0.50mm
	CEFS217A	CEFS217B	CEFS217C	CEFS217D	CEFS217E	1.00mm
	CEFS218A	CEFS218B	CEFS218C	CEFS218D	CEFS218E	2.00mm




Bayonet Monopolar Angle Up Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM219A	CEFM219B	CEFM219C	CEFM219D	CEFM219E	0.50mm
	CEFM220A	CEFM220B	CEFM220C	CEFM220D	CEFM220E	1.00mm
	CEFM221A	CEFM221B	CEFM221C	CEFM221D	CEFM221E	2.00mm
Single-Use	CEFS219A	CEFS219B	CEFS219C	CEFS219D	CEFS219E	0.50mm
	CEFS220A	CEFS220B	CEFS220C	CEFS220D	CEFS220E	1.00mm
	CEFS221A	CEFS221B	CEFS221C	CEFS221D	CEFS221E	2.00mm



Insulated Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM222A	CEFM222B	CEFM222C	CEFM222D	CEFM222E	0.50mm
	CEFM223A	CEFM223B	CEFM223C	CEFM223D	CEFM223E	1.00mm
	CEFM224A	CEFM224B	CEFM224C	CEFM224D	CEFM224E	2.00mm
Single-Use	CEFS222A	CEFS222B	CEFS222C	CEFS222D	CEFS222E	0.50mm
	CEFS223A	CEFS223B	CEFS223C	CEFS223D	CEFS223E	1.00mm
	CEFS224A	CEFS224B	CEFS224C	CEFS224D	CEFS224E	2.00mm



INSULATED MONOPOLAR FORCEPS

BANANA PLUG CONNECTION




UMAX: 0.5KVp



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM225W	CEFM225A	CEFM225B	CEFM225C	0.50mm
	CEFM226W	CEFM226A	CEFM226B	CEFM226C	1.00mm
	CEFM227W	CEFM227A	CEFM227B	CEFM227C	2.00mm
Single-Use	CEFS225W	CEFS225A	CEFS225B	CEFS225C	0.50mm
	CEFS226W	CEFS226A	CEFS226B	CEFS226C	1.00mm
	CEFS227W	CEFS227A	CEFS227B	CEFS227C	2.00mm




Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM234A	CEFM234B	CEFM234C	CEFM234D	CEFM234E	0.50mm
	CEFM235A	CEFM235B	CEFM235C	CEFM235D	CEFM235E	1.00mm
	CEFM236A	CEFM236B	CEFM236C	CEFM236D	CEFM236E	2.00mm
Single-Use	CEFS234A	CEFS234B	CEFS234C	CEFS234D	CEFS234E	0.50mm
	CEFS235A	CEFS235B	CEFS235C	CEFS235D	CEFS235E	1.00mm
	CEFS236A	CEFS236B	CEFS236C	CEFS236D	CEFS236E	2.00mm

McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM228A	CEFM228B	0.50mm
	CEFM229A	CEFM229B	1.00mm
	CEFM230A	CEFM230B	2.00mm
Single-Use	CEFS228A	CEFS228B	0.50mm
	CEFS229A	CEFS229B	1.00mm
	CEFS230A	CEFS230B	2.00mm




Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM237A	CEFM237B	CEFM237C	0.50mm
	CEFM238A	CEFM238B	CEFM238C	1.00mm
	CEFM239A	CEFM239B	CEFM239C	2.00mm
Single-Use	CEFS237A	CEFS237B	CEFS237C	0.50mm
	CEFS238A	CEFS238B	CEFS238C	1.00mm
	CEFS239A	CEFS239B	CEFS239C	2.00mm

Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM231A	CEFM231B	CEFM231C	0.50mm
	CEFM232A	CEFM232B	CEFM232C	1.00mm
	CEFM233A	CEFM233B	CEFM233C	2.00mm
Single-Use	CEFS231A	CEFS231B	CEFS231C	0.50mm
	CEFS232A	CEFS232B	CEFS232C	1.00mm
	CEFS233A	CEFS233B	CEFS233C	2.00mm



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM240A	CEFM240B	CEFM240C	CEFM240D	CEFM240E	0.50mm
	CEFM241A	CEFM241B	CEFM241C	CEFM241D	CEFM241E	1.00mm
	CEFM242A	CEFM242B	CEFM242C	CEFM242D	CEFM242E	2.00mm
Single-Use	CEFS240A	CEFS240B	CEFS240C	CEFS240D	CEFS240E	0.50mm
	CEFS241A	CEFS241B	CEFS241C	CEFS241D	CEFS241E	1.00mm
	CEFS242A	CEFS242B	CEFS242C	CEFS242D	CEFS242E	2.00mm



INSULATED MONOPOLAR FORCEPS

BANANA PLUG CONNECTION




UMAX: 0.5KVp



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM225W	CEFM225A	CEFM225B	CEFM225C	0.50mm
	CEFM226W	CEFM226A	CEFM226B	CEFM226C	1.00mm
	CEFM227W	CEFM227A	CEFM227B	CEFM227C	2.00mm
Single-Use	CEFS225W	CEFS225A	CEFS225B	CEFS225C	0.50mm
	CEFS226W	CEFS226A	CEFS226B	CEFS226C	1.00mm
	CEFS227W	CEFS227A	CEFS227B	CEFS227C	2.00mm




Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM234A	CEFM234B	CEFM234C	CEFM234D	CEFM234E	0.50mm
	CEFM235A	CEFM235B	CEFM235C	CEFM235D	CEFM235E	1.00mm
	CEFM236A	CEFM236B	CEFM236C	CEFM236D	CEFM236E	2.00mm
Single-Use	CEFS234A	CEFS234B	CEFS234C	CEFS234D	CEFS234E	0.50mm
	CEFS235A	CEFS235B	CEFS235C	CEFS235D	CEFS235E	1.00mm
	CEFS236A	CEFS236B	CEFS236C	CEFS236D	CEFS236E	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM228A	CEFM228B	0.50mm
	CEFM229A	CEFM229B	1.00mm
	CEFM230A	CEFM230B	2.00mm
Single-Use	CEFS228A	CEFS228B	0.50mm
	CEFS229A	CEFS229B	1.00mm
	CEFS230A	CEFS230B	2.00mm




Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM237A	CEFM237B	CEFM237C	0.50mm
	CEFM238A	CEFM238B	CEFM238C	1.00mm
	CEFM239A	CEFM239B	CEFM239C	2.00mm
Single-Use	CEFS237A	CEFS237B	CEFS237C	0.50mm
	CEFS238A	CEFS238B	CEFS238C	1.00mm
	CEFS239A	CEFS239B	CEFS239C	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM231A	CEFM231B	CEFM231C	0.50mm
	CEFM232A	CEFM232B	CEFM232C	1.00mm
	CEFM233A	CEFM233B	CEFM233C	2.00mm
Single-Use	CEFS231A	CEFS231B	CEFS231C	0.50mm
	CEFS232A	CEFS232B	CEFS232C	1.00mm
	CEFS233A	CEFS233B	CEFS233C	2.00mm



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM240A	CEFM240B	CEFM240C	CEFM240D	CEFM240E	0.50mm
	CEFM241A	CEFM241B	CEFM241C	CEFM241D	CEFM241E	1.00mm
	CEFM242A	CEFM242B	CEFM242C	CEFM242D	CEFM242E	2.00mm
Single-Use	CEFS240A	CEFS240B	CEFS240C	CEFS240D	CEFS240E	0.50mm
	CEFS241A	CEFS241B	CEFS241C	CEFS241D	CEFS241E	1.00mm
	CEFS242A	CEFS242B	CEFS242C	CEFS242D	CEFS242E	2.00mm



INSULATED MONOPOLAR FORCEPS

BANANA PLUG CONNECTION




UMAX: 0.5KVP



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM243A	CEFM243B	CEFM243C	CEFM243D	CEFM243E	0.50mm
	CEFM244A	CEFM244B	CEFM244C	CEFM244D	CEFM244E	1.00mm
	CEFM245A	CEFM245B	CEFM245C	CEFM245D	CEFM245E	2.00mm
Single-Use	CEFS243A	CEFS243B	CEFS243C	CEFS243D	CEFS243E	0.50mm
	CEFS244A	CEFS244B	CEFS244C	CEFS244D	CEFS244E	1.00mm
	CEFS245A	CEFS245B	CEFS245C	CEFS245D	CEFS245E	2.00mm




FOR ORDER
Please mention tip code along with product code



Regiminton Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM246C	CEFM246D	CEFM246E	0.50mm
	CEFM247C	CEFM247D	CEFM247E	1.00mm
	CEFM248C	CEFM248D	CEFM248E	2.00mm
Single-Use	CEFS246C	CEFS246D	CEFS246E	0.50mm
	CEFS247C	CEFS247D	CEFS247E	1.00mm
	CEFS248C	CEFS248D	CEFS248E	2.00mm



FOR ORDER
Please mention tip code along with product code



Monopolar Grasping Holder 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	CEFM248A	CEFM248B	CEFM248C	CEFM248D	CEFM248E
	CEFM249A	CEFM249B	CEFM249C	CEFM249D	CEFM249E
	CEFM248A	CEFM248B	CEFM248C	CEFM248D	CEFM248E
Single-Use	CEFS248A	CEFS248B	CEFS248C	CEFS248D	CEFS248E
	CEFS249A	CEFS249B	CEFS249C	CEFS249D	CEFS249E
	CEFS248A	CEFS248B	CEFS248C	CEFS248D	CEFS248E



Monopolar Dressing Scissors 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	CEFM250A	CEFM250B	CEFM250C	CEFM250D	CEFM250E
	CEFM251A	CEFM251B	CEFM251C	CEFM251D	CEFM251E
	CEFM250A	CEFM250B	CEFM250C	CEFM250D	CEFM250E
Single-Use	CEFS250A	CEFS250B	CEFS250C	CEFS250D	CEFS250E
	CEFS251A	CEFS251B	CEFS251C	CEFS251D	CEFS251E
	CEFS250A	CEFS250B	CEFS250C	CEFS250D	CEFS250E



Monopolar Riche Forceps 4.8mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	CEFM252A	CEFM252B	CEFM252C	CEFM252D	CEFM252E
	CEFM253A	CEFM253B	CEFM253C	CEFM253D	CEFM253E
	CEFM252A	CEFM252B	CEFM252C	CEFM252D	CEFM252E
Single-Use	CEFS252A	CEFS252B	CEFS252C	CEFS252D	CEFS252E
	CEFS253A	CEFS253B	CEFS253C	CEFS253D	CEFS253E
	CEFS252A	CEFS252B	CEFS252C	CEFS252D	CEFS252E

INSULATED MONOPOLAR FORCEPS

WITHOUT CONNECTOR



UMAX: 0.5KVp



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM272A	CEFM272B	CEFM272C	CEFM272D	CEFM272E	0.50mm
	CEFM273A	CEFM273B	CEFM273C	CEFM273D	CEFM273E	1.00mm
	CEFM274A	CEFM274B	CEFM274C	CEFM274D	CEFM274E	2.00mm
Single-Use	CEFS272A	CEFS272B	CEFS272C	CEFS272D	CEFS272E	0.50mm
	CEFS273A	CEFS273B	CEFS273C	CEFS273D	CEFS273E	1.00mm
	CEFS274A	CEFS274B	CEFS274C	CEFS274D	CEFS274E	2.00mm



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM278W	CEFM278A	CEFM278B	CEFM278C	0.50mm
	CEFM279W	CEFM279A	CEFM279B	CEFM279C	1.00mm
	CEFM280W	CEFM280A	CEFM280B	CEFM280C	2.00mm
Single-Use	CEFS278W	CEFS278A	CEFS278B	CEFS278C	0.50mm
	CEFS279W	CEFS279A	CEFS279B	CEFS279C	1.00mm
	CEFS280W	CEFS280A	CEFS280B	CEFS280C	2.00mm



Regimont Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM275C	CEFM275D	CEFM275E	0.50mm
	CEFM276C	CEFM276D	CEFM276E	1.00mm
	CEFM277C	CEFM277D	CEFM277E	2.00mm
Single-Use	CEFS275C	CEFS275D	CEFS275E	0.50mm
	CEFS276C	CEFS276D	CEFS276E	1.00mm
	CEFS277C	CEFS277D	CEFS277E	2.00mm



McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM281A	CEFM281B	0.50mm
	CEFM282A	CEFM282B	1.00mm
	CEFM283A	CEFM283B	2.00mm
Single-Use	CEFS281A	CEFS281B	0.50mm
	CEFS282A	CEFS282B	1.00mm
	CEFS283A	CEFS283B	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM284A	CEFM284B	CEFM284C	0.50mm
	CEFM285A	CEFM285B	CEFM285C	1.00mm
	CEFM286A	CEFM286B	CEFM286C	2.00mm
Single-Use	CEFS284A	CEFS284B	CEFS284C	0.50mm
	CEFS285A	CEFS285B	CEFS285C	1.00mm
	CEFS286A	CEFS286B	CEFS286C	2.00mm



INSULATED MONOPOLAR FORCEPS

FEMALE PLUG CONNECTION




UMAX: 0.5KVp



Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM287A CEFM288A CEFM289A	CEFM287B CEFM288B CEFM289B	CEFM287C CEFM288C CEFM289C	CEFM287D CEFM288D CEFM289D	CEFM287E CEFM288E CEFM289E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS287A CEFS288A CEFS289A	CEFS287B CEFS288B CEFS289B	CEFS287C CEFS288C CEFS289C	CEFS287D CEFS288D CEFS289D	CEFS287E CEFS288E CEFS289E	0.50mm 1.00mm 2.00mm




FOR ORDER
Please mention tip code along with product code



Gerald Straight 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM290A CEFM291A CEFM292A	CEFM290B CEFM291B CEFM292B	CEFM290C CEFM291C CEFM292C	0.50mm 1.00mm 2.00mm
Single-Use	CEFS290A CEFS291A CEFS292A	CEFS290B CEFS291B CEFS292B	CEFS290C CEFS291C CEFS292C	0.50mm 1.00mm 2.00mm




FOR ORDER
Please mention tip code along with product code



Dressing s Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM293A CEFM294A CEFM295A	CEFM293B CEFM294B CEFM295B	CEFM293C CEFM294C CEFM295C	CEFM293D CEFM294D CEFM295D	CEFM293E CEFM294E CEFM295E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS293A CEFS294A CEFS295A	CEFS293B CEFS294B CEFS295B	CEFS293C CEFS294C CEFS295C	CEFS293D CEFS294D CEFS295D	CEFS293E CEFS294E CEFS295E	0.50mm 1.00mm 2.00mm



FOR ORDER
Please mention tip code along with product code



Bayonet Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM296A CEFM297A CEFM298A	CEFM296B CEFM297B CEFM298B	CEFM296C CEFM297C CEFM298C	CEFM296D CEFM297D CEFM298D	CEFM296E CEFM297E CEFM298E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS296A CEFS297A CEFS298A	CEFS296B CEFS297B CEFS298B	CEFS296C CEFS297C CEFS298C	CEFS296D CEFS297D CEFS298D	CEFS296E CEFS297E CEFS298E	0.50mm 1.00mm 2.00mm



Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM299A CEFM300A CEFM301A	CEFM299B CEFM300B CEFM301B	CEFM299C CEFM300C CEFM301C	CEFM299D CEFM300D CEFM301D	CEFM299E CEFM300E CEFM301E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS299A CEFS300A CEFS301A	CEFS299B CEFS300B CEFS301B	CEFS299C CEFS300C CEFS301C	CEFS299D CEFS300D CEFS301D	CEFS299E CEFS300E CEFS301E	0.50mm 1.00mm 2.00mm



Regimont Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM302C CEFM303C CEFM304C	CEFM302D CEFM303D CEFM304D	CEFM302E CEFM303E CEFM304E	0.50mm 1.00mm 2.00mm
Single-Use	CEFS302C CEFS303C CEFS304C	CEFS302D CEFS303D CEFS304D	CEFS302E CEFS303E CEFS304E	0.50mm 1.00mm 2.00mm

INSULATED MONOPOLAR FORCEPS

4.8MM PLUG CONNECTION




UMAX: 0.5KVp



Adson Monopolar Forceps

	121mm	155mm	180mm	200mm	
Multi-Use	CEFM305W	CEFM305A	CEFM305B	CEFM305C	0.50mm
	CEFM306W	CEFM306A	CEFM306B	CEFM306C	1.00mm
	CEFM307W	CEFM307A	CEFM307B	CEFM307C	2.00mm
Single-Use	CEFS305W	CEFS305A	CEFS305B	CEFS305C	0.50mm
	CEFS306W	CEFS306A	CEFS306B	CEFS306C	1.00mm
	CEFS307W	CEFS307A	CEFS307B	CEFS307C	2.00mm




Dressing Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM314A	CEFM314B	CEFM314C	CEFM314D	CEFM314E	0.50mm
	CEFM315A	CEFM315B	CEFM315C	CEFM315D	CEFM315E	1.00mm
	CEFM316A	CEFM316B	CEFM316C	CEFM316D	CEFM316E	2.00mm
Single-Use	CEFS314A	CEFS314B	CEFS314C	CEFS314D	CEFS314E	0.50mm
	CEFS315A	CEFS315B	CEFS315C	CEFS315D	CEFS315E	1.00mm
	CEFS316A	CEFS316B	CEFS316C	CEFS316D	CEFS316E	2.00mm




McIndoe Monopolar Forceps

	155mm	180mm	
Multi-Use	CEFM308A	CEFM308B	0.50mm
	CEFM309A	CEFM309B	1.00mm
	CEFM310A	CEFM310B	2.00mm
Single-Use	CEFS308A	CEFS308B	0.50mm
	CEFS309A	CEFS309B	1.00mm
	CEFS310A	CEFS310B	2.00mm



Gillies 2/1 Tip Monopolar Forceps

	155mm	180mm	200mm	
Multi-Use	CEFM311A	CEFM311B	CEFM311C	0.50mm
	CEFM312A	CEFM312B	CEFM312C	1.00mm
	CEFM313A	CEFM313B	CEFM313C	2.00mm
Single-Use	CEFS311A	CEFS311B	CEFS311C	0.50mm
	CEFS312A	CEFS312B	CEFS312C	1.00mm
	CEFS313A	CEFS313B	CEFS313C	2.00mm



INSULATED MONOPOLAR FORCEPS

4.8MM PLUG CONNECTION




UMAX: 0.5KVp




Bayonet Angled Tip Monopolar Forceps

	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CEFM323A	CEFM323B	CEFM323C	CEFM323D	CEFM323E	0.50mm
	CEFM324A	CEFM324B	CEFM324C	CEFM324D	CEFM324E	1.00mm
	CEFM325A	CEFM325B	CEFM325C	CEFM325D	CEFM325E	2.00mm
Single-Use	CEFS323A	CEFS323B	CEFS323C	CEFS323D	CEFS323E	0.50mm
	CEFS324A	CEFS324B	CEFS324C	CEFS324D	CEFS324E	1.00mm
	CEFS325A	CEFS325B	CEFS325C	CEFS325D	CEFS325E	2.00mm



Regimonton Monopolar Forceps

	200mm	220mm	240mm	
Multi-Use	CEFM326C	CEFM326D	CEFM326E	0.50mm
	CEFM327C	CEFM327D	CEFM327E	1.00mm
	CEFM328C	CEFM328D	CEFM328E	2.00mm
Single-Use	CEFS326C	CEFS326D	CEFS326E	0.50mm
	CEFS327C	CEFS327D	CEFS327E	1.00mm
	CEFS328C	CEFS328D	CEFS328E	2.00mm



Monopolar Grasping Holder 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM328A	CEFM328B	CEFM328C	CEFM328D	CEFM328E
	Angled CEFM329A	CEFM329B	CEFM329C	CEFM329D	CEFM329E
Single-Use	Straight CEFS328A	CEFS328B	CEFS328C	CEFS328D	CEFS328E
	Angled CEFS329A	CEFS329B	CEFS329C	CEFS329D	CEFS329E



Monopolar Dressing Scissors 4mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM330A	CEFM330B	CEFM330C	CEFM330D	CEFM330E
	Angled CEFM331A	CEFM331B	CEFM331C	CEFM331D	CEFM331E
Single-Use	Straight CEFS330A	CEFS330B	CEFS330C	CEFS330D	CEFS330E
	Angled CEFS331A	CEFS331B	CEFS331C	CEFS331D	CEFS331E



Monopolar Riche Forceps 4.8mm

	155mm	180mm	200mm	220mm	240mm
Multi-Use	Straight CEFM332A	CEFM332B	CEFM332C	CEFM332D	CEFM332E
	Angled CEFM333A	CEFM333B	CEFM333C	CEFM333D	CEFM333E
Single-Use	Straight CEFS332A	CEFS332B	CEFS332C	CEFS332D	CEFS332E
	Angled CEFS333A	CEFS333B	CEFS333C	CEFS333D	CEFS333E





02

C-TECH BIPOLAR INSTRUMENTS

PRECISE. INSULATED. RELIABLE.

C-TECH
BIPOLAR INSTRUMENTS
SEALER



UMAX: 0.5KVP



Bipolar Vessel Sealer Stainless Steel

Curve

Angle

4.0mm Male Plug



BOX OF 20 UNITS

4.0mm Female Plug



BOX OF 20 UNITS

2.4mm Ethicon Plug



BOX OF 20 UNITS

	165mm	190mm	215mm	240mm
Curve	CPFM101K	CPFM101L	CPFM101M	CPFM101N
Angle	CPFM102K	CPFM102K	CPFM102M	CPFM102N
Single-Use	CPFS101K	CPFS101L	CPFS101M	CPFS101N
Multi-Use	CPFS102K	CPFS102L	CPFS102M	CPFS102N

	165mm	190mm	215mm	240mm
Curve	CPFM103K	CPFM103L	CPFM103M	CPFM103N
Angle	CPFM104K	CPFM104L	CPFM104M	CPFM104N
Single-Use	CPFS103K	CPFS103L	CPFS103M	CPFS103N
Multi-Use	CPFS104K	CPFS104L	CPFS104M	CPFS104N

	165mm	190mm	215mm	240mm
Curve	CPFM105K	CPFM105L	CPFM105M	CPFM105N
Angle	CPFM106K	CPFM106L	CPFM106M	CPFM106N
Single-Use	CPFS105K	CPFS105L	CPFS105M	CPFS105N
Multi-Use	CPFS106K	CPFS106L	CPFS106M	CPFS106N



C-TECH NON-STICK
BIPOLAR INSTRUMENTS
SEALER



UMAX: 0.5KVP



Non-Stick Bipolar Vessel Sealer Stainless Steel

4.0mm Male Plug




BOX OF 20 UNITS

4.0mm Female Plug



BOX OF 20 UNITS

2.4mm Ethicon Plug



BOX OF 20 UNITS

	165mm	190mm	215mm	240mm
Single-Use	Curve CPFM107K	CPFM107L	CPFM107M	CPFM107N
Multi-Use	Angle CPFM108K	CPFM108L	CPFM108M	CPFM108N
Single-Use	Curve CPFS107K	CPFS107L	CPFS107M	CPFS107N
Multi-Use	Angle CPFS108K	CPFS108L	CPFS108M	CPFS108N

	165mm	190mm	215mm	240mm
Single-Use	Curve CPFM109K	CPFM109L	CPFM109M	CPFM109N
Multi-Use	Angle CPFM110K	CPFM110L	CPFM110M	CPFM110N
Single-Use	Curve CPFS109K	CPFS109L	CPFS109M	CPFS109N
Multi-Use	Angle CPFS110K	CPFS110L	CPFS110M	CPFS110N


	165mm	190mm	215mm	240mm
Single-Use	Curve CPFM111K	CPFM111L	CPFM111M	CPFM111N
Multi-Use	Angle CPFM112K	CPFM112L	CPFM112M	CPFM112N
Single-Use	Curve CPFS111K	CPFS111L	CPFS111M	CPFS111N
Multi-Use	Angle CPFS112K	CPFS112L	CPFS112M	CPFS112N



C-TECH NON-STICK
BIPOLAR INSTRUMENTS
SCISSOR



UMAX: 0.5KVP



Bipolar Cut-Coag Scissor -Ceramics

4.0mm Male Plug



BOX OF 20 UNITS

4.0mm Female Plug



BOX OF 20 UNITS

2.4mm Ethicon Plug



BOX OF 20 UNITS

Single-Use / Multi-Use	180mm	200mm	220mm
	CPFM113P	CPFM113Q	CPFM113R
	CPFS113P	CPFS113Q	CPFS113R

Single-Use / Multi-Use	180mm	200mm	220mm
	CPFM114P	CPFM114Q	CPFM114R
	CPFS114P	CPFS114Q	CPFS114R

Single-Use / Multi-Use	180mm	200mm	220mm
	CPFM115P	CPFM115Q	CPFM115R
	CPFS115P	CPFS115Q	CPFS115R



C-TECH BIPOLAR INSTRUMENTS

FLAT PLUG CONNECTION



UMAX: 0.5KVP



102mm		115mm		
Multi-Use	CPFM116T	CPFM116V	0.50mm	
	CPFM117T	CPFM117V	1.00mm	
	CPFM118T	CPFM118V	2.00mm	
Single-Use	CPFS116T	CPFS116V	0.50mm	
	CPFS117T	CPFS117V	1.00mm	
	CPFS118T	CPFS118V	2.00mm	

Box of 20 Units



102mm		115mm		
Multi-Use	CPFM119T	CPFM119V	0.50mm	
	CPFM120T	CPFM120V	1.00mm	
	CPFM121T	CPFM121V	2.00mm	
Single-Use	CPFS119T	CPFS119V	0.50mm	
	CPFS120T	CPFS120V	1.00mm	
	CPFS121T	CPFS121V	2.00mm	

Box of 20 Units



140mm		155mm		180mm		200mm		220mm		240mm		
Multi-Use	CPFM122X	CPFM122A	CPFM122B	CPFM122C	CPFM122D	CPFM122P	0.50mm					
	CPFM123X	CPFM123A	CPFM123B	CPFM123C	CPFM123D	CPFM123P	1.00mm					
	CPFM124X	CPFM124A	CPFM124B	CPFM124C	CPFM124D	CPFM124P	2.00mm					
Single-Use	CPFS122X	CPFS122A	CPFS122B	CPFS122C	CPFS122D	CPFS122P	0.50mm					
	CPFS123X	CPFS123A	CPFS123B	CPFS123C	CPFS123D	CPFS123P	1.00mm					
	CPFS124X	CPFS124A	CPFS124B	CPFS124C	CPFS124D	CPFS124P	2.00mm					

Box of 20 Units



102mm		115mm		
Multi-Use	CPFM125T	CPFM125V	0.50mm	
	CPFM126T	CPFM126V	1.00mm	
	CPFM127T	CPFM127V	2.00mm	
Single-Use	CPFS125T	CPFS125V	0.50mm	
	CPFS126T	CPFS126V	1.00mm	
	CPFS127T	CPFS127V	2.00mm	

Box of 20 Units



102mm		115mm		
Multi-Use	CPFM128T	CPFM128V	0.50mm	
	CPFM129T	CPFM129V	1.00mm	
	CPFM130T	CPFM130V	2.00mm	
Single-Use	CPFS128T	CPFS128V	0.50mm	
	CPFS129T	CPFS129V	1.00mm	
	CPFS130T	CPFS130V	2.00mm	

Box of 20 Units



140mm		155mm		180mm		200mm		220mm		240mm		
Multi-Use	CPFM131X	CPFM131A	CPFM131B	CPFM131C	CPFM131D	CPFM131P	0.50mm					
	CPFM132X	CPFM132A	CPFM132B	CPFM132C	CPFM132D	CPFM132P	1.00mm					
	CPFM133X	CPFM133A	CPFM133B	CPFM133C	CPFM133D	CPFM133P	2.00mm					
Single-Use	CPFS131X	CPFS131A	CPFS131B	CPFS131C	CPFS131D	CPFS131P	0.50mm					
	CPFS132X	CPFS132A	CPFS132B	CPFS132C	CPFS132D	CPFS132P	1.00mm					
	CPFS133X	CPFS133A	CPFS133B	CPFS133C	CPFS133D	CPFS133P	2.00mm					

Box of 20 Units



C-TECH
BIPOLAR INSTRUMENTS

TWO PIN CONNECTION



UMAX: 0.5KVP



	102mm	115mm	
Multi-Use	CPFM134T	CPFM134V	0.50mm
	CPFM135T	CPFM135V	1.00mm
	CPFM136T	CPFM136V	2.00mm
Single-Use	CPFS134T	CPFS134V	0.50mm
	CPFS135T	CPFS135V	1.00mm
	CPFS136T	CPFS136V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM137T	CPFM137V	0.50mm
	CPFM138T	CPFM138V	1.00mm
	CPFM139T	CPFM139V	2.00mm
Single-Use	CPFS137T	CPFS137V	0.50mm
	CPFS138T	CPFS138V	1.00mm
	CPFS139T	CPFS139V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM140X	CPFM140A	CPFM140B	CPFM140C	CPFM140D	CPFM140P	0.50mm
	CPFM141X	CPFM141A	CPFM141B	CPFM141C	CPFM141D	CPFM141P	1.00mm
	CPFM142X	CPFM142A	CPFM142B	CPFM142C	CPFM142D	CPFM142P	2.00mm
Single-Use	CPFS140X	CPFS140A	CPFS140B	CPFS140C	CPFS140D	CPFS140P	0.50mm
	CPFS141X	CPFS141A	CPFS141B	CPFS141C	CPFS141D	CPFS141P	1.00mm
	CPFS142X	CPFS142A	CPFS142B	CPFS142C	CPFS142D	CPFS142P	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM143T	CPFM143V	0.50mm
	CPFM144T	CPFM144V	1.00mm
	CPFM145T	CPFM145V	2.00mm
Single-Use	CPFS143T	CPFS143V	0.50mm
	CPFS144T	CPFS144V	1.00mm
	CPFS145T	CPFS145V	2.00mm

Box of 20 Units



	102mm	115mm	
Multi-Use	CPFM146T	CPFM146V	0.50mm
	CPFM147T	CPFM147V	1.00mm
	CPFM148T	CPFM148V	2.00mm
Single-Use	CPFS146T	CPFS146V	0.50mm
	CPFS147T	CPFS147V	1.00mm
	CPFS148T	CPFS148V	2.00mm

Box of 20 Units



	140mm	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM149X	CPFM149A	CPFM149B	CPFM149C	CPFM149D	CPFM149P	0.50mm
	CPFM150X	CPFM150A	CPFM150B	CPFM150C	CPFM150D	CPFM150P	1.00mm
	CPFM151X	CPFM151A	CPFM151B	CPFM151C	CPFM151D	CPFM151P	2.00mm
Single-Use	CPFS149X	CPFS149A	CPFS149B	CPFS149C	CPFS149D	CPFS149P	0.50mm
	CPFS150X	CPFS150A	CPFS150B	CPFS150C	CPFS150D	CPFS150P	1.00mm
	CPFS151X	CPFS151A	CPFS151B	CPFS151C	CPFS151D	CPFS151P	2.00mm

Box of 20 Units



C-TECH
BIPOLAR INSTRUMENTS

FLAT PLUG CONNECTION



UMAX: 0.5KVP



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM152A	CPFM152B	CPFM152C	CPFM152D	CPFM152P	0.50mm
	CPFM153A	CPFM153B	CPFM153C	CPFM153D	CPFM153P	1.00mm
	CPFM154A	CPFM154B	CPFM154C	CPFM154D	CPFM154P	2.00mm
Single-Use	CPFS152A	CPFS152B	CPFS152C	CPFS152D	CPFS152P	0.50mm
	CPFS153A	CPFS153B	CPFS153C	CPFS153D	CPFS153P	1.00mm
	CPFS154A	CPFS154B	CPFS154C	CPFS154D	CPFS154P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM155A	CPFM155B	CPFM155C	CPFM155D	CPFM155P	0.50mm
	CPFM156A	CPFM156B	CPFM156C	CPFM156D	CPFM156P	1.00mm
	CPFM157A	CPFM157B	CPFM157C	CPFM157D	CPFM157P	2.00mm
Single-Use	CPFS155A	CPFS155B	CPFS155C	CPFS155D	CPFS155P	0.50mm
	CPFS156A	CPFS156B	CPFS156C	CPFS156D	CPFS156P	1.00mm
	CPFS157A	CPFS157B	CPFS157C	CPFS157D	CPFS157P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM158A	CPFM158B	CPFM158C	CPFM158D	CPFM158P	0.50mm
	CPFM159A	CPFM159B	CPFM159C	CPFM159D	CPFM159P	1.00mm
	CPFM160A	CPFM160B	CPFM160C	CPFM160D	CPFM160P	2.00mm
Single-Use	CPFS158A	CPFS158B	CPFS158C	CPFS158D	CPFS158P	0.50mm
	CPFS159A	CPFS159B	CPFS159C	CPFS159D	CPFS159P	1.00mm
	CPFS160A	CPFS160B	CPFS160C	CPFS160D	CPFS160P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM161A	CPFM161B	CPFM161C	CPFM161D	CPFM161P	0.50mm
	CPFM162A	CPFM162B	CPFM162C	CPFM162D	CPFM162P	1.00mm
	CPFM163A	CPFM163B	CPFM163C	CPFM163D	CPFM163P	2.00mm
Single-Use	CPFS161A	CPFS161B	CPFS161C	CPFS161D	CPFS161P	0.50mm
	CPFS162A	CPFS162B	CPFS162C	CPFS162D	CPFS162P	1.00mm
	CPFS163A	CPFS163B	CPFS163C	CPFS163D	CPFS163P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM164A	CPFM164B	CPFM164C	CPFM164D	CPFM164P	0.50mm
	CPFM165A	CPFM165B	CPFM165C	CPFM165D	CPFM165P	1.00mm
	CPFM166A	CPFM166B	CPFM166C	CPFM166D	CPFM166P	2.00mm
Single-Use	CPFS164A	CPFS164B	CPFS164C	CPFS164D	CPFS164P	0.50mm
	CPFS165A	CPFS165B	CPFS165C	CPFS165D	CPFS165P	1.00mm
	CPFS166A	CPFS166B	CPFS166C	CPFS166D	CPFS166P	2.00mm

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

TWO PIN CONNECTION



UMAX: 0.5KVp



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM167A	CPFM167B	CPFM167C	CPFM167D	CPFM167P	0.50mm
	CPFM168A	CPFM168B	CPFM168C	CPFM168D	CPFM168P	1.00mm
	CPFM169A	CPFM169B	CPFM169C	CPFM169D	CPFM169P	2.00mm
Single-Use	CPFS167A	CPFS167B	CPFS167C	CPFS167D	CPFS167P	0.50mm
	CPFS168A	CPFS168B	CPFS168C	CPFS168D	CPFS168P	1.00mm
	CPFS169A	CPFS169B	CPFS169C	CPFS169D	CPFS169P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM170A	CPFM170B	CPFM170C	CPFM170D	CPFM170P	0.50mm
	CPFM171A	CPFM171B	CPFM171C	CPFM171D	CPFM171P	1.00mm
	CPFM172A	CPFM172B	CPFM172C	CPFM172D	CPFM172P	2.00mm
Single-Use	CPFS170A	CPFS170B	CPFS170C	CPFS170D	CPFS170P	0.50mm
	CPFS171A	CPFS171B	CPFS171C	CPFS171D	CPFS171P	1.00mm
	CPFS172A	CPFS172B	CPFS172C	CPFS172D	CPFS172P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM173A	CPFM173B	CPFM173C	CPFM173D	CPFM173P	0.50mm
	CPFM174A	CPFM174B	CPFM174C	CPFM174D	CPFM174P	1.00mm
	CPFM175A	CPFM175B	CPFM175C	CPFM175D	CPFM175P	2.00mm
Single-Use	CPFS173A	CPFS173B	CPFS173C	CPFS173D	CPFS173P	0.50mm
	CPFS174A	CPFS174B	CPFS174C	CPFS174D	CPFS174P	1.00mm
	CPFS175A	CPFS175B	CPFS175C	CPFS175D	CPFS175P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM176A	CPFM176B	CPFM176C	CPFM176D	CPFM176P	0.50mm
	CPFM177A	CPFM177B	CPFM177C	CPFM177D	CPFM177P	1.00mm
	CPFM178A	CPFM178B	CPFM178C	CPFM178D	CPFM178P	2.00mm
Single-Use	CPFS176A	CPFS176B	CPFS176C	CPFS176D	CPFS176P	0.50mm
	CPFS177A	CPFS177B	CPFS177C	CPFS177D	CPFS177P	1.00mm
	CPFS178A	CPFS178B	CPFS178C	CPFS178D	CPFS178P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM179A	CPFM179B	CPFM179C	CPFM179D	CPFM179P	0.50mm
	CPFM180A	CPFM180B	CPFM180C	CPFM180D	CPFM180P	1.00mm
	CPFM181A	CPFM181B	CPFM181C	CPFM181D	CPFM181P	2.00mm
Single-Use	CPFS179A	CPFS179B	CPFS179C	CPFS179D	CPFS179P	0.50mm
	CPFS180A	CPFS180B	CPFS180C	CPFS180D	CPFS180P	1.00mm
	CPFS181A	CPFS181B	CPFS181C	CPFS181D	CPFS181P	2.00mm

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

NON-STICK IRRIGATION FORCEPS



UMAX: 0.5KVP



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM182A	CPFM182B	CPFM182C	CPFM182D	CPFM182P	0.50mm
	CPFM183A	CPFM183B	CPFM183C	CPFM183D	CPFM183P	1.00mm
	CPFM184A	CPFM184B	CPFM184C	CPFM184D	CPFM184P	2.00mm
Single-Use	CPFS182A	CPFS182B	CPFS182C	CPFS182D	CPFS182P	0.50mm
	CPFS183A	CPFS183B	CPFS183C	CPFS183D	CPFS183P	1.00mm
	CPFS184A	CPFS184B	CPFS184C	CPFS184D	CPFS184P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM188A	CPFM188B	CPFM188C	CPFM188D	CPFM188P	0.50mm
	CPFM189A	CPFM189B	CPFM189C	CPFM189D	CPFM189P	1.00mm
	CPFM190A	CPFM190B	CPFM190C	CPFM190D	CPFM190P	2.00mm
Single-Use	CPFS188A	CPFS188B	CPFS188C	CPFS188D	CPFS188P	0.50mm
	CPFS189A	CPFS189B	CPFS189C	CPFS189D	CPFS189P	1.00mm
	CPFS190A	CPFS190B	CPFS190C	CPFS190D	CPFS190P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM185A	CPFM185B	CPFM185C	CPFM185D	CPFM185P	0.50mm
	CPFM186A	CPFM186B	CPFM186C	CPFM186D	CPFM186P	1.00mm
	CPFM187A	CPFM187B	CPFM187C	CPFM187D	CPFM187P	2.00mm
Single-Use	CPFS185A	CPFS185B	CPFS185C	CPFS185D	CPFS185P	0.50mm
	CPFS186A	CPFS186B	CPFS186C	CPFS186D	CPFS186P	1.00mm
	CPFS187A	CPFS187B	CPFS187C	CPFS187D	CPFS187P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM191A	CPFM191B	CPFM191C	CPFM191D	CPFM191P	0.50mm
	CPFM192A	CPFM192B	CPFM192C	CPFM192D	CPFM192P	1.00mm
	CPFM193A	CPFM193B	CPFM193C	CPFM193D	CPFM193P	2.00mm
Single-Use	CPFS191A	CPFS191B	CPFS191C	CPFS191D	CPFS191P	0.50mm
	CPFS192A	CPFS192B	CPFS192C	CPFS192D	CPFS192P	1.00mm
	CPFS193A	CPFS193B	CPFS193C	CPFS193D	CPFS193P	2.00mm

Box of 20 Units

C-TECH
BIPOLAR INSTRUMENTS

TITANIUM FORCEPS



UMAX: 0.5KVp



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM194A	CPFM194B	CPFM194C	CPFM194D	CPFM194P	0.50mm
	CPFM195A	CPFM195B	CPFM195C	CPFM195D	CPFM195P	1.00mm
	CPFM196A	CPFM196B	CPFM196C	CPFM196D	CPFM196P	2.00mm
Single-Use	CPFS194A	CPFS194B	CPFS194C	CPFS194D	CPFS194P	0.50mm
	CPFS195A	CPFS195B	CPFS195C	CPFS195D	CPFS195P	1.00mm
	CPFS196A	CPFS196B	CPFS196C	CPFS196D	CPFS196P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM197A	CPFM197B	CPFM197C	CPFM197D	CPFM197P	0.50mm
	CPFM198A	CPFM198B	CPFM198C	CPFM198D	CPFM198P	1.00mm
	CPFM199A	CPFM199B	CPFM199C	CPFM199D	CPFM199P	2.00mm
Single-Use	CPFS197A	CPFS197B	CPFS197C	CPFS197D	CPFS197P	0.50mm
	CPFS198A	CPFS198B	CPFS198C	CPFS198D	CPFS198P	1.00mm
	CPFS199A	CPFS199B	CPFS199C	CPFS199D	CPFS199P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM200A	CPFM200B	CPFM200C	CPFM200D	CPFM200P	0.50mm
	CPFM201A	CPFM201B	CPFM201C	CPFM201D	CPFM201P	1.00mm
	CPFM202A	CPFM202B	CPFM202C	CPFM202D	CPFM202P	2.00mm
Single-Use	CPFS200A	CPFS200B	CPFS200C	CPFS200D	CPFS200P	0.50mm
	CPFS201A	CPFS201B	CPFS201C	CPFS201D	CPFS201P	1.00mm
	CPFS202A	CPFS202B	CPFS202C	CPFS202D	CPFS202P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM203A	CPFM203B	CPFM203C	CPFM203D	CPFM203P	0.50mm
	CPFM204A	CPFM204B	CPFM204C	CPFM204D	CPFM204P	1.00mm
	CPFM205A	CPFM205B	CPFM205C	CPFM205D	CPFM205P	2.00mm
Single-Use	CPFS203A	CPFS203B	CPFS203C	CPFS203D	CPFS203P	0.50mm
	CPFS204A	CPFS204B	CPFS204C	CPFS204D	CPFS204P	1.00mm
	CPFS205A	CPFS205B	CPFS205C	CPFS205D	CPFS205P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM206A	CPFM206B	CPFM206C	CPFM206D	CPFM206P	0.50mm
	CPFM207A	CPFM207B	CPFM207C	CPFM207D	CPFM207P	1.00mm
	CPFM208A	CPFM208B	CPFM208C	CPFM208D	CPFM208P	2.00mm
Single-Use	CPFS206A	CPFS206B	CPFS206C	CPFS206D	CPFS206P	0.50mm
	CPFS207A	CPFS207B	CPFS207C	CPFS207D	CPFS207P	1.00mm
	CPFS208A	CPFS208B	CPFS208C	CPFS208D	CPFS208P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM209A	CPFM209B	CPFM209C	CPFM209D	CPFM209P	0.50mm
	CPFM210A	CPFM210B	CPFM210C	CPFM210D	CPFM210P	1.00mm
	CPFM211A	CPFM211B	CPFM211C	CPFM211D	CPFM211P	2.00mm
Single-Use	CPFS209A	CPFS209B	CPFS209C	CPFS209D	CPFS209P	0.50mm
	CPFS210A	CPFS210B	CPFS210C	CPFS210D	CPFS210P	1.00mm
	CPFS211A	CPFS211B	CPFS211C	CPFS211D	CPFS211P	2.00mm

Box of 20 Units



C-TECH
BIPOLAR INSTRUMENTS

ELITE NON-STICK FORECEPS



UMAX: 0.5KVP



102mm		115mm	
Multi-Use	CPFM212T	CPFM212V	0.50mm
	CPFM213T	CPFM213V	1.00mm
	CPFM214T	CPFM214V	2.00mm
Single-Use	CPFS212T	CPFS212V	0.50mm
	CPFS213T	CPFS213V	1.00mm
	CPFS214T	CPFS214V	2.00mm

Box of 20 Units



140mm		155mm		180mm		200mm		220mm		240mm		300mm	
Multi-Use	CPFM218X	CPFM218A	CPFM218B	CPFM218C	CPFM218D	CPFM218P	CPFM218F	0.50mm					
	CPFM219X	CPFM219A	CPFM219B	CPFM219C	CPFM219D	CPFM219P	CPFM219F	1.00mm					
	CPFM220X	CPFM220A	CPFM220B	CPFM220C	CPFM220D	CPFM220P	CPFM220F	2.00mm					
Single-Use	CPFS218X	CPFS218A	CPFS218B	CPFS218C	CPFS218D	CPFS218P	CPFS218F	0.50mm					
	CPFS219X	CPFS219A	CPFS219B	CPFS219C	CPFS219D	CPFS219P	CPFS219F	1.00mm					
	CPFS220X	CPFS220A	CPFS220B	CPFS220C	CPFS220D	CPFS220P	CPFS220F	2.00mm					

Box of 20 Units



102mm		115mm	
Multi-Use	CPFM215T	CPFM215V	0.50mm
	CPFM216T	CPFM216V	1.00mm
	CPFM217T	CPFM217V	2.00mm
Single-Use	CPFS215T	CPFS215V	0.50mm
	CPFS216T	CPFS216V	1.00mm
	CPFS217T	CPFS217V	2.00mm

Box of 20 Units



140mm		155mm		180mm		200mm		220mm		240mm		300mm	
Multi-Use	CPFM221X	CPFM221A	CPFM221B	CPFM221C	CPFM221D	CPFM221P	CPFM221F	0.50mm					
	CPFM222X	CPFM222A	CPFM222B	CPFM222C	CPFM222D	CPFM222P	CPFM222F	1.00mm					
	CPFM223X	CPFM223A	CPFM223B	CPFM223C	CPFM223D	CPFM223P	CPFM223F	2.00mm					
Single-Use	CPFS221X	CPFS221A	CPFS221B	CPFS221C	CPFS221D	CPFS221P	CPFS221F	0.50mm					
	CPFS222X	CPFS222A	CPFS222B	CPFS222C	CPFS222D	CPFS222P	CPFS222F	1.00mm					
	CPFS223X	CPFS223A	CPFS223B	CPFS223C	CPFS223D	CPFS223P	CPFS223F	2.00mm					

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

ELITE NON-STICK FORECEPS



UMAX: 0.5KVP



102mm		115mm	
Multi-Use	CPFM224T	CPFM224V	0.50mm
	CPFM225T	CPFM225V	1.00mm
	CPFM226T	CPFM226V	2.00mm
Single-Use	CPFS224T	CPFS224V	0.50mm
	CPFS225T	CPFS225V	1.00mm
	CPFS226T	CPFS226V	2.00mm

Box of 20 Units



140mm		155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM230X	CPFM230A	CPFM230B	CPFM230C	CPFM230D	CPFM230P	CPFM230F	0.50mm
	CPFM231X	CPFM231A	CPFM231B	CPFM231C	CPFM231D	CPFM231P	CPFM231F	1.00mm
	CPFM232X	CPFM232A	CPFM232B	CPFM232C	CPFM232D	CPFM232P	CPFM232F	2.00mm
Single-Use	CPFS230X	CPFS230A	CPFS230B	CPFS230C	CPFS230D	CPFS230P	CPFS230F	0.50mm
	CPFS231X	CPFS231A	CPFS231B	CPFS231C	CPFS231D	CPFS231P	CPFS231F	1.00mm
	CPFS232X	CPFS232A	CPFS232B	CPFS232C	CPFS232D	CPFS232P	CPFS232F	2.00mm

Box of 20 Units



102mm		115mm	
Multi-Use	CPFM227T	CPFM227V	0.50mm
	CPFM228T	CPFM228V	1.00mm
	CPFM229T	CPFM229V	2.00mm
Single-Use	CPFS227T	CPFS227V	0.50mm
	CPFS228T	CPFS228V	1.00mm
	CPFS229T	CPFS229V	2.00mm

Box of 20 Units



140mm		155mm	180mm	200mm	220mm	240mm	300mm	
Multi-Use	CPFM233X	CPFM233A	CPFM233B	CPFM233C	CPFM233D	CPFM233P	CPFM233F	0.50mm
	CPFM234X	CPFM234A	CPFM234B	CPFM234C	CPFM234D	CPFM234P	CPFM234F	1.00mm
	CPFM235X	CPFM235A	CPFM235B	CPFM235C	CPFM235D	CPFM235P	CPFM235F	2.00mm
Single-Use	CPFS233X	CPFS233A	CPFS233B	CPFS233C	CPFS233D	CPFS233P	CPFS233F	0.50mm
	CPFS234X	CPFS234A	CPFS234B	CPFS234C	CPFS234D	CPFS234P	CPFS234F	1.00mm
	CPFS235X	CPFS235A	CPFS235B	CPFS235C	CPFS235D	CPFS235P	CPFS235F	2.00mm

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

ELITE NON-STICK FORECEPS



UMAX: 0.5KVp



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM236A	CPFM236B	CPFM236C	CPFM236D	CPFM236P	0.50mm
	CPFM237A	CPFM237B	CPFM237C	CPFM237D	CPFM237P	1.00mm
	CPFM238A	CPFM238B	CPFM238C	CPFM238D	CPFM238P	2.00mm
Single-Use	CPFS236A	CPFS236B	CPFS236C	CPFS236D	CPFS236P	0.50mm
	CPFS237A	CPFS237B	CPFS237C	CPFS237D	CPFS237P	1.00mm
	CPFS238A	CPFS238B	CPFS238C	CPFS238D	CPFS238P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM239A	CPFM239B	CPFM239C	CPFM239D	CPFM239P	0.50mm
	CPFM240A	CPFM240B	CPFM240C	CPFM240D	CPFM240P	1.00mm
	CPFM241A	CPFM241B	CPFM241C	CPFM241D	CPFM241P	2.00mm
Single-Use	CPFS239A	CPFS239B	CPFS239C	CPFS239D	CPFS239P	0.50mm
	CPFS240A	CPFS240B	CPFS240C	CPFS240D	CPFS240P	1.00mm
	CPFS241A	CPFS241B	CPFS241C	CPFS241D	CPFS241P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM242A	CPFM242B	CPFM242C	CPFM242D	CPFM242P	0.50mm
	CPFM243A	CPFM243B	CPFM243C	CPFM243D	CPFM243P	1.00mm
	CPFM245A	CPFM245B	CPFM245C	CPFM245D	CPFM245P	2.00mm
Single-Use	CPFS242A	CPFS242B	CPFS242C	CPFS242D	CPFS242P	0.50mm
	CPFS243A	CPFS243B	CPFS243C	CPFS243D	CPFS243P	1.00mm
	CPFS245A	CPFS245B	CPFS245C	CPFS245D	CPFS245P	2.00mm

Box of 20 Units

C-TECH BIPOLAR INSTRUMENTS

ELITE NON-STICK FORCECEPS



UMAX: 0.5KVP



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM246A	CPFM246B	CPFM246C	CPFM246D	CPFM246P	0.50mm
	CPFM247A	CPFM247B	CPFM247C	CPFM247D	CPFM247P	1.00mm
	CPFM248A	CPFM248B	CPFM248C	CPFM248D	CPFM248P	2.00mm
Single-Use	CPFS246A	CPFS246B	CPFS246C	CPFS246D	CPFS246P	0.50mm
	CPFS247A	CPFS247B	CPFS247C	CPFS247D	CPFS247P	1.00mm
	CPFS248A	CPFS248B	CPFS248C	CPFS248D	CPFS248P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM249A	CPFM249B	CPFM249C	CPFM249D	CPFM249P	0.50mm
	CPFM250A	CPFM250B	CPFM250C	CPFM250D	CPFM250P	1.00mm
	CPFM251A	CPFM251B	CPFM251C	CPFM251D	CPFM251P	2.00mm
Single-Use	CPFS249A	CPFS249B	CPFS249C	CPFS249D	CPFS249P	0.50mm
	CPFS250A	CPFS250B	CPFS250C	CPFS250D	CPFS250P	1.00mm
	CPFS251A	CPFS251B	CPFS251C	CPFS251D	CPFS251P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM252A	CPFM252B	CPFM252C	CPFM252D	CPFM252P	0.50mm
	CPFM253A	CPFM253B	CPFM253C	CPFM253D	CPFM253P	1.00mm
	CPFM254A	CPFM254B	CPFM254C	CPFM254D	CPFM254P	2.00mm
Single-Use	CPFS252A	CPFS252B	CPFS252C	CPFS252D	CPFS252P	0.50mm
	CPFS253A	CPFS253B	CPFS253C	CPFS253D	CPFS253P	1.00mm
	CPFS254A	CPFS254B	CPFS254C	CPFS254D	CPFS254P	2.00mm

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

ELITE TC-TIP FORECEPS



UMAX: 0.5KVP



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM255A CPFM256A CPFM257A	CPFM255B CPFM256B CPFM257B	CPFM255C CPFM256C CPFM257C	CPFM255D CPFM256D CPFM257D	CPFM255P CPFM256P CPFM257P	0.50mm 1.00mm 2.00mm
Single-Use	CPFS255A CPFS256A CPFS257A	CPFS255B CPFS256B CPFS257B	CPFS255C CPFS256C CPFS257C	CPFS255D CPFS256D CPFS257D	CPFS255P CPFS256P CPFS257P	0.50mm 1.00mm 2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM261A CPFM262A CPFM263A	CPFM261B CPFM262B CPFM263B	CPFM261C CPFM262C CPFM263C	CPFM261D CPFM262D CPFM263D	CPFM261P CPFM262P CPFM263P	0.50mm 1.00mm 2.00mm
Single-Use	CPFS261A CPFS262A CPFS263A	CPFS261B CPFS262B CPFS263B	CPFS261C CPFS262C CPFS263C	CPFS261D CPFS262D CPFS263D	CPFS261P CPFS262P CPFS263P	0.50mm 1.00mm 2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM258A CPFM259A CPFM260A	CPFM258B CPFM259B CPFM260B	CPFM258C CPFM259C CPFM260C	CPFM258D CPFM259D CPFM260D	CPFM258P CPFM259P CPFM260P	0.50mm 1.00mm 2.00mm
Single-Use	CPFS258A CPFS259A CPFS260A	CPFS258B CPFS259B CPFS260B	CPFS258C CPFS259C CPFS260C	CPFS258D CPFS259D CPFS260D	CPFS258P CPFS259P CPFS260P	0.50mm 1.00mm 2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM264A CPFM265A CPFM266A	CPFM264B CPFM265B CPFM266B	CPFM264C CPFM265C CPFM266C	CPFM264D CPFM265D CPFM266D	CPFM264P CPFM265P CPFM266P	0.50mm 1.00mm 2.00mm
Single-Use	CPFS264A CPFS265A CPFS266A	CPFS264B CPFS265B CPFS266B	CPFS264C CPFS265C CPFS266C	CPFS264D CPFS265D CPFS266D	CPFS264P CPFS265P CPFS266P	0.50mm 1.00mm 2.00mm

Box of 20 Units



C-TECH BIPOLAR INSTRUMENTS

ELITE TC-TIP FORCEPS



UMAX: 0.5KVP



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM267A	CPFM267B	CPFM267C	CPFM267D	CPFM267P	0.50mm
	CPFM268A	CPFM268B	CPFM268C	CPFM268D	CPFM268P	1.00mm
	CPFM269A	CPFM269B	CPFM269C	CPFM269D	CPFM269P	2.00mm
Single-Use	CPFS267A	CPFS267B	CPFS267C	CPFS267D	CPFS267P	0.50mm
	CPFS268A	CPFS268B	CPFS268C	CPFS268D	CPFS268P	1.00mm
	CPFS269A	CPFS269B	CPFS269C	CPFS269D	CPFS269P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM273A	CPFM273B	CPFM273C	CPFM273D	CPFM273P	0.50mm
	CPFM274A	CPFM274B	CPFM274C	CPFM274D	CPFM274P	1.00mm
	CPFM275A	CPFM275B	CPFM275C	CPFM275D	CPFM275P	2.00mm
Single-Use	CPFS273A	CPFS273B	CPFS273C	CPFS273D	CPFS273P	0.50mm
	CPFS274A	CPFS274B	CPFS274C	CPFS274D	CPFS274P	1.00mm
	CPFS275A	CPFS275B	CPFS275C	CPFS275D	CPFS275P	2.00mm

Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM270A	CPFM270B	CPFM270C	CPFM270D	CPFM270P	0.50mm
	CPFM271A	CPFM271B	CPFM271C	CPFM271D	CPFM271P	1.00mm
	CPFM272A	CPFM272B	CPFM272C	CPFM272D	CPFM272P	2.00mm
Single-Use	CPFS270A	CPFS270B	CPFS270C	CPFS270D	CPFS270P	0.50mm
	CPFS271A	CPFS271B	CPFS271C	CPFS271D	CPFS271P	1.00mm
	CPFS272A	CPFS272B	CPFS272C	CPFS272D	CPFS272P	2.00mm


Box of 20 Units



	155mm	180mm	200mm	220mm	240mm	
Multi-Use	CPFM276A	CPFM276B	CPFM276C	CPFM276D	CPFM276P	0.50mm
	CPFM277A	CPFM277B	CPFM277C	CPFM277D	CPFM277P	1.00mm
	CPFM278A	CPFM278B	CPFM278C	CPFM278D	CPFM278P	2.00mm
Single-Use	CPFS276A	CPFS276B	CPFS276C	CPFS276D	CPFS276P	0.50mm
	CPFS277A	CPFS277B	CPFS277C	CPFS277D	CPFS277P	1.00mm
	CPFS278A	CPFS278B	CPFS278C	CPFS278D	CPFS278P	2.00mm

Box of 20 Units





03

DIATHERMY
ELECTROSURGICAL
CABLES AND PLUGS

PRECISE. INSULATED. RELIABLE.

DIATHERMY ELECTROSURGICAL CABLES AND PLUGS

BIPOLAR | FLAT, TWO POOL, TWO PIN, DUAL FITTING



UMAX: 0.5KVP



Instruments Side ▶▶									
▼▼ Generator Side									
									
 SIMENS: SUTTER BM50, ERBE T-SERIES: ERBE ACC/ICC 4MM	Multi-Use	CDCM101 - 3 Meter CDCM102 - 5 Meter	CDCM103 - 3 Meter CDCM104 - 5 Meter	CDCM105 - 3 Meter CDCM106 - 5 Meter	CDCM107 - 3 Meter CDCM108 - 5 Meter	CDCM109 - 3 Meter CDCM110 - 5 Meter	CDCM111 - 3 Meter CDCM112 - 5 Meter	CDCM113 - 3 Meter CDCM114 - 5 Meter	CDCM115 - 3 Meter CDCM116 - 5 Meter
	Single-Use	CDCS101 - 3 Meter CDCS102 - 5 Meter	CDCS103 - 3 Meter CDCS104 - 5 Meter	CDCS105 - 3 Meter CDCS106 - 5 Meter	CDCS107 - 3 Meter CDCS108 - 5 Meter	CDCS109 - 3 Meter CDCS110 - 5 Meter	CDCS111 - 3 Meter CDCS112 - 5 Meter	CDCS113 - 3 Meter CDCS114 - 5 Meter	CDCS115 - 3 Meter CDCS116 - 5 Meter
 ASSCULAP: WOLF: MARTIN BERCHTOLD: SUTTER BM780	Multi-Use	CDCM117 - 3 Meter CDCM118 - 5 Meter	CDCM119 - 3 Meter CDCM120 - 5 Meter	CDCM121 - 3 Meter CDCM122 - 5 Meter	CDCM123 - 3 Meter CDCM124 - 5 Meter	CDCM125 - 3 Meter CDCM126 - 5 Meter	CDCM127 - 3 Meter CDCM128 - 5 Meter	CDCM129 - 3 Meter CDCM130 - 5 Meter	CDCM131 - 3 Meter CDCM132 - 5 Meter
	Single-Use	CDCS117 - 3 Meter CDCS118 - 5 Meter	CDCS119 - 3 Meter CDCS120 - 5 Meter	CDCS121 - 3 Meter CDCS122 - 5 Meter	CDCS123 - 3 Meter CDCS124 - 5 Meter	CDCS125 - 3 Meter CDCS126 - 5 Meter	CDCS127 - 3 Meter CDCS128 - 5 Meter	CDCS129 - 3 Meter CDCS130 - 5 Meter	CDCS131 - 3 Meter CDCS132 - 5 Meter
 STANDARD INTERNATIONAL 4MM BANANA PLUG FITTING	Multi-Use	CDCM133 - 3 Meter CDCM134 - 5 Meter	CDCM135 - 3 Meter CDCM136 - 5 Meter	CDCM137 - 3 Meter CDCM138 - 5 Meter	CDCM139 - 3 Meter CDCM140 - 5 Meter	CDCM141 - 3 Meter CDCM142 - 5 Meter	CDCM143 - 3 Meter CDCM144 - 5 Meter	CDCM145 - 3 Meter CDCM146 - 5 Meter	CDCM147 - 3 Meter CDCM148 - 5 Meter
	Single-Use	CDCS133 - 3 Meter CDCS134 - 5 Meter	CDCS135 - 3 Meter CDCS136 - 5 Meter	CDCS137 - 3 Meter CDCS138 - 5 Meter	CDCS139 - 3 Meter CDCS140 - 5 Meter	CDCS141 - 3 Meter CDCS142 - 5 Meter	CDCS143 - 3 Meter CDCS144 - 5 Meter	CDCS145 - 3 Meter CDCS146 - 5 Meter	CDCS147 - 3 Meter CDCS148 - 5 Meter
 VALLEYLAB: LAMIDEY: EMC BOWA, ETHICON	Multi-Use	CDCM149 - 3 Meter CDCM150 - 5 Meter	CDCM151 - 3 Meter CDCM152 - 5 Meter	CDCM153 - 3 Meter CDCM154 - 5 Meter	CDCM155 - 3 Meter CDCM156 - 5 Meter	CDCM157 - 3 Meter CDCM158 - 5 Meter	CDCM159 - 3 Meter CDCM160 - 5 Meter	CDCM161 - 3 Meter CDCM162 - 5 Meter	CDCM163 - 3 Meter CDCM164 - 5 Meter
	Single-Use	CDCS149 - 3 Meter CDCS150 - 5 Meter	CDCS151 - 3 Meter CDCS152 - 5 Meter	CDCS153 - 3 Meter CDCS154 - 5 Meter	CDCS155 - 3 Meter CDCS156 - 5 Meter	CDCS157 - 3 Meter CDCS158 - 5 Meter	CDCS159 - 3 Meter CDCS160 - 5 Meter	CDCS161 - 3 Meter CDCS162 - 5 Meter	CDCS163 - 3 Meter CDCS164 - 5 Meter
 EMC: DOLLEY, ELLMAN	Multi-Use	CDCM165 - 3 Meter CDCM166 - 5 Meter	CDCM167 - 3 Meter CDCM168 - 5 Meter	CDCM169 - 3 Meter CDCM170 - 5 Meter	CDCM171 - 3 Meter CDCM172 - 5 Meter	CDCM173 - 3 Meter CDCM174 - 5 Meter	CDCM175 - 3 Meter CDCM176 - 5 Meter	CDCM177 - 3 Meter CDCM178 - 5 Meter	CDCM179 - 3 Meter CDCM180 - 5 Meter
	Single-Use	CDCS165 - 3 Meter CDCS166 - 5 Meter	CDCS167 - 3 Meter CDCS168 - 5 Meter	CDCS169 - 3 Meter CDCS170 - 5 Meter	CDCS171 - 3 Meter CDCS172 - 5 Meter	CDCS173 - 3 Meter CDCS174 - 5 Meter	CDCS175 - 3 Meter CDCS176 - 5 Meter	CDCS177 - 3 Meter CDCS178 - 5 Meter	CDCS179 - 3 Meter CDCS180 - 5 Meter
 VALLEYLAB (AUTOBIPOLAR MODELS ONLY)	Multi-Use	CDCM181 - 3 Meter CDCM182 - 5 Meter	CDCM183 - 3 Meter CDCM184 - 5 Meter	CDCM185 - 3 Meter CDCM186 - 5 Meter	CDCM187 - 3 Meter CDCM188 - 5 Meter	CDCM189 - 3 Meter CDCM190 - 5 Meter	CDCM191 - 3 Meter CDCM192 - 5 Meter	CDCM193 - 3 Meter CDCM194 - 5 Meter	CDCM195 - 3 Meter CDCM196 - 5 Meter
	Single-Use	CDCS181 - 3 Meter CDCS182 - 5 Meter	CDCS183 - 3 Meter CDCS184 - 5 Meter	CDCS185 - 3 Meter CDCS186 - 5 Meter	CDCS187 - 3 Meter CDCS188 - 5 Meter	CDCS189 - 3 Meter CDCS190 - 5 Meter	CDCS191 - 3 Meter CDCS192 - 5 Meter	CDCS193 - 3 Meter CDCS194 - 5 Meter	CDCS195 - 3 Meter CDCS196 - 5 Meter

DIATHERMY ELECTROSURGICAL CABLES AND PLUGS

BIPOLAR | FLAT, TWO POOL, TWO PIN, DUAL FITTING



AUTOCALVABLE
REUSABLE



MAX TEMPERATURE
RESISTANCE



MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp



Instruments Side ▶	Generator Side ▼▼			
	4mm Female Plugs	4mm Male Plugs Fitting	2.4mm Female Ethicon	Block Connector Fitting
 SIMENS: SUTTER BM50, ERBE T-SERIES: ERBE ACC/ICC 4MM	CDCM236 - 3 Meter CDCM237 - 5 Meter	CDCM238 - 3 Meter CDCM239 - 5 Meter	CDCM240 - 3 Meter CDCM241 - 5 Meter	CDCM242 - 3 Meter CDCM243 - 5 Meter
	CDCS236 - 3 Meter CDCS237 - 5 Meter	CDCS238 - 3 Meter CDCS239 - 5 Meter	CDCS240 - 3 Meter CDCS241 - 5 Meter	CDCS242 - 3 Meter CDCS243 - 5 Meter
 ASSCULAP: WOLF: MARTIN BERCHTOLD: SUTTER BM780	CDCM244 - 3 Meter CDCM245 - 5 Meter	CDCM246 - 3 Meter CDCM247 - 5 Meter	CDCM248 - 3 Meter CDCM249 - 5 Meter	CDCM250 - 3 Meter CDCM251 - 5 Meter
	CDCS244 - 3 Meter CDCS245 - 5 Meter	CDCS246 - 3 Meter CDCS247 - 5 Meter	CDCS248 - 3 Meter CDCS249 - 5 Meter	CDCS250 - 3 Meter CDCS251 - 5 Meter
 STANDARD INTERNATIONAL 4MM BANANA PLUG FITTING	CDCM252 - 3 Meter CDCM253 - 5 Meter	CDCM254 - 3 Meter CDCM255 - 5 Meter	CDCM256 - 3 Meter CDCM257 - 5 Meter	CDCM258 - 3 Meter CDCM259 - 5 Meter
	CDCS252 - 3 Meter CDCS253 - 5 Meter	CDCS254 - 3 Meter CDCS255 - 5 Meter	CDCS256 - 3 Meter CDCS257 - 5 Meter	CDCS258 - 3 Meter CDCS259 - 5 Meter
 VALLEYLAB: LAMIDEY: EMC BOWA, ETHICON	CDCM260 - 3 Meter CDCM261 - 5 Meter	CDCM262 - 3 Meter CDCM263 - 5 Meter	CDCM264 - 3 Meter CDCM265 - 5 Meter	CDCM266 - 3 Meter CDCM267 - 5 Meter
	CDCS260 - 3 Meter CDCS261 - 5 Meter	CDCS262 - 3 Meter CDCS263 - 5 Meter	CDCS264 - 3 Meter CDCS265 - 5 Meter	CDCS266 - 3 Meter CDCS267 - 5 Meter
 EMC: DOLLEY, ELLMAN	CDCM268 - 3 Meter CDCM269 - 5 Meter	CDCM270 - 3 Meter CDCM271 - 5 Meter	CDCM272 - 3 Meter CDCM273 - 5 Meter	CDCM274 - 3 Meter CDCM275 - 5 Meter
	CDCS268 - 3 Meter CDCS269 - 5 Meter	CDCS270 - 3 Meter CDCS271 - 5 Meter	CDCS272 - 3 Meter CDCS273 - 5 Meter	CDCS274 - 3 Meter CDCS275 - 5 Meter
 VALLEYLAB (AUTOBIPOLAR MODELS ONLY)	CDCM276 - 3 Meter CDCM277 - 5 Meter	CDCM278 - 3 Meter CDCM279 - 5 Meter	CDCM280 - 3 Meter CDCM281 - 5 Meter	CDCM282 - 3 Meter CDCM283 - 5 Meter
	CDCS276 - 3 Meter CDCS277 - 5 Meter	CDCS278 - 3 Meter CDCS279 - 5 Meter	CDCS280 - 3 Meter CDCS281 - 5 Meter	CDCS282 - 3 Meter CDCS283 - 5 Meter

▼▼ Generator Side

▼▼ Instruments Side

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM284 - 3 Meter
CDCM285 - 5 Meter

Single-Use

CDCS284 - 3 Meter
CDCS285 - 5 Meter

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM286 - 3 Meter
CDCM287 - 5 Meter

Single-Use

CDCS286 - 3 Meter
CDCS287 - 5 Meter

Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Conmed 30mm



Multi-Use

CDCM288 - 3 Meter
CDCM289 - 5 Meter

Single-Use

CDCS288 - 3 Meter
CDCS289 - 5 Meter

4mm banana to 6mm Alsa Excel connector adapter for Alsa Excelgenerator



Multi-Use

CDCM290 - 3 Meter
CDCM291 - 5 Meter

Single-Use

CDCS290 - 3 Meter
CDCS291 - 5 Meter

DIATHERMY ELECTROSURGICAL CABLES AND PLUGS

MONOPOLAR | INTERNATIONAL FITTING PLUG, ADAPTER



UMAX: 0.5KVp



Instruments Side ▶▶ ▼▼ Generator Side	   			
	Multi-Use	Single-Use	Multi-Use	Single-Use
 MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CDCM197 - 3 Meter CDCM198 - 5 Meter	CDCM199 - 3 Meter CDCM200 - 5 Meter	CDCM201 - 3 Meter CDCM202 - 5 Meter	CDCM203 - 3 Meter CDCM204 - 5 Meter
	CDCS197 - 3 Meter CDCS198 - 5 Meter	CDCS199 - 3 Meter CDCS200 - 5 Meter	CDCS201 - 3 Meter CDCS202 - 5 Meter	CDCS203 - 3 Meter CDCS204 - 5 Meter
 ERBE ACC/ICC	CDCM205 - 3 Meter CDCM206 - 5 Meter	CDCM207 - 3 Meter CDCM208 - 5 Meter	CDCM209 - 3 Meter CDCM210 - 5 Meter	CDCM211 - 3 Meter CDCM212 - 5 Meter
	CDCS205 - 3 Meter CDCS206 - 5 Meter	CDCS207 - 3 Meter CDCS208 - 5 Meter	CDCS209 - 3 Meter CDCS210 - 5 Meter	CDCS211 - 3 Meter CDCS212 - 5 Meter
 VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CDCM213 - 3 Meter CDCM214 - 5 Meter	CDCM215 - 3 Meter CDCM216 - 5 Meter	CDCM217 - 3 Meter CDCM218 - 5 Meter	CDCM219 - 3 Meter CDCM220 - 5 Meter
	CDCS213 - 3 Meter CDCS214 - 5 Meter	CDCS215 - 3 Meter CDCS216 - 5 Meter	CDCS217 - 3 Meter CDCS218 - 5 Meter	CDCS219 - 3 Meter CDCS220 - 5 Meter
 MARTIN: BERCHTOLD:	CDCM221 - 3 Meter CDCM222 - 5 Meter	CDCM223 - 3 Meter CDCM224 - 5 Meter	CDCM225 - 3 Meter CDCM226 - 5 Meter	CDCM227 - 3 Meter CDCM228 - 5 Meter
	CDCS221 - 3 Meter CDCS222 - 5 Meter	CDCS223 - 3 Meter CDCS224 - 5 Meter	CDCS225 - 3 Meter CDCS226 - 5 Meter	CDCS227 - 3 Meter CDCS228 - 5 Meter
 INTERNATIONAL 3 PIN PLUG DIATHERMY	CDCM229 - 3 Meter CDCM230 - 5 Meter	N/A	N/A	CDCM231 - 3 Meter CDCM232 - 5 Meter
	CDCS229 - 3 Meter CDCS230 - 5 Meter	N/A	N/A	CDCS231 - 3 Meter CDCS232 - 5 Meter

Connecting Finger switch handles with 3 pin plugs, for ERBE T-Series generators.



Connecting Finger switch handles with 3 pin plugs, for Martin, Integra generators





04

NON-CONDUCTIVE GYNE INSTRUMENTS

PRECISE. INSULATED. RELIABLE.

NON CONDUCTIVE GYNE INSTRUMENTS

FORCEPS, SPECULUM, RETRACTORS, HOOKS



UMAX: 0.5KVp



Pederson Vaginal Speculum
with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM101	SMALL	CGIS101
CGIM102	MEDIUM	CGIS102
CGIM103	LARGE	CGIS103

BOX OF 1 UNITS



Graves Vaginal Speculum
with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM104	SMALL	CGIS104
CGIM105	MEDIUM	CGIS105
CGIM106	LARGE	CGIS106

BOX OF 1 UNITS



Cusco Speculum
with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM107	SMALL	CGIS107
CGIM108	MEDIUM	CGIS108
CGIM109	LARGE	CGIS109


BOX OF 1 UNITS



Collin Speculum
with Smoke Evacuation Tube

Multi-Use		Single-Use
CGIM110	SMALL	CGIS110
CGIM111	MEDIUM	CGIS111
CGIM112	LARGE	CGIS112

BOX OF 1 UNITS



Lateral Vaginal Retractor

Multi-Use		Single-Use
CGIM113	SMALL	CGIS113
CGIM114	MEDIUM	CGIS114
CGIM115	LARGE	CGIS115

BOX OF 1 UNITS



Kogan Endocervical

Multi-Use		Single-Use
CGIM116	SMALL	CGIS116
CGIM117	MEDIUM	CGIS117
CGIM118	LARGE	CGIS118

BOX OF 1 UNITS




NON CONDUCTIVE GYNE INSTRUMENTS

Forceps, Speculum, Retractors, Hooks



UMAX: 0.5KVp



NON CONDUCTIVE GYNE INSTRUMENTS

Forceps, Speculum, Retractors, Hooks



UMAX: 0.5KVP



Dressing Forcep

Multi-Use		Single-Use
CGIM137	155mm	CGIS137
CGIM138	180mm	CGIS138
CGIM139	210mm	CGIS139
CGIM140	250mm	CGIS140


BOX OF 20 UNITS



Tissue Forcep

Multi-Use		Single-Use
CGIM141	155mm	CGIS141
CGIM142	180mm	CGIS142
CGIM143	210mm	CGIS143
CGIM144	250mm	CGIS144


BOX OF 20 UNITS



Iris Hook

Multi-Use		Single-Use
CGIM145	155mm	CGIS145
CGIM146	180mm	CGIS146
CGIM147	210mm	CGIS147
CGIM148	250mm	CGIS148


BOX OF 20 UNITS



Stragih Hook

Multi-Use		Single-Use
CGIM149	155mm	CGIS149
CGIM150	180mm	CGIS150
CGIM151	210mm	CGIS151


BOX OF 20 UNITS



Stragih Hook 2

Multi-Use		Single-Use
CGIM153	155mm	CGIS153
CGIM154	180mm	CGIS154
CGIM155	210mm	CGIS155

BOX OF 20 UNITS



Stragih Hook 3

Multi-Use		Single-Use
CGIM157	155mm	CGIS157
CGIM158	180mm	CGIS158
CGIM159	210mm	CGIS159

BOX OF 20 UNITS





05

DIATHERMY CAUTERY AND ELECTRODES

PRECISE. INSULATED. RELIABLE.

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

HAND CONTROL CAUTERY PENCIL



UMAX: 0.5KVP



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use		Single-Use
CDPM101	3M	CDPS101
CDPM102	5M	CDPS102



Martin, ERBE T-A Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM103	3M	CDPS103
CDPM104	5M	CDPS104



Valleylab, Bowa, ERBE ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM105	3M	CDPS105
CDPM106	5M	CDPS106



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM107	3M	CDPS107
CDPM108	5M	CDPS108



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use		Single-Use
CDPM117	3M	CDPS117
CDPM118	5M	CDPS118



Martin, ERBE T-A Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM119	3M	CDPS119
CDPM120	5M	CDPS120



Valleylab, Bowa, ERBE ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM121	3M	CDPS121
CDPM122	5M	CDPS122



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM123	3M	CDPS123
CDPM124	5M	CDPS124



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use		Single-Use
CDPM109	3M	CDPS109
CDPM110	5M	CDPS110



Martin, ERBE T-A Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM111	3M	CDPS111
CDPM112	5M	CDPS112



Valleylab, Bowa, ERBE ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM113	3M	CDPS113
CDPM114	5M	CDPS114



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM115	3M	CDPS115
CDPM116	5M	CDPS116



BOX OF 20 UNITS



4mm ERBE T-Series

Multi-Use		Single-Use
CDPM125	3M	CDPS125
CDPM126	5M	CDPS126



Martin, ERBE T-A Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM127	3M	CDPS127
CDPM128	5M	CDPS128



Valleylab, Bowa, ERBE ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM129	3M	CDPS129
CDPM130	5M	CDPS130



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM131	3M	CDPS131
CDPM132	5M	CDPS132

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

FOOT CONTROL CAUTERY PENCIL



UMAX: 0.5KVp



BOX OF 20 UNITS



4mm International, Bowa, ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM133	3M	CDPS133
CDPM134	5M	CDPS134



8mm ERBE T-Series
Valleylab, ERBE Int'l

Multi-Use		Single-Use
CDPM135	3M	CDPS135
CDPM136	5M	CDPS136



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM137	3M	CDPS137
CDPM138	5M	CDPS138



Martin, ERBE T-Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM139	3M	CDPS139
CDPM140	5M	CDPS140



BOX OF 20 UNITS



4mm International, Bowa, ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM149	3M	CDPS149
CDPM150	5M	CDPS150



8mm ERBE T-Series
Valleylab, ERBE Int'l

Multi-Use		Single-Use
CDPM151	3M	CDPS151
CDPM152	5M	CDPS152



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM153	3M	CDPS153
CDPM154	5M	CDPS154



Martin, ERBE T-Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM155	3M	CDPS155
CDPM156	5M	CDPS156



BOX OF 20 UNITS



4mm International, Bowa, ICC
ERBE Vio, Eschmann, Conmed

Multi-Use		Single-Use
CDPM141	3M	CDPS141
CDPM142	5M	CDPS142



8mm ERBE T-Series
Valleylab, ERBE Int'l

Multi-Use		Single-Use
CDPM143	3M	CDPS143
CDPM144	5M	CDPS144



5mm ERBE ACC Series Standard
ICC Series Vio, Karl Storz

Multi-Use		Single-Use
CDPM145	3M	CDPS145
CDPM146	5M	CDPS146



Martin, ERBE T-Series, Integra
Berchtold, Aesculap

Multi-Use		Single-Use
CDPM147	3M	CDPS147
CDPM148	5M	CDPS148

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

FOOT CONTROL CAUTERY PENCIL



UMAX: 0.5KVp



BOX OF 20 UNITS



INTERNATIONAL 3 PIN PLUG DIATHERMY



BOX OF 20 UNITS



INTERNATIONAL 3 PIN PLUG DIATHERMY



BOX OF 20 UNITS



INTERNATIONAL 3 PIN PLUG DIATHERMY



DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH NON-STICK ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS Multi-Use CDEM190	Single-Use CDES190		4mm Shaft		BOX OF 50 UNITS Multi-Use CDEM191	Single-Use CDES191
Multi-Use CDEM191	Single-Use CDES192		4mm Shaft		Multi-Use CDEM193	Single-Use CDES193
Multi-Use CDEM194	Single-Use CDES194		4mm Shaft		Multi-Use CDEM195	Single-Use CDES195
Multi-Use CDEM196	Single-Use CDES196		4mm Shaft		Multi-Use CDEM197	Single-Use CDES197
Multi-Use CDEM198	Single-Use CDES198		4mm Shaft		Multi-Use CDEM199	Single-Use CDES199
Multi-Use CDEM200	Single-Use CDES200		2.4mm Shaft		Multi-Use CDEM201	Single-Use CDES201
Multi-Use CDEM202	Single-Use CDES202		2.4mm Shaft		Multi-Use CDEM203	Single-Use CDES203

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH NON-STICK ELECTRODES



UMAX: 0.5KVP







BOX OF 50 UNITS Multi-Use CDEM204	Single-Use CDES204	 Wire Loop electrodes, Straight	5 mm	2.4mm Shaft	 Wire Loop electrodes, Straight	10 mm	BOX OF 50 UNITS Multi-Use CDEM205	Single-Use CDES205
Multi-Use CDEM206	Single-Use CDES206	 Wire Loop electrodes, Straight	15 mm	2.4mm Shaft	 Ball Electrodes, Straight	2 mm	Multi-Use CDEM207	Single-Use CDES207
Multi-Use CDEM208	Single-Use CDES208	 Ball Electrodes, Straight	4 mm	2.4mm Shaft	 Ball Electrodes, Straight	6 mm	Multi-Use CDEM209	Single-Use CDES1209
Multi-Use CDEM210	Single-Use CDES210	 Needle electrodes, Straight	0.5 mm	2.38mm Shaft	 Blade Electrodes, Straight	2.5x20mm	Multi-Use CDEM211	Single-Use CDES211
Multi-Use CDEM212	Single-Use CDES212	 Needle electrodes, Straight	0.3mm	2.38mm Shaft	 Needle Electrodes, Straight	0.5 mm	Multi-Use CDEM213	Single-Use CDES213
Multi-Use CDEM214	Single-Use CDES214	 Needle electrodes, Straight	0.3 mm	2.38mm Shaft	 Needle electrodes, Straight	0.5 mm	Multi-Use CDEM215	Single-Use CDES215
Multi-Use CDEM216	Single-Use CDES216	 Needle electrodes, Straight	0.5 mm	2.38mm Shaft	 Ball Electrodes, Straight	2 mm	Multi-Use CDEM217	Single-Use CDES217

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES
C-TECH NON-STICK ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM218 Single-Use CDES218		 Needle Electrodes, Insulated 45° Angled	1.6mm Shaft	 Needle Electrodes, Insulated 45° Angled	BOX OF 50 UNITS Multi-Use CDEM219 Single-Use CDES219
Multi-Use CDEM220 Single-Use CDES220		 Needle Electrodes, Insulated 45° Angled	1.6mm Shaft	 Needle Electrodes, Insulated 45° Angled	Multi-Use CDEM221 Single-Use CDES221
Multi-Use CDEM222 Single-Use CNES222		 Needle Electrodes, Insulated 90° Angled	1.6mm Shaft	 Wire Loop electrodes, Straight Triangle	Multi-Use CDEM223 Single-Use CDES223
Multi-Use CDEM224 Single-Use CDES224		 Wire Loop electrodes	1.6mm Shaft	 Wire Loop electrodes	Multi-Use CDEM225 Single-Use CDES225
Multi-Use CDEM226 Single-Use CDES226		 Wire Loop electrodes	1.6mm Shaft		

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS Multi-Use CDEM232	Single-Use CDES232		Knife Electrodes	4mm Shaft		BOX OF 50 UNITS Multi-Use CDEM233	Single-Use CDES233
Multi-Use CDEM234	Single-Use CDES234		Lancet-Cutting Electrodes	4mm Shaft		Multi-Use CDEM235	Single-Use CDES235
Multi-Use CDEM236	Single-Use CDES236		Spatula Electrodes, Oval Shape	4mm Shaft		Multi-Use CDEM237	Single-Use CDES237
Multi-Use CDEM238	Single-Use CDES238		Knife Electrodes	4mm Shaft		Multi-Use CDEM239	Single-Use CDES239
Multi-Use CDEM240	Single-Use CDES240		Knife Electrode, Lancet	4mm Shaft		Multi-Use CDEM241	Single-Use CDES241
Multi-Use CDEM242	Single-Use CDES242		Needle Electrodes, Insulated 10MM WL	4mm Shaft		Multi-Use CDEM243	Single-Use CDES243
Multi-Use CDEM244	Single-Use CDES244		Needle Electrode	4mm Shaft		Multi-Use CDEM245	Single-Use CDES245

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS Multi-Use CDEM246	Single-Use CDES246		4mm Shaft		BOX OF 50 UNITS Multi-Use CDEM247	Single-Use CDES247
Multi-Use CDEM248	Single-Use CDES248		4mm Shaft		Multi-Use CDEM249	Single-Use CDES249
Multi-Use CDEM250	Single-Use CDES250		4mm Shaft		Multi-Use CDEM251	Single-Use CDES251
Multi-Use CDEM252	Single-Use CDES252		4mm Shaft		Multi-Use CDEM253	Single-Use CDES253
Multi-Use CDEM254	Single-Use CDES254		4mm Shaft		Multi-Use CDEM255	Single-Use CDES255
Multi-Use CDEM256	Single-Use CDES256		4mm Shaft		Multi-Use CDEM257	Single-Use CDES257
Multi-Use CDEM258	Single-Use CDES258		4mm Shaft		Multi-Use CDEM259	Single-Use CDES259

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS Multi-Use CDEM260	Single-Use CDES260	 Wire Loop Electrode, Triangular	4mm Shaft	 Wire Loop Electrode, Triangular	BOX OF 50 UNITS Multi-Use CDEM261	Single-Use CDES261
Multi-Use CDEM262	Single-Use CDES262	 Wire Loop Electrode, Rhombic	4mm Shaft	 Wire Loop Electrode, Rhombic	Multi-Use CDEM263	Single-Use CDES263
Multi-Use CDEM264	Single-Use CDES264	 Wire Loop Electrode, Circular	4mm Shaft	 Wire Loop Electrode, Circular	Multi-Use CDEM265	Single-Use CDES265
Multi-Use CDEM266	Single-Use CDES266	 Needle Electrode	4mm Shaft	 Ball Electrode	Multi-Use CDEM267	Single-Use CDES267
Multi-Use CDEM268	Single-Use CDES268	 Ball Electrode	4mm Shaft	 Ball Electrode	Multi-Use CDEM269	Single-Use CDES269
Multi-Use CDEM270	Single-Use CDES270	 Ball Electrode	4mm Shaft	 Spatula Electrode	Multi-Use CDEM271	Single-Use CDES271

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR STEEL ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM272	Single-Use CDES272	 2mm Ball electrode	4mm Shaft	 4mm Ball electrode	BOX OF 50 UNITS Multi-Use CDEM273	Single-Use CDES273
Multi-Use CDEM274	Single-Use CDES274	 6mm Ball electrode	4mm Shaft	 5mm Wire loop electrode	Multi-Use CDEM275	Single-Use CDES275
Multi-Use CDEM276	Single-Use CDES276	 8mm Wire loop electrode	4mm Shaft	 10mm Wire loop electrode	Multi-Use CDEM277	Single-Use CDES277
Multi-Use CDEM278	Single-Use CDES278	 5mm Ribbon loop electrode	4mm Shaft	 10mm 15mm Ribbon loop electrode	Multi-Use CDEM279	Single-Use CDES279
Multi-Use CDEM280	Single-Use CDES280	 15mm Ribbon loop electrode	4mm Shaft	 1.6mm, 19mm Lancet electrode, Stright	Multi-Use CDEM281	Single-Use CDES281
Multi-Use CDEM282	Single-Use CDES282	 1.6x19mm LANCET ELECTRODE ANGULAR	4mm Shaft	 3.3X25mm Knife electrode	Multi-Use CDEM283	Single-Use CDES283
Multi-Use CDEM284	Single-Use CDES284	 3.3X24mm Knife electrode	4mm Shaft	 22mm, 30mm Needle electrode	Multi-Use CDEM285	Single-Use CDES285

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM286	Single-Use CDES286		4mm	2.4mm Shaft		2mm	BOX OF 50 UNITS Multi-Use CDEM287	Single-Use CDES287
Multi-Use CDEM288	Single-Use CDES288		4mm	2.4mm Shaft		6mm	Multi-Use CDEM289	Single-Use CDES289
Multi-Use CDEM290	Single-Use CDES290		2mm	2.4mm Shaft		4mm	Multi-Use CDEM291	Single-Use CDES291
Multi-Use CDEM292	Single-Use CDES292		6mm	2.4mm Shaft		2mm	Multi-Use CDEM293	Single-Use CDES293
Multi-Use CDEM294	Single-Use CDES294		2.5X20mm	2.4mm Shaft		2.5X20mm	Multi-Use CDEM295	Single-Use CDES295
Multi-Use CDEM296	Single-Use CDES296		0.8X20mm	2.4mm Shaft		0.5X20mm	Multi-Use CDEM297	Single-Use CDES297
Multi-Use CDEM298	Single-Use CDES298		0.5X20mm	2.4mm Shaft		0.5X20mm	Multi-Use CDEM299	Single-Use CDES299

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM300	Single-Use CDES300	 Needle electrode with ceramic isolation & Tungsten needle	2.4mm Shaft	 Needle electrode	BOX OF 50 UNITS Multi-Use CDEM301	Single-Use CDES301
Multi-Use CDEM302	Single-Use CDES302	 Wire loop electrode , Triangular	2.4mm Shaft	 Wire loop electrode , Rhombic	Multi-Use CDEM303	Single-Use CDES303
Multi-Use CDEM304	Single-Use CDES304	 Spatula electrode	2.4mm Shaft	 Wire loop electrode , Triangular	Multi-Use CDEM305	Single-Use CDES305
Multi-Use CDEM306	Single-Use CDES306	 Wire loop electrode	2.4mm Shaft	 Wire loop electrode	Multi-Use CDEM307	Single-Use CDES307
Multi-Use CDEM308	Single-Use CDES308	 Ribbon loop electrode	2.4mm Shaft	 Ribbon loop electrode	Multi-Use CDEM309	Single-Use CDES309
Multi-Use CDEM310	Single-Use CDES310	 Ribbon loop electrode	2.4mm Shaft	 Ribbon loop electrode	Multi-Use CDEM311	Single-Use CDES311
Multi-Use CDEM312	Single-Use CDES312	 Ball electrode	2.4mm Shaft	 Ball electrode	Multi-Use CDEM313	Single-Use CDES313

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM314	Single-Use CDES314		Electrode coagulating and fluguration	1.6mm Shaft		Wire loop electrode 3mm	BOX OF 50 UNITS Multi-Use CDEM315	Single-Use CDES315
Multi-Use CDEM316	Single-Use CDES316		Wire loop electrode	1.6mm Shaft		Wire loop electrode	Multi-Use CDEM317	Single-Use CDES317
Multi-Use CDEM318	Single-Use CDES318		Needle electrode wired WL 5mm	1.6mm Shaft		Needle electrode wired, Angled 45° WL 5mm	Multi-Use CDEM319	Single-Use CDES319
Multi-Use CDEM320	Single-Use CDES320		Coagulation Needle electrode, Angled 45° WL 5mm	1.6mm Shaft		Ball electrode, Angled 45°	Multi-Use CDEM321	Single-Use CDES321
Multi-Use CDEM322	Single-Use CDES322		Ball electrode	1.6mm Shaft		Ball electrode	Multi-Use CDEM323	Single-Use CDES323
Multi-Use CDEM324	Single-Use CDES324		Ball electrode	1.6mm Shaft		Ball electrode	Multi-Use CDEM325	Single-Use CDES325
Multi-Use CDEM326	Single-Use CDES326		Ball electrode, Angled 45°	2.38mm Shaft		Wire loop electrode, Triangular	Multi-Use CDEM327	Single-Use CDES327

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVP



Multi-Use CDEM328	Single-Use CNES328	 Wire loop electrode, Triangular	2.38mm Shaft	 Wire loop electrode	Multi-Use CDEM329	Single-Use CNES329
Multi-Use CDEM330	Single-Use CNES330	 Wire loop electrode, Rohimbic	2.38mm Shaft	 Wire loop electrode	Multi-Use CDEM331	Single-Use CNES331
Multi-Use CDEM332	Single-Use CNES332	 Wire loop electrode, Rohimbic	2.38mm Shaft	 Ball electrode	Multi-Use CDEM333	Single-Use CNES333
Multi-Use CDEM334	Single-Use CNES334	 Spatula electrode	2.38mm Shaft	 Needle electrodes, Straight	Multi-Use CDEM335	Single-Use CNES335

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM336	Single-Use CDES336		
		Wire loop electrode	
Multi-Use CDEM337	Single-Use CDES337		
		Wire loop electrode	
Multi-Use CDEM338	Single-Use CDES338		
		Wire loop electrode	
Multi-Use CDEM339	Single-Use CDES339		
		Wire loop electrode	
Multi-Use CDEM340	Single-Use CDES340		
		Ball electrode	
Multi-Use CDEM341	Single-Use CDES341		
		Ball electrode	
Multi-Use CDEM342	Single-Use CNES342		
		Needle electrode	

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS		140mm		
Multi-Use CDEM343	Single-Use CDES343		4mm Shaft	Needle electrode, Angled
		140mm		0.8X12mm
Multi-Use CDEM344	Single-Use CDES344		4mm Shaft	Knife electrode
		140mm		2.4X10mm
Multi-Use CDEM345	Single-Use CDES345			
		140mm		3.3X25mm
Multi-Use CDEM346	Single-Use CDES346			
		140mm		3.3X20mm
Multi-Use CDEM347	Single-Use CDES347		4mm Shaft	Knife Electrode
		140mm		2.4X20mm
Multi-Use CDEM348	Single-Use CDES348			
		140mm		16X20mm
Multi-Use CDEM349	Single-Use CNES349		4mm Shaft	Lancet Electrode

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS		140mm	2.2X20mm
Multi-Use CDEM350	Single-Use CDES350		
		4mm Shaft	
		100mm	2X16mm
Multi-Use CDEM351	Single-Use CDES351		
			electrode
		160mm	10X10mm
Multi-Use CDEM352	Single-Use CDES352		
		160mm	15X15mm
Multi-Use CDEM353	Single-Use CDES353		
		160mm	15X20mm
Multi-Use CDEM354	Single-Use CDES354		
		4mm Shaft	Loop Electrode
		160mm	20X20mm
Multi-Use CDEM355	Single-Use CDES355		
			BOX OF 50 UNITS
		160mm	25mm
Multi-Use CDEM356	Single-Use CDES356		

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS		160mm	10mm	BOX OF 50 UNITS
Multi-Use CDEM357	Single-Use CDES357			
		4mm Shaft	Loop Electrode	
Multi-Use CDEM358	Single-Use CDES358	160mm	15mm	
		4mm Shaft	Loop Electrode	
Multi-Use CDEM39	Single-Use CDES359	160mm	20mm	
Multi-Use CDEM3560	Single-Use CDES360	160mm	25mm	
Multi-Use CDEM361	Single-Use CDES3561	160mm	30mm	
		4mm Shaft	Loop Electrode	

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS		160mm		10X25mm		
Multi-Use CDEM362	Single-Use CDES362			4mm Shaft	Loop Electrode	BOX OF 50 UNITS
		160mm		15X25mm		
Multi-Use CDEM363	Single-Use CDES363			4mm Shaft	Loop Electrode	
		160mm		20X25mm		
Multi-Use CDEM364	Single-Use CDES364			4mm Shaft	Loop Electrode	
		140mm		1X10mm		
Multi-Use CDEM365	Single-Use CDES365			4mm Shaft	Coagulation Needle electrodes	
		140mm		4mm		
Multi-Use CDEM366	Single-Use CDES366			4mm Shaft	Coagulation Ball electrodes	

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES EXTENSIONS



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM370 Single-Use CDES370		100mm, 150mm 4mm Extension To Insert 2.4mm Electrode Shaft With Cage 4mm Shaft:	
Multi-Use CDEM371	Single-Use CDES371	4mm 2.4mm 4mm Extension To Insert 1.6mm Electrode 4mm Shaft	4mm 1.6mm 4mm Extension To Insert 2.4mm Electrode Multi-Use CDEM372 Single-Use CDES372
Multi-Use CDEM373	Single-Use CDES373	100mm, 150mm 2.4mm Extension To Insert 4mm Electrode Shaft With Cage 4mm Shaft	
Multi-Use CDEM374	Single-Use CDES374	100mm, 150mm 2.4mm Extension To Insert 1.6mm Electrode Shaft With Cage	
Multi-Use CDEM375	Single-Use CDES375	140mm 2x1.6mm 2.4x20mm Spatula Electrode 2.4mm Shaft	
Multi-Use CDEM376	Single-Use CDES376	140mm 2.4x10mm 3.3x25mm Knife Electrode 2.4mm Shaft	
Multi-Use CDEM377	Single-Use CDES377	140mm 0.8x12mm 0.5x20mm Needle Electrode 2.4mm Shaft	

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM378	Single-Use CDES378		
		2.4mm Shaft	Ball electrode
Multi-Use CDEM379	Single-Use CDES379		
		2.4mm Shaft	Ball electrode
Multi-Use CDEM380	Single-Use CDES380		
		2.4mm Shaft	Ball electrode
Multi-Use CDEM381	Single-Use CDES381		
		2.4mm Shaft	Ball electrode
Multi-Use CDEM382	Single-Use CDES382		
		2.4mm Shaft	Loop Electrode
Multi-Use CDEM383	Single-Use CDES383		
		2.4mm Shaft	Loop Electrode
Multi-Use CDEM384	Single-Use CDES384		
		2.4mm Shaft	Loop Electrode

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM385	Single-Use CDES385	160mm 2.4mmShaft Loop Electrode 10x10mm	
Multi-Use CDEM386	Single-Use CDES386	160mm 2.4mmShaft Loop Electrode 15x15mm	
Multi-Use CDEM387	Single-Use CDES387	160mm 2.4mmShaft Loop Electrode 15x20mm	
Multi-Use CDEM388	Single-Use CDES388	160mm 2.4mmShaft Loop Electrode 20x20mm	
Multi-Use CDEM389	Single-Use CDES389	160mm 2.4mmShaft Loop Electrode 20x25mm	
Multi-Use CDEM390	Single-Use CDES390	160mm 2.4mmShaft Loop Electrode 20x25mm	



DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MONOPOLAR ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM391Single-Use CDES391		160mm	10mm	2.4mmShaft	Loop Electrode
Multi-Use CDEM392Single-Use CDES392		160mm	15mm	2.4mmShaft	Loop Electrode
Multi-Use CDEM393Single-Use CDES393		160mm	20mm	2.4mmShaft	
Multi-Use CDEM394Single-Use CDES394		160mm	25mm	2.4mmShaft	



DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MICRODISSECTION NEEDLE ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS Multi-Use CDEM398 Single-Use CDES398				BOX OF 50 UNITS Multi-Use CDEM399 Single-Use CDES399	
Multi-Use CDEM400 Single-Use CDES400				Multi-Use CDEM401 Single-Use CDES401	
Multi-Use CDEM402 Single-Use CDES402					
Multi-Use CDEM403 Single-Use CDES403					
Multi-Use CDEM404 Single-Use CDES404				Multi-Use CDEM405 Single-Use CDES405	
Multi-Use CDEM406 Single-Use CDES406					
Multi-Use CDEM407 Single-Use CDES407					

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ACTIVE MICRODISSECTION NEEDLE ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS		
Multi-Use CDEM408	Single-Use CDES408	 Microdissection Needle Electrode, Flexible, Tungsten
Multi-Use CDEM409	Single-Use CDES409	 Microdissection Needle Electrode Flexible, Tungsten 4mmShaft
Multi-Use CDEM410	Single-Use CDES410	 Microdissection Needle Electrode Flexible, Tungsten 4mmShaft
Multi-Use CDEM411	Single-Use CDES411	 Microdissection Needle Electrode Flexible, Tungsten 4mmShaft
Multi-Use CDEM412	Single-Use CNES412	 Microdissection Needle Electrode Flexible, Tungsten 4mmShaft
Multi-Use CDEM413	Single-Use CDES413	 Microdissection Needle Electrode, Flexible, 45°, Tungsten

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH ARTHROSCOPIC ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS

Multi-Use
CDEM416
Single-Use
CDES416



Multi-Use
CDEM417
Single-Use
CDES417



Multi-Use
CDEM418
Single-Use
CDES418



Multi-Use
CDEM419
Single-Use
CDES419



Multi-Use
CDEM420
Single-Use
CNES420



Multi-Use
CDEM421
Single-Use
CDES421



Multi-Use
CDEM422
Single-Use
CDES422



Multi-Use
CDEM423
Single-Use
CDES423




Multi-Use
CDEM424
Single-Use
CDES424



BOX OF 50 UNITS

Multi-Use
CDEM425
Single-Use
CDES425




Multi-Use
CDEM426
Single-Use
CDES426




Multi-Use
CDEM428
Single-Use
CNES428



Multi-Use
CDEM430
Single-Use
CDES430



Multi-Use
CDEM432
Single-Use
CDES432



Multi-Use
CDEM427
Single-Use
CDES427



Multi-Use
CDEM429
Single-Use
CDES429



Multi-Use
CDEM431
Single-Use
CDES431



Multi-Use
CDEM433
Single-Use
CDES433

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ARTHROSCOPIC ELECTRODES



UMAX: 0.5KVp



BOX OF 50 UNITS		
Multi-Use CDEM434	Single-Use CDES434	Needle Electrode 1 '90° Angular With Handle 2.4mmShaft
Multi-Use CDEM435	Single-Use CDES435	Needle Electrode 90° Angular With Handle 2.4mmShaft
Multi-Use CDEM436	Single-Use CNES436	Knife Electrode 45° Angled Length 115 Mm 2.4mmShaft
Multi-Use CDEM413	Single-Use CDES413	Button Electrode With Handle 90° Angled Length 115 Mm 2.4mmShaft
Multi-Use CDEM437	Single-Use CDES437	Ball Electrode With Handle 2.4mmShaft

DIATHERMY CAUTERY
HAND PIECE & ELECTRODES

C-TECH ARTHROSCOPIC ELECTRODES



UMAX: 0.5KVP



BOX OF 50 UNITS		
Multi-Use CDEM438	Single-Use CDES438	Needle Electrode Angular with handle 4mmShaft
Multi-Use CDEM439	Single-Use CDES439	Knife Electrode 45° Angled 4mmShaft
Multi-Use CDEM440	Single-Use CDES440	Button Electrode with handle 90° Angled 4mmShaft
Multi-Use CDEM441	Single-Use CDES441	Ball Electrode with handle 4mmShaft
Multi-Use CDEM442	Single-Use CDES442	Needle Electrode 90° Angular, with handle 4mmShaft

DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH BIPOLAR ENT ELECTRODES




UMAX: 0.5KVP




Bipolar ENT Electrodes Steel

BOX OF 50 UNITS



Bipolar ENT Electrodes Teflon

BOX OF 50 UNITS



DIATHERMY CAUTERY HAND PIECE & ELECTRODES

C-TECH LAPROSCOPIC ELECTRODES



UMAX: 0.5KVP



Laprosopic Electrodes without Suction

BOX OF 50 UNITS



Multi-Use CDEM454
Single-Use CDES454



Multi-Use CDEM455
Single-Use CDES455



Multi-Use CDEM456
Single-Use CNES456



Multi-Use CDEM457
Single-Use CDES457



Multi-Use CDEM458
Single-Use CDES458

Laprosopic Electrodes with Suction

BOX OF 50 UNITS



Multi-Use CDEM459
Single-Use CDES459



Multi-Use CDEM460
Single-Use CDES460



Multi-Use CDEM461
Single-Use CNES461



Multi-Use CDEM462
Single-Use CDES462



Multi-Use CDEM463
Single-Use CDES463



06

ELECTROSURGICAL ELECTRODE PAD

PRECISE. INSULATED. RELIABLE.

ELECTROSURGICAL DIATHERMY ELECTRODE PADS

CLASSIC, UPGRADED



UMAX: 0.5KVP



Classic Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Adult Plate Dimension 260x150x3.3mm




Classic Silicon Rubber Patient Return Plate – Without Cable



Adult Plate Dimension
260x150x3.3mm
CDP103

Child Child Plate Dimension
200x80x3.3mm
CDP104

Upgraded – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Upgraded – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



ELECTROSURGICAL DIATHERMY ELECTRODE PADS

PREMIUM, WITH CABLE



UMAX: 0.5KVp

Premium – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Bipolar Adult Plate Dimension 200x150x3.3mm



ERBE, Berchtold, Tekno Series
Generators
CDP109



Valleylab, Conmed, BOWA Series
Generators
CDP110

Premium – Silicon Rubber Patient Return Plate – Available in 2 Different Connectors



Monopolar Adult Plate Dimension 200x150x3.3mm




ERBE, Berchtold, Tekno Series
Generators



Valleylab, Conmed, BOWA Series
Generators
CDP112

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



ERBE, Berchtold, Tekno Series
Generators

CDP113

Single-Use Monopolar
Patient Return Plates - Child

CDP114

Single-Use Monopolar
Patient Return Plates - Adult




Valleylab, Conmed, BOWA Series
Generators

CDP115

CDP116

Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



ERBE, Berchtold, Tekno Series
Generators

CDP117

Single-Use Monopolar
Patient Return Plates - Child

CDP118

Single-Use Monopolar
Patient Return Plates - Adult



Valleylab, Conmed, BOWA Series
Generators

CDP119

CDP120

ELECTROSURGICAL DIATHERMY
ELECTRODE PADS

WITHOUT CABLE, CABLE, TIP CLEANER



UMAX: 0.5KVp



Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



Single-Use Monopolar
Patient Return Plates - Child

CDP121

Single-Use Monopolar
Patient Return Plates - Adult

CDP122




ERBE, Berchtold, Tekno Series
Generators

Valleylab, Conmed, BOWA Series
Generators

Multi-Use	CDP125	High Frequency Jack for	CDP127
Single-Use	CDP126	Split/Round Bipolar Plates	CDP128
Multi-Use	CDP129	High Frequency Jack for	CDP131
Single-Use	CDP130	Split/Round Bipolar Plates	CDP132



Disposable Patient Return Plate With Cable – Available in 2 Different Connectors



Single-Use Monopolar
Patient Return Plates - Child

CDP123

Single-Use Monopolar
Patient Return Plates - Adult

CDP124

Tip Cleaner Pads – Available in 3 Different Sizes



Single-Use Tip Cleaner Pad

CDP133	45x45mm
CDP134	48x48mm
CDP135	50x50mm



07

SINGLE USE
FORCEPS
INTEGRATED CABLE

PRECISE. INSULATED. RELIABLE.

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-BLUE

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp

Regular Single Use Mcpherson Bipolar Forceps with Cable-3m

		102mm	108mm	115mm
 International 4mm Banana Plug Fitting	CSUB101T	CSUB101U	CSUB101V	 0.50mm
	CSUB102T	CSUB102U	CSUB102V	 1.00mm
	CSUB103T	CSUB103U	CSUB103V	 2.00mm
 BOWA, Limdey, EMC Valleylab, Ethicon	CSUB104T	CSUB104U	CSUB104V	 0.50mm
	CSUB105T	CSUB105U	CSUB105V	 1.00mm
	CSUB106T	CSUB106U	CSUB106V	 2.00mm
 ERBE 4mm, Siemens Shutter BM50	CSUB107T	CSUB107U	CSUB107V	 0.50mm
	CSUB108T	CSUB108U	CSUB108V	 1.00mm
	CSUB109T	CSUB109U	CSUB109V	 2.00mm
 Aesculap, Berchtold, Wolf, Martin	CSUB110T	CSUB110U	CSUB110V	 0.50mm
	CSUB111T	CSUB111U	CSUB111V	 1.00mm
	CSUB112T	CSUB112U	CSUB112V	 2.00mm
 Eichmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB113T	CSUB113U	CSUB113V	 0.50mm
	CSUB114T	CSUB114U	CSUB114V	 1.00mm
	CSUB115T	CSUB115U	CSUB115V	 2.00mm

Box of 20 Units

Regular Single Use Mcpherson Bipolar Forceps Angled with Cable-3m

					
		102mm	108mm	115mm	
	mm	CSUB116T	CSUB116U	CSUB116V	 0.50mm
	mm	CSUB117T	CSUB117U	CSUB117V	 1.00mm
		CSUB118T	CSUB118U	CSUB118V	 2.00mm
	22mm	CSUB119T	CSUB119U	CSUB119V	 0.50mm
		CSUB120T	CSUB120U	CSUB120V	 1.00mm
		CSUB121T	CSUB121U	CSUB121V	 2.00mm
	22mm	CSUB122T	CSUB122U	CSUB122V	 0.50mm
		CSUB123T	CSUB123U	CSUB123V	 1.00mm
		CSUB124T	CSUB124U	CSUB124V	 2.00mm
	22mm	CSUB125T	CSUB125U	CSUB125V	 0.50mm
		CSUB126T	CSUB126U	CSUB126V	 1.00mm
		CSUB127T	CSUB127U	CSUB127V	 2.00mm
		CSUB128T	CSUB128U	CSUB128V	 0.50mm
		CSUB129T	CSUB129U	CSUB129V	 1.00mm
		CSUB130T	CSUB130U	CSUB130V	 2.00mm

Box of 20 Units

Regular Single Use Jeweller Bipolar Forceps with Cable-3m

				
102mm		115mm	125mm	
	CSUB131T	CSUB131V	CSUB131Y	 0.50mm
International 4mm Banana Plug Fitting	CSUB132T	CSUB132V	CSUB132Y	 1.00mm
	CSUB133T	CSUB133V	CSUB133Y	 2.00mm
	CSUB134T	CSUB134V	CSUB134Y	 0.50mm
BOWA, Limdey, EMC Valleylab, Ethicon	CSUB135T	CSUB135V	CSUB135Y	 1.00mm
	CSUB136T	CSUB136V	CSUB136Y	 2.00mm
	CSUB137T	CSUB137V	CSUB137Y	 0.50mm
ERBE 4mm, Siemens Shutter BM50	CSUB138T	CSUB138V	CSUB138Y	 1.00mm
	CSUB139T	CSUB139V	CSUB139Y	 2.00mm
	CSUB140T	CSUB140V	CSUB140Y	 0.50mm
Aesculap, Berchtold, Wolf, Martin	CSUB141T	CSUB141V	CSUB141Y	 1.00mm
	CSUB142T	CSUB142V	CSUB142Y	 2.00mm
	CSUB143T	CSUB143V	CSUB143Y	 0.50mm
Eichmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUB144T	CSUB144V	CSUB144Y	 1.00mm
	CSUB145T	CSUB145V	CSUB145Y	 2.00mm

Box of 20 Units

Regular Single Use Jeweller Bipolar Forceps Angled with Cable-3m

102mm			115mm			125mm					
	CSUB146T	CSUB146V	CSUB146Y		0.50mm		CSUB147T	CSUB147V	CSUB147Y		1.00mm
	CSUB147T	CSUB147V	CSUB147Y		2.00mm		CSUB148T	CSUB148V	CSUB148Y		2.00mm
	CSUB148T	CSUB148V	CSUB148Y		2.00mm						
	CSUB149T	CSUB149V	CSUB149Y		0.50mm		CSUB150T	CSUB150V	CSUB150Y		1.00mm
	CSUB150T	CSUB150V	CSUB150Y		2.00mm		CSUB151T	CSUB151V	CSUB151Y		2.00mm
	CSUB151T	CSUB151V	CSUB151Y		2.00mm						
	CSUB152T	CSUB152V	CSUB152Y		0.50mm		CSUB153T	CSUB153V	CSUB153Y		1.00mm
	CSUB153T	CSUB153V	CSUB153Y		2.00mm		CSUB154T	CSUB154V	CSUB154Y		2.00mm
	CSUB154T	CSUB154V	CSUB154Y		2.00mm						
	CSUB155T	CSUB155V	CSUB155Y		0.50mm		CSUB156T	CSUB156V	CSUB156Y		1.00mm
	CSUB156T	CSUB156V	CSUB156Y		2.00mm		CSUB157T	CSUB157V	CSUB157Y		2.00mm
	CSUB157T	CSUB157V	CSUB157Y		2.00mm						
	CSUB158T	CSUB158V	CSUB158Y		0.50mm		CSUB159T	CSUB159V	CSUB159Y		1.00mm
	CSUB159T	CSUB159V	CSUB159Y		2.00mm		CSUB160T	CSUB160V	CSUB160Y		2.00mm
	CSUB160T	CSUB160V	CSUB160Y		2.00mm						

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-BLUE

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp

Regular Single Use Adson Bipolar Forceps with Cable-3m

	121mm	155mm	180mm	200mm	
International 4mm Banana Plug Fitting	CSUB161W CSUB162W CSUB163W	CSUB161A CSUB162A CSUB163A	CSUB161B CSUB162B CSUB163B	CSUB161C CSUB162C CSUB163C	0.50mm 1.00mm 2.00mm
BOWA, Lumidey, EMC Valleylab, Ethicon	CSUB164W CSUB165W CSUB166W	CSUB164A CSUB165A CSUB166A	CSUB164B CSUB165B CSUB166B	CSUB164C CSUB165C CSUB166C	0.50mm 1.00mm 2.00mm
ERBE 4mm, Siemens Shutter BMSO	CSUB167W CSUB168W CSUB169W	CSUB167A CSUB168A CSUB169A	CSUB167B CSUB168B CSUB169B	CSUB167C CSUB168C CSUB169C	0.50mm 1.00mm 2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB170W CSUB171W CSUB172W	CSUB170A CSUB171A CSUB172A	CSUB170B CSUB171B CSUB172B	CSUB170C CSUB171C CSUB172C	0.50mm 1.00mm 2.00mm
Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Cornead 30mm	CSUB173W CSUB174W CSUB175W	CSUB173A CSUB174A CSUB175A	CSUB173B CSUB174B CSUB175B	CSUB173C CSUB174C CSUB175C	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps with Cable-3m

	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB191X CSUB192X CSUB193X	CSUB191Z CSUB192Z CSUB193Z	0.50mm 1.00mm 2.00mm
BOWA, Lumidey, EMC Valleylab, Ethicon	CSUB194X CSUB195X CSUB196X	CSUB194Z CSUB195Z CSUB196Z	0.50mm 1.00mm 2.00mm
ERBE 4mm, Siemens Shutter BMSO	CSUB197X CSUB198X CSUB199X	CSUB197Z CSUB198Z CSUB199Z	0.50mm 1.00mm 2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB200X CSUB201X CSUB202X	CSUB200Z CSUB201Z CSUB202Z	0.50mm 1.00mm 2.00mm
Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Cornead 30mm	CSUB203X CSUB204X CSUB205X	CSUB203Z CSUB204Z CSUB205Z	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps with Cable-3m

	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB176X CSUB177X CSUB178X	CSUB176Z CSUB177Z CSUB178Z	0.50mm 1.00mm 2.00mm
BOWA, Lumidey, EMC Valleylab, Ethicon	CSUB179X CSUB180X CSUB181X	CSUB179Z CSUB180Z CSUB181Z	0.50mm 1.00mm 2.00mm
ERBE 4mm, Siemens Shutter BMSO	CSUB182X CSUB183X CSUB184X	CSUB182Z CSUB183Z CSUB184Z	0.50mm 1.00mm 2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB185X CSUB186X CSUB187X	CSUB185Z CSUB186Z CSUB187Z	0.50mm 1.00mm 2.00mm
Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Cornead 30mm	CSUB188X CSUB189X CSUB190X	CSUB188Z CSUB189Z CSUB190Z	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps Angled with Cable-3m

	140mm	165mm	
International 4mm Banana Plug Fitting	CSUB206X CSUB207X CSUB208X	CSUB206Z CSUB207Z CSUB208Z	0.50mm 1.00mm 2.00mm
BOWA, Lumidey, EMC Valleylab, Ethicon	CSUB209X CSUB210X CSUB211X	CSUB209Z CSUB210Z CSUB211Z	0.50mm 1.00mm 2.00mm
ERBE 4mm, Siemens Shutter BMSO	CSUB212X CSUB213X CSUB214X	CSUB212Z CSUB213Z CSUB214Z	0.50mm 1.00mm 2.00mm
Aesculap, Berchtold, Wolf, Martin	CSUB215X CSUB216X CSUB217A	CSUB215Z CSUB216Z CSUB217Z	0.50mm 1.00mm 2.00mm
Eschmann 18mm, Ellman 22mm, Vallylab 28mm, Cornead 30mm	CSUB218X CSUB219X CSUB220X	CSUB218Z CSUB219Z CSUB220Z	0.50mm 1.00mm 2.00mm

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-BLUE

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp

Regular Single Use Straight Bipolar Forceps with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB221A	CSUB221B	CSUB221C	CSUB221D	CSUB221E	0.50mm
	CSUB222A	CSUB222B	CSUB222C	CSUB222D	CSUB222E	1.00mm
	CSUB223A	CSUB223B	CSUB223C	CSUB223D	CSUB223E	2.00mm
BOWA Lumsley, EMC Valleylab, Ethicon	CSUB224A	CSUB224B	CSUB224C	CSUB224D	CSUB224E	0.50mm
	CSUB225A	CSUB225B	CSUB225C	CSUB225D	CSUB225E	1.00mm
	CSUB226A	CSUB226B	CSUB226C	CSUB226D	CSUB226E	2.00mm
ERBE 4mm, Siemens Shutter BMS0	CSUB227A	CSUB227B	CSUB227C	CSUB227D	CSUB227E	0.50mm
	CSUB228A	CSUB228B	CSUB228C	CSUB228D	CSUB228E	1.00mm
	CSUB229A	CSUB229B	CSUB229C	CSUB229D	CSUB229E	2.00mm
Aesculap, Bercitold, Wolf, Martin	CSUB230A	CSUB230B	CSUB230C	CSUB230D	CSUB230E	0.50mm
	CSUB231A	CSUB231B	CSUB231C	CSUB231D	CSUB231E	1.00mm
	CSUB232A	CSUB232B	CSUB232C	CSUB232D	CSUB232E	2.00mm
Eschmann 18mm, Elman 22mm, Valleylab 28mm, Conmed 30mm	CSUB233A	CSUB233B	CSUB233C	CSUB233D	CSUB233E	0.50mm
	CSUB234A	CSUB234B	CSUB234C	CSUB234D	CSUB234E	1.00mm
	CSUB235A	CSUB235B	CSUB235C	CSUB235D	CSUB235E	2.00mm

Box of 20 Units

Regular Single Use Straight Bipolar Forceps Angled with Cable-3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB236A	CSUB236B	CSUB236C	CSUB236D	CSUB236E	0.50mm
	CSUB237A	CSUB237B	CSUB237C	CSUB237D	CSUB237E	1.00mm
	CSUB238A	CSUB238B	CSUB238C	CSUB238D	CSUB238E	2.00mm
BOWA Lumsley, EMC Valleylab, Ethicon	CSUB239A	CSUB239B	CSUB239C	CSUB239D	CSUB239E	0.50mm
	CSUB240A	CSUB240B	CSUB240C	CSUB240D	CSUB240E	1.00mm
	CSUB241A	CSUB241B	CSUB241C	CSUB241D	CSUB241E	2.00mm
ERBE 4mm, Siemens Shutter BMS0	CSUB242A	CSUB242B	CSUB242C	CSUB242D	CSUB242E	0.50mm
	CSUB243A	CSUB243B	CSUB243C	CSUB243D	CSUB243E	1.00mm
	CSUB244A	CSUB244B	CSUB244C	CSUB244D	CSUB244E	2.00mm
Aesculap, Bercitold, Wolf, Martin	CSUB245A	CSUB245B	CSUB245C	CSUB245D	CSUB245E	0.50mm
	CSUB246A	CSUB246B	CSUB246C	CSUB246D	CSUB246E	1.00mm
	CSUB247A	CSUB247B	CSUB247C	CSUB247D	CSUB247E	2.00mm
Eschmann 18mm, Elman 22mm, Valleylab 28mm, Conmed 30mm	CSUB248A	CSUB248B	CSUB248C	CSUB248D	CSUB248E	0.50mm
	CSUB249A	CSUB249B	CSUB249C	CSUB249D	CSUB249E	1.00mm
	CSUB250A	CSUB250B	CSUB250C	CSUB250D	CSUB250E	2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB251A	CSUB251B	CSUB251C	CSUB251D	CSUB251E	0.50mm
	CSUB252A	CSUB252B	CSUB252C	CSUB252D	CSUB252E	1.00mm
	CSUB253A	CSUB253B	CSUB253C	CSUB253D	CSUB253E	2.00mm
BOWA Lumsley, EMC Valleylab, Ethicon	CSUB254A	CSUB254B	CSUB254C	CSUB254D	CSUB254E	0.50mm
	CSUB255A	CSUB255B	CSUB255C	CSUB255D	CSUB255E	1.00mm
	CSUB256A	CSUB256B	CSUB256C	CSUB256D	CSUB256E	2.00mm
ERBE 4mm, Siemens Shutter BMS0	CSUB257A	CSUB257B	CSUB257C	CSUB257D	CSUB257E	0.50mm
	CSUB258A	CSUB258B	CSUB258C	CSUB258D	CSUB258E	1.00mm
	CSUB259A	CSUB259B	CSUB259C	CSUB259D	CSUB259E	2.00mm
Aesculap, Bercitold, Wolf, Martin	CSUB260A	CSUB260B	CSUB260C	CSUB260D	CSUB260E	0.50mm
	CSUB261A	CSUB261B	CSUB261C	CSUB261D	CSUB261E	1.00mm
	CSUB262A	CSUB262B	CSUB262C	CSUB262D	CSUB262E	2.00mm
Eschmann 18mm, Elman 22mm, Valleylab 28mm, Conmed 30mm	CSUB263A	CSUB263B	CSUB263C	CSUB263D	CSUB263E	0.50mm
	CSUB264A	CSUB264B	CSUB264C	CSUB264D	CSUB264E	1.00mm
	CSUB265A	CSUB265B	CSUB265C	CSUB265D	CSUB265E	2.00mm

Box of 20 Units

Regular Single Use Bayonet Bipolar Forceps Angled with Cable -3m



	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUB266A	CSUB266B	CSUB266C	CSUB266D	CSUB266E	0.50mm
	CSUB267A	CSUB267B	CSUB267C	CSUB267D	CSUB267E	1.00mm
	CSUB268A	CSUB268B	CSUB268C	CSUB268D	CSUB268E	2.00mm
BOWA Lumsley, EMC Valleylab, Ethicon	CSUB269A	CSUB269B	CSUB269C	CSUB269D	CSUB269E	0.50mm
	CSUB270A	CSUB270B	CSUB270C	CSUB270D	CSUB270E	1.00mm
	CSUB271A	CSUB271B	CSUB271C	CSUB271D	CSUB271E	2.00mm
ERBE 4mm, Siemens Shutter BMS0	CSUB272A	CSUB272B	CSUB272C	CSUB272D	CSUB272E	0.50mm
	CSUB273A	CSUB273B	CSUB273C	CSUB273D	CSUB273E	1.00mm
	CSUB274A	CSUB274B	CSUB274C	CSUB274D	CSUB274E	2.00mm
Aesculap, Bercitold, Wolf, Martin	CSUB275A	CSUB275B	CSUB275C	CSUB275D	CSUB275E	0.50mm
	CSUB276A	CSUB276B	CSUB276C	CSUB276D	CSUB276E	1.00mm
	CSUB277A	CSUB277B	CSUB277C	CSUB277D	CSUB277E	2.00mm
Eschmann 18mm, Elman 22mm, Valleylab 28mm, Conmed 30mm	CSUB278A	CSUB278B	CSUB278C	CSUB278D	CSUB278E	0.50mm
	CSUB279A	CSUB279B	CSUB279C	CSUB279D	CSUB279E	1.00mm
	CSUB280A	CSUB280B	CSUB280C	CSUB280D	CSUB280E	2.00mm

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-PUMPKIN

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVP

Non-Stick Single Use Mcpherson Bipolar Forceps with Cable-3m

102mm			108mm			115mm		
	CSUP281T	0.50mm		CSUP281U	0.50mm		CSUP281V	0.50mm
	CSUP282T	1.00mm		CSUP282U	1.00mm		CSUP282V	1.00mm
	CSUP283T	2.00mm		CSUP283U	2.00mm		CSUP283V	2.00mm
	CSUP284T	0.50mm		CSUP284U	0.50mm		CSUP284V	0.50mm
	CSUP285T	1.00mm		CSUP285U	1.00mm		CSUP285V	1.00mm
	CSUP286T	2.00mm		CSUP286U	2.00mm		CSUP286V	2.00mm
	CSUP287T	0.50mm		CSUP287U	0.50mm		CSUP287V	0.50mm
	CSUP288T	1.00mm		CSUP288U	1.00mm		CSUP288V	1.00mm
	CSUP289T	2.00mm		CSUP289U	2.00mm		CSUP289V	2.00mm
	CSUP290T	0.50mm		CSUP290U	0.50mm		CSUP290V	0.50mm
	CSUP291T	1.00mm		CSUP291U	1.00mm		CSUP291V	1.00mm
	CSUP292T	2.00mm		CSUP292U	2.00mm		CSUP292V	2.00mm
	CSUP293T	0.50mm		CSUP293U	0.50mm		CSUP293V	0.50mm
	CSUP294T	1.00mm		CSUP294U	1.00mm		CSUP294V	1.00mm
	CSUP295T	2.00mm		CSUP295U	2.00mm		CSUP295V	2.00mm

Box of 20 Units

Non-Stick Single Use Jeweller Bipolar Forceps with Cable-3m

102mm			115mm			125mm		
	CSUB311T	0.50mm		CSUB311V	0.50mm		CSUB311Y	0.50mm
	CSUB312T	1.00mm		CSUB312V	1.00mm		CSUB312Y	1.00mm
	CSUB313T	2.00mm		CSUB313V	2.00mm		CSUB313Y	2.00mm
	CSUB314T	0.50mm		CSUB314V	0.50mm		CSUB314Y	0.50mm
	CSUB315T	1.00mm		CSUB315V	1.00mm		CSUB315Y	1.00mm
	CSUB316T	2.00mm		CSUB316V	2.00mm		CSUB316Y	2.00mm
	CSUB317T	0.50mm		CSUB317V	0.50mm		CSUB317Y	0.50mm
	CSUB318T	1.00mm		CSUB318V	1.00mm		CSUB318Y	1.00mm
	CSUB319T	2.00mm		CSUB319V	2.00mm		CSUB319Y	2.00mm
	CSUB320T	0.50mm		CSUB320V	0.50mm		CSUB320Y	0.50mm
	CSUB321T	1.00mm		CSUB321V	1.00mm		CSUB321Y	1.00mm
	CSUB322T	2.00mm		CSUB322V	2.00mm		CSUB322Y	2.00mm
	CSUB323T	0.50mm		CSUB323V	0.50mm		CSUB323Y	0.50mm
	CSUB324T	1.00mm		CSUB324V	1.00mm		CSUB324Y	1.00mm
	CSUB325T	2.00mm		CSUB325V	2.00mm		CSUB325Y	2.00mm

Box of 20 Units

Non-Stick Single Use Mcpherson Bipolar Forceps Angled with Cable-3m

102mm			108mm			115mm		
	CSUP296T	0.50mm		CSUP296U	0.50mm		CSUP296V	0.50mm
	CSUP297T	1.00mm		CSUP297U	1.00mm		CSUP297V	1.00mm
	CSUP298T	2.00mm		CSUP298U	2.00mm		CSUP298V	2.00mm
	CSUP299T	0.50mm		CSUP299U	0.50mm		CSUP299V	0.50mm
	CSUP300T	1.00mm		CSUP300U	1.00mm		CSUP300V	1.00mm
	CSUP301T	2.00mm		CSUP301U	2.00mm		CSUP301V	2.00mm
	CSUP302T	0.50mm		CSUP302U	0.50mm		CSUP302V	0.50mm
	CSUP303T	1.00mm		CSUP303U	1.00mm		CSUP303V	1.00mm
	CSUP304T	2.00mm		CSUP304U	2.00mm		CSUP304V	2.00mm
	CSUP305T	0.50mm		CSUP305U	0.50mm		CSUP305V	0.50mm
	CSUP306T	1.00mm		CSUP306U	1.00mm		CSUP306V	1.00mm
	CSUP307T	2.00mm		CSUP307U	2.00mm		CSUP307V	2.00mm
	CSUP308T	0.50mm		CSUP308U	0.50mm		CSUP308V	0.50mm
	CSUP309T	1.00mm		CSUP309U	1.00mm		CSUP309V	1.00mm
	CSUP310T	2.00mm		CSUP310U	2.00mm		CSUP310V	2.00mm

Box of 20 Units

Non-Stick Single Use Jeweller Bipolar Forceps Angled with Cable-3m

102mm			115mm			125mm		
	CSUP326T	0.50mm		CSUP326V	0.50mm		CSUP326Y	0.50mm
	CSUP327T	1.00mm		CSUP327V	1.00mm		CSUP327Y	1.00mm
	CSUP328T	2.00mm		CSUP328V	2.00mm		CSUP328Y	2.00mm
	CSUP329T	0.50mm		CSUP329V	0.50mm		CSUP329Y	0.50mm
	CSUP330T	1.00mm		CSUP330V	1.00mm		CSUP330Y	1.00mm
	CSUP331T	2.00mm		CSUP331V	2.00mm		CSUP331Y	2.00mm
	CSUP332T	0.50mm		CSUP332V	0.50mm		CSUP332Y	0.50mm
	CSUP333T	1.00mm		CSUP333V	1.00mm		CSUP333Y	1.00mm
	CSUP334T	2.00mm		CSUP334V	2.00mm		CSUP334Y	2.00mm
	CSUP335T	0.50mm		CSUP335V	0.50mm		CSUP335Y	0.50mm
	CSUP336T	1.00mm		CSUP336V	1.00mm		CSUP336Y	1.00mm
	CSUP337T	2.00mm		CSUP337V	2.00mm		CSUP337Y	2.00mm
	CSUP338T	0.50mm		CSUP338V	0.50mm		CSUP338Y	0.50mm
	CSUP339T	1.00mm		CSUP339V	1.00mm		CSUP339Y	1.00mm
	CSUP340T	2.00mm		CSUP340V	2.00mm		CSUP340Y	2.00mm

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-PUMPKIN

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp

Non-Stick Single Use Adson Bipolar Forceps with Cable-3m

121mm	155mm	180mm	200mm	
 CUSP341W CUSP342W CUSP343W	 CUSP344W CUSP345W CUSP346W	 CUSP347W CUSP348W CUSP349W	 CUSP350W CUSP351W CUSP352W	 CUSP353W CUSP354W CUSP355W
CUSP341A CUSP342A CUSP343A	CUSP344A CUSP345A CUSP346A	CUSP347A CUSP348A CUSP349A	CUSP350A CUSP351A CUSP352A	CUSP353A CUSP354A CUSP355A
CUSP341B CUSP342B CUSP343B	CUSP344B CUSP345B CUSP346B	CUSP347B CUSP348B CUSP349B	CUSP350B CUSP351B CUSP352B	CUSP353B CUSP354B CUSP355B
CUSP341C CUSP342C CUSP343C	CUSP344C CUSP345C CUSP346C	CUSP347C CUSP348C CUSP349C	CUSP350C CUSP351C CUSP352C	CUSP353C CUSP354C CUSP355C
0.50mm	0.50mm	0.50mm	0.50mm	0.50mm
1.00mm	1.00mm	1.00mm	1.00mm	1.00mm
2.00mm	2.00mm	2.00mm	2.00mm	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps with Cable-3m

140mm	165mm	
 CSUP356X CSUP357X CSUP358X	 CSUP359X CSUP360X CSUP361X	 CSUP362X CSUP363X CSUP364X
 CSUP365X CSUP366X CSUP367X	 CSUP368X CSUP369X CSUP370X	CSUP356Z CSUP357Z CSUP358Z
CSUP359Z CSUP360Z CSUP361Z	CSUP362Z CSUP363Z CSUP364Z	CSUP365Z CSUP366Z CSUP367Z
0.50mm	0.50mm	0.50mm
1.00mm	1.00mm	1.00mm
2.00mm	2.00mm	2.00mm

Box of 20 Units

Non-Stick Single Use Bayonet Bipolar Forceps with Cable-3m

140mm	165mm	
 CSUP371X CSUP372X CSUP373X	 CSUP374X CSUP375X CSUP376X	 CSUP377X CSUP378X CSUP379X
 CSUP380X CSUP381X CSUP382X	 CSUP383X CSUP384X CSUP385X	CSUP371Z CSUP372Z CSUP373Z
CSUP374Z CSUP375Z CSUP376Z	CSUP377Z CSUP378Z CSUP379Z	CSUP380Z CSUP381Z CSUP382Z
0.50mm	0.50mm	0.50mm
1.00mm	1.00mm	1.00mm
2.00mm	2.00mm	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps Angled with Cable-3m

140mm	165mm	
 CSUP386X CSUP387X CSUP388X	 CSUP389X CSUP390X CSUP391X	 CSUP392X CSUP393X CSUP394X
 CSUP395X CSUP396X CSUP397X	 CSUP398X CSUP399X CSUP400X	CSUP386Z CSUP387Z CSUP388Z
CSUP389Z CSUP390Z CSUP391Z	CSUP392Z CSUP393Z CSUP394Z	CSUP395Z CSUP396Z CSUP397Z
0.50mm	0.50mm	0.50mm
1.00mm	1.00mm	1.00mm
2.00mm	2.00mm	2.00mm

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-PUMPKIN

STERILE EO
ETHYLENE OXIDE
STERILIZED

SINGLE USE

max
Vp
MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp



Non-Stick Single Use Straight Bipolar Forceps with Cable-3m

	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP401A	CSUP401B	CSUP401C	CSUP401D	CSUP401E	0.50mm
	CSUP402A	CSUP402B	CSUP402C	CSUP402D	CSUP402E	1.00mm
	CSUP403A	CSUP403B	CSUP403C	CSUP403D	CSUP403E	2.00mm
BOWA, Lummed, EMC, Valleylab, Ethicon	CSUP404A	CSUP404B	CSUP404C	CSUP404D	CSUP404E	0.50mm
	CSUP405A	CSUP405B	CSUP405C	CSUP405D	CSUP405E	1.00mm
	CSUP406A	CSUP406B	CSUP406C	CSUP406D	CSUP406E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP407A	CSUP407B	CSUP407C	CSUP407D	CSUP407E	0.50mm
	CSUP408A	CSUP408B	CSUP408C	CSUP408D	CSUP408E	1.00mm
	CSUP409A	CSUP409B	CSUP409C	CSUP409D	CSUP409E	2.00mm
Aesculap, Becthield, Wolf, Martin	CSUP410A	CSUP410B	CSUP410C	CSUP410D	CSUP410E	0.50mm
	CSUP411A	CSUP411B	CSUP411C	CSUP411D	CSUP411E	1.00mm
	CSUP412A	CSUP412B	CSUP412C	CSUP412D	CSUP412E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP413A	CSUP413B	CSUP413C	CSUP413D	CSUP413E	0.50mm
	CSUP414A	CSUP414B	CSUP414C	CSUP414D	CSUP414E	1.00mm
	CSUP415A	CSUP415B	CSUP415C	CSUP415D	CSUP415E	2.00mm

Box of 20 Units

Non-Stick Single Use Straight Bipolar Forceps Angled with Cable-3m

	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP416A	CSUP416B	CSUP416C	CSUP416D	CSUP416E	0.50mm
	CSUP417A	CSUP417B	CSUP417C	CSUP417D	CSUP417E	1.00mm
	CSUP418A	CSUP418B	CSUP418C	CSUP418D	CSUP418E	2.00mm
BOWA, Lummed, EMC, Valleylab, Ethicon	CSUP419A	CSUP419B	CSUP419C	CSUP419D	CSUP419E	0.50mm
	CSUP420A	CSUP420B	CSUP420C	CSUP420D	CSUP420E	1.00mm
	CSUP421A	CSUP421B	CSUP421C	CSUP421D	CSUP421E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP422A	CSUP422B	CSUP422C	CSUP422D	CSUP422E	0.50mm
	CSUP423A	CSUP423B	CSUP423C	CSUP423D	CSUP423E	1.00mm
	CSUP424A	CSUP424B	CSUP424C	CSUP424D	CSUP424E	2.00mm
Aesculap, Becthield, Wolf, Martin	CSUP425A	CSUP425B	CSUP425C	CSUP425D	CSUP425E	0.50mm
	CSUP426A	CSUP426B	CSUP426C	CSUP426D	CSUP426E	1.00mm
	CSUP427A	CSUP427B	CSUP427C	CSUP427D	CSUP427E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP428A	CSUP428B	CSUP428C	CSUP428D	CSUP428E	0.50mm
	CSUP429A	CSUP429B	CSUP429C	CSUP429D	CSUP429E	1.00mm
	CSUP430A	CSUP430B	CSUP430C	CSUP430D	CSUP430E	2.00mm

Box of 20 Units

Non-Stick Single Use Bayonet Bipolar Forceps with Cable -3m

	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP431A	CSUP431B	CSUP431C	CSUP431D	CSUP431E	0.50mm
	CSUP432A	CSUP432B	CSUP432C	CSUP432D	CSUP432E	1.00mm
	CSUP433A	CSUP433B	CSUP433C	CSUP433D	CSUP433E	2.00mm
BOWA, Lummed, EMC, Valleylab, Ethicon	CSUP434A	CSUP434B	CSUP434C	CSUP434D	CSUP434E	0.50mm
	CSUP435A	CSUP435B	CSUP435C	CSUP435D	CSUP435E	1.00mm
	CSUP436A	CSUP436B	CSUP436C	CSUP436D	CSUP436E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP437A	CSUP437B	CSUP437C	CSUP437D	CSUP437E	0.50mm
	CSUP438A	CSUP438B	CSUP438C	CSUP438D	CSUP438E	1.00mm
	CSUP439A	CSUP439B	CSUP439C	CSUP439D	CSUP439E	2.00mm
Aesculap, Becthield, Wolf, Martin	CSUP440A	CSUP440B	CSUP440C	CSUP440D	CSUP440E	0.50mm
	CSUP441A	CSUP441B	CSUP441C	CSUP441D	CSUP441E	1.00mm
	CSUP442A	CSUP442B	CSUP442C	CSUP442D	CSUP442E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP443A	CSUP443B	CSUP443C	CSUP443D	CSUP443E	0.50mm
	CSUP444A	CSUP444B	CSUP444C	CSUP444D	CSUP444E	1.00mm
	CSUP445A	CSUP445B	CSUP445C	CSUP445D	CSUP445E	2.00mm

Box of 20 Units


Non-Stick Single Use Bayonet Bipolar Forceps Angled with Cable -3m

	155mm	180mm	200mm	220mm	240mm	
International 4mm Banana Plug Fitting	CSUP446A	CSUP446B	CSUP446C	CSUP446D	CSUP446E	0.50mm
	CSUP447A	CSUP447B	CSUP447C	CSUP447D	CSUP447E	1.00mm
	CSUP448A	CSUP448B	CSUP448C	CSUP448D	CSUP448E	2.00mm
BOWA, Lummed, EMC, Valleylab, Ethicon	CSUP449A	CSUP449B	CSUP449C	CSUP449D	CSUP449E	0.50mm
	CSUP450A	CSUP450B	CSUP450C	CSUP450D	CSUP450E	1.00mm
	CSUP451A	CSUP451B	CSUP451C	CSUP451D	CSUP451E	2.00mm
ERBE 4mm, Siemens Shutter BMS50	CSUP452A	CSUP452B	CSUP452C	CSUP452D	CSUP452E	0.50mm
	CSUP453A	CSUP453B	CSUP453C	CSUP453D	CSUP453E	1.00mm
	CSUP454A	CSUP454B	CSUP454C	CSUP454D	CSUP454E	2.00mm
Aesculap, Becthield, Wolf, Martin	CSUP455A	CSUP455B	CSUP455C	CSUP455D	CSUP455E	0.50mm
	CSUP456A	CSUP456B	CSUP456C	CSUP456D	CSUP456E	1.00mm
	CSUP457A	CSUP457B	CSUP457C	CSUP457D	CSUP457E	2.00mm
Eschmann 18mm, Ellman 22mm, Valleylab 28mm, Conmed 30mm	CSUP458A	CSUP458B	CSUP458C	CSUP458D	CSUP458E	0.50mm
	CSUP459A	CSUP459B	CSUP459C	CSUP459D	CSUP459E	1.00mm
	CSUP460A	CSUP460B	CSUP460C	CSUP460D	CSUP460E	2.00mm

Box of 20 Units

SU BIPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-GREEN



UMAX: 0.5KVP



Irrigation and With Suction Single Use Bayonet Bipolar Forceps – Attached Wire



	Irrigation	With Suction	Description and Length	Wire Length
	CIR-1227	CWS-1327	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1228	CWS-1328	Length: 18.00cm	
	CIR-1229	CWS-1329	Length: 20.00cm	
	CIR-1230	CWS-1330	Length: 22.00cm	
	CIR-1231	CWS-1331	Length: 24.00cm	
	CIR-1227A	CWS-1327A	Length: 5.50cm	3M = 3meter 5M = 5meter
	CIR-1228A	CWS-1328A	Length: 18.00cm	
	CIR-1229A	CWS-1329A	Length: 20.00cm	
	CIR-1230A	CWS-1330A	Length: 22.00cm	
	CIR-1231A	CWS-1331A	Length: 24.00cm	
	CIR-1227B	CWS-1327B	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1228B	CWS-1328B	Length: 18.00cm	
	CIR-1229B	CWS-1329B	Length: 20.00cm	
	CIR-1230B	CWS-1330B	Length: 22.00cm	
	CIR-1231B	CWS-1331B	Length: 24.00cm	
	CIR-1227C	CWS-1327C	Length: 5.50cm	3M = 3meter 5M = 5meter
	CIR-1228C	CWS-1328C	Length: 18.00cm	
	CIR-1229C	CWS-1329C	Length: 20.00cm	
	CIR-1230C	CWS-1330C	Length: 22.00cm	
	CIR-1231C	CWS-1331C	Length: 24.00cm	
	CIR-1227D	CWS-1327D	Length: 5.50cm	3M = 3meter 5M = 5meter
	CIR-1228D	CWS-1328D	Length: 18.00cm	
	CIR-1229D	CWS-1329D	Length: 20.00cm	
	CIR-1230D	CWS-1330D	Length: 22.00cm	
	CIR-1231D	CWS-1331D	Length: 24.00cm	

Importance Note: Always mention tip sizes while preparing order for these items



Irrigation and With Suction Single Use Bayonet Bipolar Forceps Angled – Attached Wire



	Irrigation	With Suction	Description and Length	Wire Length
	CIR-1232	CWS-1332	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1233	CWS-1333	Length: 18.00cm	
	CIR-1234	CWS-1334	Length: 20.00cm	
	CIR-1235	CWS-1335	Length: 22.00cm	
	CIR-1236	CWS-1336	Length: 24.00cm	
	CIR-1232A	CWS-1332A	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1233A	CWS-1333A	Length: 18.00cm	
	CIR-1234A	CWS-1334A	Length: 20.00cm	
	CIR-1235A	CWS-1335A	Length: 22.00cm	
	CIR-1236A	CWS-1336A	Length: 24.00cm	
	CIR-1232B	CWS-1332B	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1233B	CWS-1333B	Length: 18.00cm	
	CIR-1234B	CWS-1334B	Length: 20.00cm	
	CIR-1235B	CWS-1335B	Length: 22.00cm	
	CIR-1236B	CWS-1336B	Length: 24.00cm	
	CIR-1232C	CWS-1332C	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1233C	CWS-1333C	Length: 18.00cm	
	CIR-1234C	CWS-1334C	Length: 20.00cm	
	CIR-1235C	CWS-1335C	Length: 22.00cm	
	CIR-1236C	CWS-1336C	Length: 24.00cm	
	CIR-1232D	CWS-1332D	Length: 15.50cm	3M = 3meter 5M = 5meter
	CIR-1233D	CWS-1333D	Length: 18.00cm	
	CIR-1234D	CWS-1334D	Length: 20.00cm	
	CIR-1235D	CWS-1335D	Length: 22.00cm	
	CIR-1236D	CWS-1336D	Length: 24.00cm	

Importance Note: Always mention tip sizes while preparing order for these items



SU MONPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-BLUE



ETHYLENE OXIDE
STERILIZED



SINGLE USE



MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp

Regular Single Use Adson Monopolar Forceps with Cable-3m

	121mm	155mm	180mm	200mm	220mm	
	MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB461W CSUB462W CSUB463W	CSUB461A CSUB462A CSUB463A	CSUB461B CSUB462B CSUB463B	CSUB461C CSUB462C CSUB463C	CSUB461D CSUB462D CSUB463D
		CSUB464W CSUB465W CSUB466W	CSUB464A CSUB465A CSUB466A	CSUB464B CSUB465B CSUB466B	CSUB464C CSUB465C CSUB466C	CSUB464D CSUB465D CSUB466D
	VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD					

Box of 20 Units

Regular Single Use Narrow serration Dissecting Monopolar Forceps with Cable-3m

155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB473A CSUB474A CSUB475A	CSUB473B CSUB474B CSUB475B	CSUB473C CSUB474C CSUB475C	CSUB473D CSUB474D CSUB475D	CSUB473E CSUB474E CSUB475E
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB476A CSUB477A CSUB478A	CSUB476B CSUB477B CSUB478B	CSUB476C CSUB477C CSUB478C	CSUB476D CSUB477D CSUB478D	CSUB476E CSUB477E CSUB478E

Box of 20 Units

Regular Single Use McIndoe Monopolar Forceps with Cable-3m

155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB467A CSUB468A CSUB469A	CSUB467B CSUB468B CSUB469B	CSUB467C CSUB468C CSUB469C	CSUB467D CSUB468D CSUB469D	CSUB467E CSUB468E CSUB469E
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB470A CSUB471A CSUB472A	CSUB470B CSUB471B CSUB472B	CSUB470C CSUB471C CSUB472C	CSUB470D CSUB471D CSUB472D	CSUB470E CSUB471E CSUB472E

Box of 20 Units

Regular Single Use Wide serration Monopolar Forceps with Cable-3m

155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB479A CSUB480A CSUB481A	CSUB479B CSUB480B CSUB481B	CSUB479C CSUB480C CSUB481C	CSUB479D CSUB480D CSUB481D	CSUB479E CSUB480E CSUB481E
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB482A CSUB483A CSUB484A	CSUB482B CSUB483B CSUB484B	CSUB482C CSUB483C CSUB484C	CSUB482D CSUB483D CSUB484D	CSUB482E CSUB483E CSUB484E

Box of 20 Units

Regular Single Use McIndoe Monopolar Forceps with Cable-3m

155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB467A CSUB468A CSUB469A	CSUB467B CSUB468B CSUB469B	CSUB467C CSUB468C CSUB469C	CSUB467D CSUB468D CSUB469D	CSUB467E CSUB468E CSUB469E
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB470A CSUB471A CSUB472A	CSUB470B CSUB471B CSUB472B	CSUB470C CSUB471C CSUB472C	CSUB470D CSUB471D CSUB472D	CSUB470E CSUB471E CSUB472E

Box of 20 Units

Regular Single Use Bayonet Monopolar Forceps with Cable-3m

155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB485A CSUB486A CSUB487A	CSUB485B CSUB486B CSUB487B	CSUB485C CSUB486C CSUB487C	CSUB485D CSUB486D CSUB487D	CSUB485E CSUB486E CSUB487E
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB488A CSUB489A CSUB490A	CSUB488B CSUB489B CSUB490B	CSUB488C CSUB489C CSUB490C	CSUB488D CSUB489D CSUB490D	CSUB488E CSUB489E CSUB490E

Box of 20 Units

SU MONPOLAR FORCEPS INTEGRATED CABLE

REGULAR INSULATED FORCEPS-BLUE



ETHYLENE OXIDE
STERILIZED



SINGLE USE



MAX VOLTAGE
PROTECTION

UMAX: 0.5KVp



Regular Single Use Bayonet Monopolar Forceps Angled Down tip with Cable-3m

	121mm	155mm	180mm	200mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB491W CSUB492W CSUB493W	CSUB491A CSUB492A CSUB493A	CSUB491B CSUB492B CSUB493B	CSUB491C CSUB492C CSUB493C	CSUB491D CSUB492D CSUB493D	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB494W CSUB495W CSUB496W	CSUB494A CSUB495A CSUB496A	CSUB494B CSUB495B CSUB496B	CSUB494C CSUB495C CSUB496C	CSUB494D CSUB495D CSUB496D	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Bayonet Monopolar Forceps Angled Up tip with Cable-3m

	155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB497A CSUB498A CSUB499A	CSUB497B CSUB498B CSUB499B	CSUB497C CSUB498C CSUB499C	CSUB497D CSUB498D CSUB499D	CSUB497E CSUB498E CSUB499E	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB500A CSUB501A CSUB502A	CSUB500B CSUB501B CSUB502B	CSUB500C CSUB501C CSUB502C	CSUB500D CSUB501D CSUB502D	CSUB500E CSUB501E CSUB502E	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Striaht Monopolar Forceps with Cable-3m

	155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB503A CSUB504A CSUB505A	CSUB503B CSUB504B CSUB505B	CSUB503C CSUB504C CSUB505C	CSUB503D CSUB504D CSUB505D	CSUB503E CSUB504E CSUB505E	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB506A CSUB507A CSUB508A	CSUB506B CSUB507B CSUB508B	CSUB506C CSUB507C CSUB508C	CSUB506D CSUB507D CSUB508D	CSUB506E CSUB507E CSUB508E	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Monopolar Scissor with Cable-3m

	155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB509A CSUB510A CSUB511A	CSUB509B CSUB510B CSUB511B	CSUB509C CSUB510C CSUB511C	CSUB509D CSUB510D CSUB511D	CSUB509E CSUB510E CSUB511E	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB512A CSUB513A CSUB514A	CSUB512B CSUB513B CSUB514B	CSUB512C CSUB513C CSUB514C	CSUB512D CSUB513D CSUB514D	CSUB512E CSUB513E CSUB514E	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Monopolar Artery Forceps with Cable-3m

	155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB515A CSUB516A CSUB517A	CSUB515B CSUB516B CSUB517B	CSUB515C CSUB516C CSUB517C	CSUB515D CSUB516D CSUB517D	CSUB515E CSUB516E CSUB517E	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB518A CSUB519A CSUB520A	CSUB518B CSUB519B CSUB520B	CSUB518C CSUB519C CSUB520C	CSUB518D CSUB519D CSUB520D	CSUB518E CSUB519E CSUB520E	0.50mm 1.00mm 2.00mm

Box of 20 Units

Regular Single Use Regminton Monopolar Forceps with Cable-3m

	155mm	180mm	200mm	240mm	220mm	
MARTIN: BERCHTOLD: ERBE T-SERIES, AESCULAP, SUTTER	CSUB521A CSUB522A CSUB523A	CSUB521B CSUB522B CSUB523B	CSUB521C CSUB522C CSUB523C	CSUB521D CSUB522D CSUB523D	CSUB521E CSUB522E CSUB523E	0.50mm 1.00mm 2.00mm
VALLEYLAB: CONMED: BOVIE 8MM PLUG, BOWA, BARD	CSUB524A CSUB525A CSUB526A	CSUB524B CSUB525B CSUB526B	CSUB524C CSUB525C CSUB526C	CSUB524D CSUB525D CSUB526D	CSUB524E CSUB525E CSUB526E	0.50mm 1.00mm 2.00mm

Box of 20 Units




08

ELECTROSURGICAL GENERATORS

PRECISE. INSULATED. RELIABLE.

ELECTROSURGICAL GENERATORS

C-TECH SMART ELECTROSURGICAL UNITS



UMAX: 0.5KVp

CEG101

Power Supply:

110V, 50/60Hz or 220V, 50/60Hz

Dimension:

Product: 330(L) x 280(W) x 100(H) mm

Package: 440(L) x 390(W) x 260(H) mm

Smart Electrosurgical Unit

Compact design for both monopolar & bipolar.
Max 4 working modes: Pure/Blend, Standard/Bipolar.
Fast trouble-shooting by default Error alarm system
REM (Return Electrode Monitoring) system for highest safety; Excellent operation performance with minimal sticking & carbonization. Different & adjustable operation sounds in case of misoperation; Wide clinical application, Suitable for different surgery

Electrosurgical Unit 150WATTS

Weight:

Net Weight: 4.7KG

Gross Weight: 8.9KG




Easy-to-Use

Convenient Operation by Finger-suit Buttons



Error Code

Simple Error Code for Fast Trouble-shooting



High Safety

Rem (Return Electrode Monitoring) system to avoid any burnt or damage to the surgeon and patient.



Wide Application

It can be connected to endoscopes for more clinical applications.

Standard Accessory

Optional



ELECTROSURGICAL GENERATORS

C-TECH SMART ELECTROSURGICAL UNITS



UMAX: 0.5KVp

CEG102

Power Supply:

110V, 50/60Hz or 220V, 50/60Hz

Dimension:

Product: 421(L) x 357(W) x 162(H) mm
Package: 520(L) x 445(W) x 305(H) mm

Electrosurgical Unit 400 WATTS

Weight:

Net Weight: 8.7KG
Gross Weight: 14.5KG

High-frequency Electrosurgical Unit

Easy-to-use buttons for simple operation
High power output 400W for more operation selection
Up to 7 work modes: Pure Blend 1 Blend 2, Standard Spray, Bipolar 1 Bipolar 2
Increment of 5W micro power control system to ensure optimal effect Minimal sticking and carbonization, reduce the surgery damage
Fast trouble-shooting by default error alarm system
Auto-terminate of power output after 5 seconds to avoid misoperation Support endoscopes connection and underwater cutting (TURP etc) OEM projects



Multiple Functions

Clear Color Indicators for Work



Easy-to-Use


Convenient Operation by Finger-suit Buttons



Error Code

Simple Error Code for Fast Trouble-shooting

Mode & Alarm



High Safety

Rem (Return Electrode Monitoring) system to avoid any burn or damage to the surgeon and patient.



Wide Application

It can be connected to endoscopes for more clinical applications.

Standard Accessory

Optional





COSIL[®]
INSTRUMENTS COMPANY

COSIL[®]
INSTRUMENTS COMPANY

We Provide
24/7 Customer Support

+92 331 732 7777

sales@cosilinstruments.com

www.cosilinstruments.com

CE
0197



EU-MDR ✓
Transit in Process